

# INSTRUCTION MANUAL

## GLAD PCR Assay Kit

Catalog No. K009, K010/1, K010/3, K010/5, K010/32

### Highlights

- Study of DNA methylation in human and mammals
- Determination of presence and quantity of R(5mC)GY site in the point of interest in analyzed DNA by comparison with the calibration line
- Assay based on the primers and TaqMan probe designed on customer's request
- All reagents for control experiments and formation of the calibration line are included
- Assay performed in under 4 hours including DNA hydrolysis, adapter ligation and PCR
- **NO BISULFITE CONVERSION**

### Contents

Product Contents.....	2
Product description.....	3
Protocol.....	4
Appendix 1: Protocol of the control experiments .....	5
Appendix 2: Protocol of the calibration line formation for determination of the number of DNA copies carrying A(5mC)GT site in regulation region of URB1 gene [4] .....	6
Appendix 3: Hybrid primers sequences .....	8
References .....	8
Ordering information.....	9
Services .....	9

### Web page

<http://sibenzyme.com/products/kits>



# Product Contents

GLAD PCR Assay Kit (200 reactions)		K009	K010/1 , K010/3, K010/5, K010/32
1	dH <sub>2</sub> O	4 mL	4 mL
2	10X SE TMN Buffer	150 µL	150 µL
3	DMSO	35 µL	35 µL
4	BSA, 10 mg/mL	70 µL	70 µL
5	GlaI ( 20 u/µL)	10 µL	10 µL
6	Universal adapter, double-stranded (10 µM)	110 µL	110 µL
7	ATP, 10 mM	110 µL	110 µL
8	DNA ligase T4 concentrated (2000 u/µL)	12 µL	12 µL
9	10X SE GLAD-Mg Buffer	430 µL	430 µL
10	MgCl <sub>2</sub> , 50 mM	120 µL	120 µL
11	dNTP Mix, 10 mM each	90 µL	90 µL
12	SP Taq DNA Polymerase, 5 u/µL	40µL	40 µL
13	Control DNA Raji, 18ng/ µL	10 µL	10 µL
14	Control DNA HeLa, 18ng/ µL	10 µL	10 µL
15	Control DNA λ, 18ng/ µL	10 µL	10 µL
16	Control URB1 mix (primers + TaqMan probe), 10 µM each	40 µL	40 µL
17	Control CEBPD mix (primers + TaqMan probe), 10 µM each	15 µL	15 µL
18	Hybrid primers, (10 µM)	—	N vials, 170 µL each

Storage conditions: -20°C

**Note - Integrity of the kit components is guaranteed for one year from a date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide maximal performance and reliability.**

SibEnzyme GLAD PCR Assay kit contains all reagents for analysis of any R(5mC)GY site in human or mammalian genomes except TaqMan Probe and primers. The sequences of these oligonucleotides should correspond to the DNA region nearby the analyzed R(5mC)GY site (see page 4 for protocol).

GLAD PCR Assay of two R(5mC)GY sites in human genome are used as an example of the method application. Mixes of TaqMan Probes and primers for analysis of only these two R(5mC)GY sites are included:

1. "Control URB1 mix (primers + TaqMan probe)" (p.16) is provided for analysis of A(5mC)GT site in regulation region of URB1 gene. Position of this site (according to a human genome assembly GRCh38/hg38) is 32334291-32334294 in chromosome 21.
2. "Control CEBPD mix (primers + TaqMan probe)" (p.17) is provided for analysis of G(5mC)GC site in regulation region of CEBPD gene. Position of this site (according to a human genome assembly GRCh38/hg38) is 47738502-47738505 in chromosome 8.

---

\* You can choose 1, 3, 5 hybrid primers or order full range of 32 hybrid primers from the list in Appendix 3


# Product description

Today an abnormal methylation of regulation regions (promoter and/or first exon) of genes was shown at initial stage of several diseases such as *cancer*, *cardiovascular disease*, *diabetes* and some others.

This abnormal de novo DNA methylation is performed by DNMT3A and DNMT3B DNA methyltransferases. These enzymes recognize and methylate site 5'-RCGY-3' with formation of 5'-R(5mC)GY-3'/3'YG(5mC)R-5'[1].



Recently discovered site-specific methyl-directed DNA-endonuclease Glal recognizes exactly this DNA sequence 5'-R(5mC)GY-3'/3'YG(5mC)R-5' and cleaves it as indicated by symbols [2].



GLAD PCR assay is based on Glal application and includes three steps [3]:

- 1) Glal hydrolysis of studied DNA
- 2) the universal adapter ligation
- 3) subsequent real-time PCR with Taqman probe.


- ❖ *One primer is designed for DNA region of interest, structure of another primer is based on an adapter sequence.*

PCR product is observed only if R(5mC)GY site is present in the point of interest in studied DNA.

## 1. DNA hydrolysis

## 2. Adapter ligation

## 3. Real-time PCR



GLAD PCR assay has been developed to determine minimal quantities of 5'-R(5mC)GY-3'/3'YG(5mC)R-5' site in presence of unmethylated DNA. Such DNA preparations are typical for clinical samples of blood and tissues. GLAD PCR assay allows to determine methylation of RCGY site of interest in human and mammalian genomes.

GLAD PCR assay is performed in one tube, takes about 3-4 hours and determines even several copies of R(5mC)GY site of interest.

In comparison with other methylation detection methods GLAD-PCR has strong advantages:

- ✓ **Simple** – 3 easy steps
- ✓ **Accessible** – requires only real time PCR-machine
- ✓ **Quick** – only 3-4 hours
- ✓ **Sensitive** – detects even several copies of selected R(5mC)GY site



# Appendix 1: Protocol of the control experiments

Below we provide a protocol of GLAD-PCR assay for 6 triplets of control *Lambda DNA*, *HeLa DNA* and *Raji DNA*:

1. Mark 6 reaction tubes (200  $\mu$ L) with letters **A,B,C,D,E,F**; mix 14  $\mu$ L H<sub>2</sub>O and 1  $\mu$ L DNA in each tube (**A, B** - *Lambda DNA*; **C, D** - *HeLa DNA*; **E, F** - *Raji DNA*).
2. Prepare DNA Hydrolysis Mix:  
38.2  $\mu$ L H<sub>2</sub>O + 32.6  $\mu$ L 10X SE TMN Buffer + 6.6  $\mu$ L DMSO + 3.3  $\mu$ L BSA; vortex, briefly centrifuge (VBC)  
Add 0.4  $\mu$ L Glal.  
*\*mix by pipette tip carefully, without bubbles*
3. Add 6.5  $\mu$ L of DNA Hydrolysis Mix to each tube.  
*\*mix by pipette tip carefully, without bubbles*
4. Incubate 30 min at 30 °C.
5. Prepare DNA Ligation Mix:  
33.6  $\mu$ L H<sub>2</sub>O + 9.5  $\mu$ L Adapter + 9.5  $\mu$ L ATP + 1.0  $\mu$ L T4 DNA Ligase; (VBC).
6. Add 8.5  $\mu$ L of DNA Ligation Mix to each tube after Glal hydrolysis (p.4); (VBC).
7. Incubate 15 min at 25 °C.
8. Prepare PCR Mix for URB1 for 3 triplets:  
57  $\mu$ L H<sub>2</sub>O + 19  $\mu$ L 10X SE GLAD Buffer + 4.8  $\mu$ L MgCl<sub>2</sub> + 3.9  $\mu$ L dNTP mix + 1.9  $\mu$ L BSA + 7.5  $\mu$ L URB1 mix + 0.6  $\mu$ L SP Taq DNA pol; (VBC).
9. Prepare PCR Mix for CEBPD for 3 triplets:  
57  $\mu$ L H<sub>2</sub>O + 19  $\mu$ L 10X SE GLAD Buffer + 4.8  $\mu$ L MgCl<sub>2</sub> + 3.9  $\mu$ L dNTP + 1.9  $\mu$ L BSA + 7.5  $\mu$ L CEBPD mix + 0.6  $\mu$ L SP Taq DNA pol; (VBC).
10. Add 30  $\mu$ L of PCR Mix for URB1 to tube A, to tube C and to tube E; (VBC).
11. Add 30  $\mu$ L of PCR Mix for CEBPD to tube B, to tube D and to tube F; (VBC).
12. Withdraw 20  $\mu$ L of the reaction mixture 3 times (triplet) from each tube and place them in PCR plate (totally 18 samples), SEAL, CENTRIFUGE, PLACE INTO THERMOCYCLER.
13. Perform real-time PCR: in case of CFX-96 Thermal Cycler Amplifier" (Bio-Rad, USA) the amplification profile is as follows:
  - 3 min at 95 °C;
  - 5 "blind" cycles:
  - 95 °C for 10 sec;
  - 63 °C for 20 sec ;
  - 72 °C for 5 sec;
  - 40 cycles:
  - 95 °C for 10 sec;
  - 63 °C for 20 sec (with fluorescence detection in FAM channel);
  - 72°C for 5 sec.

Amplification profile may differ for another PCR amplification system.

## Expected results of control experiment



Region \ DNA	<i>Raji DNA</i>	<i>Hela DNA</i>	<i>Lambda DNA</i> (negative control)
URB1	(+) Methylated	(+) Methylated	(-)
CEBPD	(+) Methylated	(-) Unmethylated	(-)

## Appendix 2: Protocol of the calibration line formation for determination of the number of DNA copies carrying A(5mC)GT site in regulation region of URB1 gene [4]

Below we provide a protocol of GLAD-PCR assay for determination of the number of DNA copies carrying A(5mC)GT site in regulation region of URB1 gene. One tube with control Lambda DNA and 5 tubes with dilutions of Raji DNA have been used to obtain the calibration line.

1. Mark 6 reaction tubes (200  $\mu$ L) with numbers 1-6.  
Add 28  $\mu$ L H<sub>2</sub>O to the vial 1, add 15  $\mu$ L H<sub>2</sub>O to vials 2, 3, 4, 5 and add 14  $\mu$ L H<sub>2</sub>O to the vial 6.  
Add 2  $\mu$ L Raji DNA to the vial 1, mix carefully.  
Withdraw 15  $\mu$ L of DNA solution from vial 1 and add this aliquot to vial 2, mix carefully.  
Withdraw 15  $\mu$ L of DNA solution from vial 2 and add this aliquot to vial 3, mix carefully.  
Withdraw 15  $\mu$ L of DNA solution from vial 3 and add this aliquot to vial 4, mix carefully.  
Withdraw 15  $\mu$ L of DNA solution from vial 4 and add this aliquot to vial 5, mix carefully.  
Withdraw 15  $\mu$ L of DNA solution from vial 5 and discard this aliquot. Add 1  $\mu$ L Lambda DNA to the vial 6, mix carefully, briefly centrifuge.
2. Prepare DNA Hydrolysis Mix:  
38.2  $\mu$ L H<sub>2</sub>O + 32.6  $\mu$ L 10X SE TMN Buffer + 6.6  $\mu$ L DMSO + 3.3  $\mu$ L BSA; after each addition mix carefully, briefly centrifuge. Add 0.4  $\mu$ L Glal, mix by the pipette tip carefully, without bubbles, briefly centrifuge.
3. Add 6.5  $\mu$ L of DNA Hydrolysis Mix to each tube, mix by pipette tip carefully, without bubbles, briefly centrifuge.
4. Incubate 30 min at 30 °C.
5. Prepare DNA Ligation Mix:  
33.6  $\mu$ L H<sub>2</sub>O + 9.5  $\mu$ L Adapter + 9.5  $\mu$ L ATP + 1.0  $\mu$ L T4 DNA Ligase;  
after each addition mix carefully, briefly centrifuge.
6. Add 8.5  $\mu$ L DNA Ligation Mix to each tube after Glal hydrolysis (p.4);  
mix carefully, briefly centrifuge.
7. Incubate 15 min at 25 °C.
8. Prepare PCR Mix for A(5mC)GT site in regulation region of URB1 gene:  
114  $\mu$ L H<sub>2</sub>O + 38  $\mu$ L 10X SE GLAD Buffer + 9.6  $\mu$ L MgCl<sub>2</sub> + 7.8  $\mu$ L dNTP mix + 3.8  $\mu$ L BSA +  
15  $\mu$ L URB1 mix+ 1,2  $\mu$ L SP Taq DNA pol;  
after each addition mix carefully, briefly centrifuge.
9. Add 30  $\mu$ L of PCR Mix for URB1 to tubes 1, 2, 3, 4, 5 and 6.  
Mix carefully, briefly centrifuge.
10. Withdraw 20  $\mu$ L of the reaction mixture 3 times (triplet) from each tube and place them in PCR plate (totally 18 samples), SEAL, CENTRIFUGE, PLACE INTO THERMOCYCLER.
11. Perform real-time PCR: in case of CFX-96 Thermal Cycler Amplifier" (Bio-Rad, USA) the amplification profile is as follows: 3 min at 95 °C; 5 "blind" cycles: 95 °C for 10 sec; 63 °C for 20 sec; 72 °C for 5 sec; 40 cycles: 95 °C for 10 sec; 63 °C for 20 sec (with fluorescence detection in FAM channel); 72°C for 5 sec. Amplification profile may differ for another PCR amplification system.

## Expected results of experiment



## Appendix 3: Hybrid primers sequences

HP1: 5' CCTGCTCTTTCATCGGCAA 3' (19)	HP17 5' CCTGCTCTTTCATCGGTAA 3' (19)
HP2: 5' CCTGCTCTTTCATCGGCAC 3' (19)	HP18 5' CCTGCTCTTTCATCGGTAC 3' (19)
HP3: 5' CCTGCTCTTTCATCGGCAG 3' (19)	HP19 5' CCTGCTCTTTCATCGGTAG 3' (19)
HP4: 5' CCTGCTCTTTCATCGGCAT 3' (19)	HP20 5' CCTGCTCTTTCATCGGTAT 3' (19)
HP5: 5' CCTGCTCTTTCATCGGCCA 3' (19)	HP21 5' CCTGCTCTTTCATCGGTCA 3' (19)
HP6: 5' CCTGCTCTTTCATCGGCC 3' (19)	HP22 5' CCTGCTCTTTCATCGGTCC 3' (19)
HP7: 5' CCTGCTCTTTCATCGGCCG 3' (19)	HP23 5' CCTGCTCTTTCATCGGTCCG 3' (19)
HP8: 5' CCTGCTCTTTCATCGGCCT 3' (19)	HP24 5' CCTGCTCTTTCATCGGTCT 3' (19)
HP9: 5' CCTGCTCTTTCATCGGCCGA 3' (19)	HP25 5' CCTGCTCTTTCATCGGTGA 3' (19)
HP10 5' CCTGCTCTTTCATCGGCCG 3' (19)	HP26 5' CCTGCTCTTTCATCGGTGC 3' (19)
HP11 5' CCTGCTCTTTCATCGGCCGG 3' (19)	HP27 5' CCTGCTCTTTCATCGGTGG 3' (19)
HP12 5' CCTGCTCTTTCATCGGCCGT 3' (19)	HP28 5' CCTGCTCTTTCATCGGTGT 3' (19)
HP13 5' CCTGCTCTTTCATCGGCTA 3' (19)	HP29 5' CCTGCTCTTTCATCGGTTA 3' (19)
HP14 5' CCTGCTCTTTCATCGGCTC 3' (19)	HP30 5' CCTGCTCTTTCATCGGTTTC 3' (19)
HP15 5' CCTGCTCTTTCATCGGCTG 3' (19)	HP31 5' CCTGCTCTTTCATCGGTTG 3' (19)
HP16 5' CCTGCTCTTTCATCGGCTT 3' (19)	HP32 5' CCTGCTCTTTCATCGGTTT 3' (19)

## References

1. Handa, V., and Jeltsch A. "Profound sequence preference of Dnmt3a and Dnmt3b mammalian DNA methyltransferases shape the human epigenome" // 2005, J. Mol. Biol. 348, 1103-1112
2. Tarasova G. V., Nayakshina T. N., Degtyarev S. Kh. "Substrate specificity of new methyl-directed DNA endonuclease Glal" // BMC Molecular Biology 2008, 9:7
3. Kuznetsov V.V., Abdurashitov M.A., Akishev A.G., Degtyarev S.Kh. "Method of determining nucleotide sequence Pu(5mC)GPy at predetermined position of continuous DNA" // Patent RU 2525710
4. A.G. Akishev et al. "GLAD-PCR assay of selected R(5mC)GY sites in URB1 and CEPBD genes in human genome" // Res J Pharm Biol Chem Sci, 8(1): pp.465-475, 2017



## Ordering information

Product description	Catalog No.
GLAD PCR Assay Kit, 200 rxns	K009
GLAD PCR Assay Kit with hybrid primers, 200 rxns, 1 hybrid primer	K010/1
GLAD PCR Assay Kit with hybrid primers, 200 rxns, 3 hybrid primers	K010/3
GLAD PCR Assay Kit with hybrid primers, 200 rxns, 5 hybrid primers	K010/5
GLAD PCR Assay Kit with hybrid primers, 200 rxns, 32 hybrid primers	K010/32

For individual sale	Catalog No.	Amount(s)
GlaI (20 units/ $\mu$ L)	E493/ E494	100/500 units
DNA ligase T4 concentrated (2000 units/ $\mu$ L)	E329/E330	50000/250000 units
10xSE GLAD Buffer	B022	1 mL
dNTP Mix, 10 mM each	N025	40 $\mu$ mol
SP Taq DNA Polymerase (5 units/ $\mu$ L)	E333/E334	200/1000 units
TaqMan probe synthesis	K009TM	1 set for 200 rxns
Taq Man (1) + Primers (2) synthesis	K009PTM	1 set for 200 rxns

## Services

Please get in touch with us ([info@sibenzyme.com](mailto:info@sibenzyme.com)) if you need our assistance with primers design for your R(5mC)GY site of interest.

Primers mix (primers + TaqMan probe) for analysis of your R(5mC)GY site of interest may be synthesized by customer or ordered at SibEnzyme separately.