

About us.

SibEnzyme Ltd. (Siberian Enzyme, SE),
is a privately owned Russian company founded in 1991.
SibEnzyme Ltd. is located in Academtown near Novosibirsk city – the capital of
Siberia and a geographic center of Russia.

SibEnzyme's primary focus is production of enzymes and related products for
molecular biology, PCR and genetic engineering. SE product line includes more than
200 enzymes, a set of perfect DNA ladders,
high quality enzymatic dNTPs and DNA preparations.

SibEnzyme Ltd. is one of the leading companies in research, development and
production of restriction endonucleases and new DNA endonucleases.

A new type of methyl-directed DNA endonucleases
has been discovered and characterized in our laboratories.
These new enzymes are a good instrument to detect DNA methylation status
in epigenetics studies and medical applications.

Information about SibEnzyme products may be found at
www.sibenzyme.com

Scientific publications with detailed study of new enzymes are available at
www.science.sibenzyme.com



New product



Mammalian Genome Qualified



New package

Contents

| | |
|---|----|
| Restriction Endonucleases..... | 3 |
| Nickases | 51 |
| Methyl-directed DNA endonucleases | 52 |
| DNA – Methyltransferases..... | 59 |
| Polymerases..... | 60 |
| DNA Modifying Enzymes | 64 |
| DNA Ladders | 67 |
| Plasmid and Phage DNAs | 71 |
| Human Genomic DNA..... | 72 |
| dNTPs(enzymatic)..... | 73 |
| dNTPs(chemical)..... | 75 |
| Substrates | 77 |
| Kits | 78 |
| Buffers composition | 80 |
| SE Buffer Activity Chart for Restriction Enzymes..... | 82 |
| Isoschizomers..... | 86 |
| Alphabetized List of SE Recognition Sequences..... | 88 |

Restriction Endonucleases

Restriction endonucleases are bacterial site-specific DNA-endonucleases that recognize and cleave with high specificity short (usually 4-8 base pairs) sequences in double stranded DNA. Type II restriction endonucleases cut both strands of substrate DNA in the presence of Mg^{2+} ions near or within the recognition sequence.

SibEnzyme Ltd. is one of the leading companies in research, development and production of restriction endonucleases. At present SibEnzyme Ltd. has **more than 200 commercially available restriction endonucleases** in the product list including many unique enzymes. 56 of these enzymes have been cloned and are purified from recombinant sources.

To simplify the reaction of DNA hydrolysis with restriction endonucleases and to set up **double digestion** reaction, a special reaction universal **Buffer “ROSE”** was designed. Over 150 enzymes reveal high activity in this buffer. The information about activity of restriction endonucleases in “ROSE” buffer is available in section “Info/SE-buffers activity”.

Fast digestion of different types of DNA (mainly plasmid DNA) may be performed with **Turbo qualified Restriction Enzymes** (marked with symbol **T**, ex: E003T/E004T). Currently **25** Turbo SE restriction enzymes are commercially available. They may be used for a short time (10-15 min) DNA digestion as well as in a standard reaction. The reaction with most of these enzymes may be performed using optimal or universal (“**ROSE**”) Buffer (both are supplied with Turbo enzymes).

Furthermore, to allow the performance of a limited number of reactions (up to 50) with minimum expenses, SE 20 most commonly used restriction endonucleases are now supplied in lesser quantities (marked with a symbol “m”, ex: E003m). The symbol indicates that these enzymes are available in **minimal** packages at optimal prices.

Restriction Endonucleases: Quality Controls

Unit Determination

One unit of restriction endonuclease is defined as the amount of enzyme required to digest 1 μ g of substrate DNA in a total reaction volume of 50 μ l in 1 hour using the optimal SE-Buffer provided.

Quality Controls

The results of all quality control assays are reported on the Certificate of analysis provided with each enzyme.

Ligation of DNA fragments

DNA fragments are produced by an excessive over-digestion of substrate DNA with each restriction endonuclease. These fragments are then ligated with T4 DNA Ligase at a 5' termini concentration of 0.1-1.0 μ M. The ligated fragments are then recut with the same restriction endonuclease. Ligation can only occur if the 3' and 5' termini are left intact, and only those molecules with a perfectly restored recognition site can be recleaved. A normal banding pattern after cleavage indicates that both the 3' and 5' termini are intact and the enzyme preparation is free of detectable exonucleases and phosphatases. An example of FbII ligation is presented in the figure.

16-hour assay for nuclease contamination

All SibEnzyme restriction endonucleases are incubated for 16 hours in optimal buffer with 1 μ g of substrate DNA in a volume of 50 μ l. The characteristic DNA fragments pattern produced by the enzyme in 1 hour is compared to the pattern produced from an excess of enzyme incubated for 16 hours. A sharp, unaltered pattern under these conditions is an indication that the enzyme preparation is free of detectable levels of nonspecific DNAses. The maximum number of units which can be incubated for 16 hours is indicated.

Assay for exonuclease and phosphatase contamination

All restriction endonucleases are incubated for 3 hours with 5'-[³²P]-labeled synthetic oligonucleotides (single-stranded and double-stranded) in a volume of 20 μ l. After incubation of these labeled oligonucleotides with an enzyme, denatured reaction products are separated on a polyacrylamide gel and then analysed by phospho-imaging. No detectable degradation of single-stranded and double-stranded oligonucleotides indicates that the enzyme preparation is free of exonuclease and phosphatase contamination.

Turbo qualified Restriction Endonucleases

Turbo qualified Restriction Endonucleases may be used for 10-15 minutes DNA digestion as well as in a standard reaction. The reaction with the most of these enzymes may be performed using optimal or universal (“ROSE”) Buffer (both are supplied with Turbo enzymes).

| No | Turbo enzyme | Prototype | Sequence | Cat.No | Digestion time (plasmid, min.) | | | Recommended buffer |
|----|--------------|--------------|-----------|--------------|--------------------------------|----|----|--------------------|
| | | | | | 5 | 10 | 15 | |
| 1 | Acc65I | Acc65I, KpnI | G↑GTACC | E003T/ E004T | | + | + | ROSE |
| 2 | AhlI | SpeI | A↑CTAGT | E173T/ E174T | | + | + | ROSE |
| 3 | AluI | AluI | AG↑CT | E015T/ E016T | | + | + | ROSE |
| 4 | ApaI | ApaI | GGGCC↑C | E019T/ E020T | | + | + | ROSE |
| 5 | AsiGI | AgeI | A↑CCGGT | E235T/ E236T | | + | + | ROSE |
| 6 | BamHI | BamHI | G↑GATCC | E021T/ E022T | + | + | + | ROSE |
| 7 | BglII | BglII | A↑GATCT | E027T/ E028T | | + | + | ROSE |
| 8 | Bsa29I | ClaI | AT↑CGAT | E205T/ E206T | + | + | + | ROSE |
| 9 | Bsp19I | NcoI | C↑CATGG | E047T/ E048T | + | + | + | ROSE |
| 10 | CciNI | NotI | GC↑GGCCGC | E203T/ E204T | | | + | ROSE |
| 11 | EcoRI | EcoRI | G↑AATTC | E057T/ E058T | + | + | + | ROSE |
| 12 | EcoRV | EcoRV | GAT↑ATC | E059T/ E060T | + | + | + | ROSE |
| 13 | FauNDI | NdeI | CA↑TATG | E009T/ E010T | + | + | + | ROSE |
| 14 | HindIII | HindIII | A↑AGCTT | E073T/ E074T | + | + | + | ROSE |
| 15 | HinFI | HinFI | G↑ANTC | E075T/ E076T | + | + | + | ROSE |
| 16 | HpaI | HpaI | GTT↑AAC | E077T/ E078T | | + | + | ROSE |
| 17 | MluI | MluI | A↑CGCGT | E085T/ E086T | | + | + | ROSE |
| 18 | PceI | StuI | AGG↑CCT | E105T/ E106T | | + | + | ROSE |
| 19 | Psp124BI | SacI | GAGCT↑C | E107T/ E108T | + | + | + | ROSE |
| 20 | PstI | PstI | CTGCA↑G | E109T/ E110T | + | + | + | ROSE |
| 21 | SaII | SaII | G↑TCGAC | E115T/ E116T | | + | + | O |
| 22 | Sfr274I | XhoI | C↑TCGAG | E125T/ E126T | | + | + | ROSE |
| 23 | SmaI | SmaI | CCC↑GGG | E177T/ E178T | | + | + | ROSE |
| 24 | SphI | SphI | GCATG↑C | E129T/ E130T | | + | + | ROSE |
| 25 | XbaI | XbaI | T↑CTAGA | E141T/ E142T | + | + | + | ROSE |

Minimal packages at optimal price.

| No | Cat. No | Enzyme | Qty u/pack. | Price EUR/ 1 pack |
|----|---------|-----------|-------------|----------------------|
| 1 | E021m | BamHI | 500 | 2 |
| 2 | E057m | EcoR I | 500 | 2 |
| 3 | E059m | EcoR V | 500 | 2 |
| 4 | E067m | Hae III | 500 | 2 |
| 5 | E073m | Hind III | 500 | 2 |
| 6 | E201m | Hind II | 250 | 2 |
| 7 | E075m | Hinf I | 500 | 2 |
| 8 | E161m | Hpa II | 250 | 2 |
| 9 | E109m | Pst I | 500 | 2 |
| 10 | E125m | Sfr274 I | 500 | 2 |
| No | Cat. No | Enzyme | Qty u/pack. | Price EUR/ 1 pack |
| 11 | E003m | Acc65 I | 500 | 5 |
| 12 | E019m | Apa I | 500 | 5 |
| 13 | E027m | Bgl II | 500 | 5 |
| 14 | E009m | FauND I | 500 | 5 |
| 15 | E077m | Hpa I | 250 | 5 |
| 16 | E069m | HspA I | 500 | 5 |
| 17 | E107m | Psp124B I | 500 | 5 |
| 18 | E177m | Sma I | 500 | 5 |
| 19 | E139m | Vsp I | 500 | 5 |
| 20 | E141m | Xba I | 500 | 5 |

| | | | | | |
|--|---|--|-------------------------------|-------|-----|
| Aat II (prototype Aat II) Isolated from an <i>E.coli</i> strain that carries the cloned Aat II gene from <i>Acetobacter acetii</i> | GACGT[^]C C[^]TGCAG | E287 E288 | 500 u.a. 2500 u.a. | | |
| Concentration: 10 000 units/ml Assayed on λ DNA Reagents Supplied with Enzyme: 10×SE Buffer Y Reaction conditions: 1×SE Buffer Y Incubate at 37°C. | Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM NaCl; 0,1 mM EDTA; 200 µg/ml BSA; 1 mM DTT; and 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 10-fold overdigestion with Aat II, about 90% of the DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of Aat II for 16 hours. Star activity: High enzyme concentration may result in star activity. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 10-25 | 25-50 | 10-25 | 25-50 | 100 |

| | | | | | |
|--|---|---|-----------------------------|-------|------|
| Abs I (prototype Abs I) Isolated from <i>Arthrobacter species</i> 7M06 | CC[^]TCGAGG GGAGCT[^]CC | E535 E536 | 50 u.a. 250 u.a. | | |
| Concentration: 500 - 1 000 units/ml Assayed on pUC19SE/DriI digest Reagents Supplied with Enzyme: 10×SE Buffer AbsI Reaction conditions: 1× SE Buffer AbsI Incubate at 37°C. | Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 0.05% Triton X-100 and 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 2-fold overdigestion with Abs I, ~90% of the DNA fragments can be ligated with T4 DNA Ligase and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of DNA with 2 units of Abs I for 16 hours. Star activity: A long incubation time may result in star activity. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 10-25 | 0 | 50-75 | 0-10 |

| | | | | | |
|---|---|--|-------------------------------|-----|--------|
| Acc16 I (prototype Mst I) Isolated from <i>Acinetobacter calcoaceticus</i> 16 | TGC[^]GCA ACG[^]CGT | E001 E002 | 200 u.a. 1000 u.a. | | |
| Concentration: 5 000-10 000 units/ml Assayed on λ DNA Reagents Supplied with Enzyme: 10×SE Buffer W Reaction conditions: 1× SE Buffer W Incubate at 37°C. | Storage buffer: 10 mM Tris-HCl (pH 7.5); 200 mM KCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; and 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 5-fold overdigestion with Acc16 I, ~90% of the DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of Acc16 I for 16 hours. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 50-75 | 75-100 | 25-50 | 100 | 75-100 |

| | | | | | |
|--|---|---|------------------------------|-------|-----|
| Acc36 I (prototype BspM I) Isolated from <i>Acinetobacter calcoaceticus</i> 36 | ACCTGC(N)₄[^] TGGACG(N)₈[^] | E289 E290 | 100 u.a. 500 u.a. | | |
| Concentration: 2000-5000 units/ml Assayed on λ DNA Reagents Supplied with Enzyme: 10×SE Buffer Y Reaction conditions: 1× SE Buffer Y Incubate at 37°C. | Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM KCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 100 µg/ml BSA; and 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 3-fold overdigestion with Acc36I, 90% of the DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 3 units of Acc36 I for 16 hours. Star activity: High enzyme concentration may result in star activity. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 25-50 | 25-50 | 50-75 | 50-75 | 100 |

| | | | | | |
|--|--|---|--------------------------------|-----|-------|
| Acc65 I (prototype Kpn I) Isolated from <i>Acinetobacter calcoaceticus</i> 65 | G[^]GTACC CCATG[^]G | E003 E004 | 1000 u.a. 5000 u.a. | | |
| Concentration: 10 000 - 30 000 units/ml Assayed on λ DNA (dcm ⁻) Reagents Supplied with Enzyme: 10×SE Buffer W Reaction conditions: 1× SE Buffer W Incubate at 37°C. Blocked by overlapping dcm-methylation (C^mCWGG): GGTACCWGG Acc65I is a neoschizomer of KpnI. | Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; and 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 10-fold overdigestion with Acc65 I, >90% of the DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of Acc65 I for 16 hours. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 10-25 | 25-50 | 75-100 | 100 | 10-25 |

| | | | | | |
|--|---|--|-------------------------------|--------|-------|
| AccB1 I (prototype HgiC I) Isolated from <i>Acinetobacter calcoaceticus</i> B1 | G[^]GYRCC CCRYG[^]G | E163 E164 | 500 u.a. 2500 u.a. | | |
| Concentration: 5 000 -10 000 units/ml Assayed on λ DNA Reagents Supplied with Enzyme: 10×SE Buffer K, BSA Reaction conditions: 1× SE Buffer K +BSA Incubate at 37°C. <i>To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.</i> | Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 10-fold overdigestion with AccB1 I, 95% of the DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of AccB1 I for 16 hours. Star activity: High enzyme concentration may result in star activity. Do not use BSA for long incubation. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 50-75 | 10-25 | 10-25 | 75-100 | 50-75 |

| | | | | | |
|--|---|--|-------------------------------|-------|-------|
| AccB7 I (prototype PflM I) Isolated from <i>Acinetobacter calcoaceticus</i> B7 | CCANNN[^]NTGG GGTN[^]NNNACC | E179 E180 | 200 u.a. 1000 u.a. | | |
| Concentration: 5 000 units/ml Assayed on λ DNA (dcm ⁻) Reagents supplied with Enzyme: 10×SE Buffer G Reaction conditions: 1× SE Buffer G Incubate at 37°C. Blocked by overlapping <i>dcm</i> -methylation (C ^m CWGG): CCANN<u>NCCTGG</u> or CCAGG<u>NNNTGG</u> | Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA, 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 5-fold overdigestion with AccB7 I, 95% of the DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 5 units of AccB7 I for 16 hours. Star activity: High enzyme concentration may result in star activity. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 10-25 | 100 | 25-50 | 50-75 | 50-75 |

| | | | | | |
|---|---|---|--------------------------------|-------|-----|
| AccBS I (prototype BsrB I) Isolated from <i>Acinetobacter calcoaceticus</i> BS | GAG[^]CGG CTC[^]GCC | E007 E008 | 1000 u.a. 5000 u.a. | | |
| Concentration: 5000 – 20 000 units/ml Assayed on: λ DNA Reagents supplied with Enzyme: 10× SE Buffer Y Reaction conditions: 1× SE Buffer Y Incubate at 37°C. | Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 5-fold overdigestion with AccBS I, 90% of the DNA fragments can be ligated. Of these, 50% can be recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of AccBS I for 16 hours. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 25-50 | 25-50 | 100 |

| | | | | | |
|---|--|---|-------------------------------|------|-----|
| Acl I (prototype Acl I) Isolated from <i>Acinetobacter calcoaceticus</i> | AA[^]CGTT TTGC[^]AA | E011 E012 | 200 u.a. 1000 u.a. | | |
| Concentration: 3 000 – 5 000 units/ml Assayed on: λ DNA Reagents Supplied with Enzyme: 10×SE Buffer Y, BSA Reaction condition: 1xSE Buffer Y + BSA Incubate at 37°C. Blocked by CpG methylation. <i>To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.</i> | Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 0,05% Triton X-100; 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 2-fold overdigestion with Acl I, 90% of the DNA fragments can be ligated and recut with this enzyme. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 2 units of Acl I for 16 hours. Do not use BSA for long incubation. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 0-10 | 0-10 | 0-10 | 0-10 | 100 |

| | | | | | |
|--|--|--|------------------------------|------|-----|
| AclW I (prototype Bin I) Isolated from <i>Acinetobacter calcoaceticus</i> W2131 | GGATC(N)₄[^] CCTAG(N)₅[^] | E211 E212 | 100 u.a. 500 u.a. | | |
| Concentration: 1000 - 3000 units/ml Assayed on λ DNA (dam ⁻) Reagents Supplied with Enzyme: 10×SE Buffer Y, BSA Reaction conditions: 1×SE Buffer Y+BSA Incubate at 37°C. Blocked by overlapping <i>dam</i> -methylation (G ^m ATC): GG^mATC <i>To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.</i> | Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM KCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; and 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 2-fold overdigestion with AclW I, about 50% of the DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 6 units of AclW I for 16 hours. Do not use BSA for long incubation. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 50-75 | 0-10 | 0-10 | 100 |

| | | | | | |
|---|--|--|------------------------------|-------|--------|
| Aco I (prototype Cfr I) Isolated from <i>Acinetobacter calcoaceticus</i> | Y[^]GGCCR RCCGG[^]Y | E499 E500 | 100 u.a. 500 u.a. | | |
| Concentration: 500 – 2 000 units/ml Assayed on: λ DNA(dam ⁻ ,dcm ⁻) Reagents Supplied with Enzyme: 10×SE Buffer G Reaction condition: 1xSE Buffer G Incubate at 37°C. Blocked by overlapping <i>dcm</i> -methylation (C ^m CWGG): CCTGGCCR. | Storage buffer: 10 mM Tris-HCl (pH 7.5); 200 mM NaCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 100 µg/ml BSA; and 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 3-fold overdigestion with Aco I, >90% of λ DNA fragments can be ligated with T4 DNA Ligase at 16°C and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 2 units of Aco I for 16 hours. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 50-75 | 100 | 50-75 | 25-50 | 75-100 |

| | | | | | |
|---|---|---|-------------------------------|-----|-------|
| Acs I (prototype Apo I) Isolated from <i>Arthrobacter citreus</i> | R[^]AATTY YTAA[^]R | E013 E014 | 500 u.a. 2500 u.a. | | |
| Concentration: 10 000 – 20 000 units/ml Assayed on λ DNA Reagents Supplied with Enzyme: 10×SE Buffer W, BSA Reaction condition: 1xSE Buffer W+BSA Incubate at 50°C. <i>To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.</i> | Storage buffer: 20 mM Tris-HCl (pH 7.5); 50 mM KCl; 0,1 mM EDTA; 10 mM 2-mercaptoethanol; 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: Yes (80° C for 20 minutes) | Ligation/recutting assay: After 20-fold overdigestion with Acs I, >95% of the DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 50 units of Acs I for 16 hours. Do not use BSA for long incubation. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 25-50 | 50-75 | 50-75 | 100 | 10-25 |

| | | | | | |
|---|--|--|-----------------------------|--------|--------|
| Acu I (prototype Eco57 I) Isolated from an <i>E.coli</i> strain that carries the cloned Acu I gene from <i>Acinetobacter calcoaceticus</i> SRW4 | CTGAAG(N)₁₆[^] GACTTC(N)₁₄[^] | E451 E452 | 50 u.a. 250 u.a. | | |
| Concentration: 1 000 units/ml Assayed on λ DNA Reagents Supplied with Enzyme: 10×SE BufferY, BSA, SAM Reaction conditions: 1×SE Buffer Y+BSA+SAM Incubate at 37°C. <i>To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml, and SAM should be added to a final concentration 0.01 mM.</i> | Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0,1 mM EDTA; 1 mM DTT; 200 µg/ml BSA; 50% glycerol. Store at -20° C. Diluent: SE Buffer A Heat inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 2-fold overdigestion with Acu I, about 80% of the DNA fragments can be ligated. Of these, 80% can be recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 2 units of Acu I for 16 hours. Star activity: High enzyme concentration may result in star activity. Do not use BSA for long incubation. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 25-50 | 50-75 | 50-75 | 75-100 | 75-100 |

| | | | | | |
|---|--|---|-------------------------------|-------------------------------|------------------|
| Afe I (prototype Eco47 III) Isolated from an <i>E.coli</i> strain that carries the cloned Afe I gene from <i>Alcaligenes faecalis</i> T2774 | AGC[^]GCT TCG[^]CGA | E213 E214 | 200 u.a. 1000 u.a. | | |
| Concentration: 10 000 and 50 000 units/ml Assayed on λ DNA (BamHI-digest) Reagents Supplied with Enzyme: 10×SE Buffer Y Reaction conditions: 1×SE Buffer Y Incubate at 37°C. | Storage buffer: 10 mM Tris-HCl (pH 7.6); 50 mM NaCl; 0,1 mM EDTA; 1 mM DTT; 200 µg/ml BSA; and 50% glycerol. Store at -20°C Diluent: SE Buffer A Heat inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 10-fold overdigestion with Afe I, > 80% of pBR322 DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of the λ DNA with 40 units of Afe I for 16 hours. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 10-25 | 25-50 | 75-100 | 75-100 | 100 |
| | | | | For high concentration | |
| | | | | E214X | 1000 u.a. |

| | | | | | |
|---|---|--|-------------------------------|-------|-----|
| Ags I (prototype Ags I) Isolated from <i>Agrococcus species 25</i> | TTS[^]AA AA[^]STT | E573 E574 | 200 u.a. 1000 u.a. | | |
| Concentration: 5000 units/ml Assayed on λ DNA Reagents Supplied with Enzyme: 10×SE Buffer Y, BSA Reaction conditions: 1×SE Buffer Y + BSA Incubate at 37°C. <i>To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.</i> | Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM KCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 5-fold overdigestion with AgsI, >90% of the DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of Ags I for 16 hours. Do not use BSA for long incubation | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 50-75 | 10-25 | 10-25 | 100 |

Ahl I (prototype Spe I) A[^]CTAGT
TGATC[^]A E173
E174 1000 u.a.
5000 u.a.

Isolated from *Alteromonas haloplanktis* SP

Concentration: 10 000 - 30 000 units/ml
Assayed on T7 DNA
Reagents Supplied with Enzyme:
 10×SE Buffer B, BSA
Reaction conditions: 1×SE Buffer B+BSA
 Incubate at 37°C.
To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM KCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA, 50% glycerol. Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
 No (65°C, 80°C for 20 minutes)

Ligation/recutting assay: After 20-fold overdigestion with Ahl I, >90% of T7 DNA fragments can be ligated and recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of T7 DNA with 40 units of Ahl I for 16 hours.
Do not use BSA for long incubation.

| | | | | | |
|-----------------------------------|-----|--------|-------|-------|--------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 100 | 75-100 | 25-50 | 25-50 | 75-100 |

Ajn I (prototype EcoR II) A[^]CCWGG
GGWCC[^] E473
E474 200 u.a.
1000 u.a.


Isolated from *Acinetobacter johnsonii* R2

Concentration: 500 – 3 000 units/ml
Assayed on λ DNA
Reagents Supplied with Enzyme:
 10×SE Buffer Y
Reaction conditions: 1×SE Buffer Y
 Incubate at 55°C.
Note: At 37°C activity is about 10% from maximum.
Not blocked by overlapping *dcm*-methylation (C^mCWGG): **CCWGG**

Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol. Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
 Yes (65° C for 20 minutes)

Ligation/recutting assay: After 3-fold overdigestion with AjnI, 95% of the DNA fragments can be ligated and recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 6 units of Ajn I for 16 hours.

| | | | | | |
|-----------------------------------|-------|-------|-------|-------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 25-50 | 10-25 | 10-25 | 25-50 | 100 |

Alu I (prototype Alu I)  AG[^]CT
TC[^]GA E015
E016 200 u.a.
1000 u.a.

Isolated from an *E.coli* strain that carries the cloned AluI gene from *Arthrobacter luteus*

Concentration: 2 000 – 5000 units/ml
Assayed on λ DNA
Reagents Supplied with Enzyme:
 10×SE Buffer Y
Reaction conditions: 1×SE Buffer Y
 Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 200 mM KCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol. Store at -20°C
Diluent: SE Buffer A
Heat inactivation:
 Yes (65° C for 20 minutes)

Ligation/recutting assay: After 3-fold overdigestion with Alu I, >90% of the DNA fragments can be ligated and recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 4 units of Alu I for 16 hours.

| | | | | | |
|-----------------------------------|--------|--------|-------|-------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 10-25 | 50-75 | 100 |

Qualified for human genome studies: http://science.sibenzyme.com/article8_article_31_1.phtml

AluB I (prototype Alu I) AG[^]CT
TC[^]GA E549
E550 200 u.a.
1000 u.a.

Isolated from *Arthrobacter luteus* B

Concentration: 5 000 units/ml
Assayed on λ DNA
Reagents Supplied with Enzyme:
 10×SE Buffer B, BSA
Reaction conditions: 1×SE Buffer B + BSA
 Incubate at 37°C.
To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.
Note: AluB I is able to cleave some methylated DNA substrates.
AluB I is an isoschizomer of Alu I
 See: http://science.sibenzyme.com/article8_article_30_1.phtml

Storage buffer: 10 mM Tris-HCl (pH 7.5); 200 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol. Store at -20°C
Diluent: SE Buffer A
Heat inactivation:
 Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with AluB I, 80% of λ DNA fragments can be ligated with T4 DNA Ligase and can be recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of AluB I for 16 hours.
Do not use BSA for long incubation

| | | | | | |
|-----------------------------------|-----|--------|-------|-------|--------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 100 | 75-100 | 10-25 | 10-25 | 75-100 |

Ama87 I (prototype Ava I) C[^]YCGRG
GRGCY[^]C E017
E018 1000 u.a.
5000 u.a.

Isolated from *Alteromonas macleodii* 87

Concentration: 10 000 - 20 000 units/ml
Assayed on λ DNA
Reagents Supplied with Enzyme:
 10×SE Buffer W, BSA
Reaction conditions: 1×SE Buffer W+BSA
 Incubate at 37°C.
To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM KH₂PO₄ (pH 7.2); 100 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 50% glycerol. Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
 Yes (65° C for 20 minutes)

Ligation/recutting assay: After 20-fold overdigestion with Ama87 I, >90% of the DNA fragments can be ligated and recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of Ama87 I for 16 hours.
Do not use BSA for long incubation.

| | | | | | |
|-----------------------------------|-------|-------|--------|-----|------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 10-25 | 50-75 | 75-100 | 100 | 0-10 |

Aox I (prototype Aox I)Isolated from *Arthrobacter oxydans* 25K**^PuG(5mC)Py
Py(5mC)GPu^****E569
E570****50 u.a.
250 u.a.**

See page 52 for more information about this enzyme.

Apa I (prototype Apa I)Isolated from *Acetobacter pasteurianus***GGGCC^C
C^CCGGG****E019
E020****1000 u.a.
5000 u.a.****Concentration:** 10 000 - 30 000 units/ml**Assayed on** λ DNA (*dcm*⁻, BamHI-digest)**Reagents Supplied with Enzyme:**

10×SE Buffer Y, BSA

Reaction conditions: 1× SE Buffer Y+BSA

Incubate at 37°C.

Blocked by overlapping *dcm*-methylation(C^mCWGG): **GGGCCCWGG**.*To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.***Storage buffer:** 10 mM Tris-HCl (pH 7.5);

200 mM KCl; 0,1 mM EDTA;

7 mM 2-mercaptoethanol; 200 µg/ml BSA;

50% glycerol. Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:****Yes** (65° C for 20 minutes)**Ligation/recutting assay:** After 20-fold overdigestion with Apa I, >90% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 40 units of Apa I for 16 hours.**Do not use** BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|------|------|-----|
| Activity in SE Buffers (% of max) | 50-75 | 25-50 | 0-10 | 0-10 | 100 |

Ars I (prototype Ars I)Isolated from *Arthrobacter* species NTS**^(N)₈GAC(N)₆TTYG(N)₁₁^
^(N)₁₃CTG(N)₆AARC(N)₆^****E575
E576****50 u.a.
250 u.a.****Concentration:** 500 units/ml**Assayed on** T7 DNA**Reagents Supplied with Enzyme:**

10×SE Buffer Y, BSA

Reaction conditions: 1×SE Buffer Y + BSA

Incubate at 30°C.

*To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.***Storage buffer:** 10 mM KH₂PO₄ (pH 7.4);

200 mM KCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 µg/ml BSA;

50% glycerol. Store at -20°C

Diluent: SE Buffer A**Heat inactivation:****Yes** (65° C for 20 minutes)**Ligation/recutting assay:** After 3-fold overdigestion with Ars I, about 70% of the DNA fragments can be ligated. Of these, 80% can be recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of DNA with 1 units of Ars I for 16 hours.**Do not use** BSA for long incubation

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|---|---|---|---|-----|
| Activity in SE Buffers (% of max) | 0 | 0 | 0 | 0 | 100 |

AsiG I (prototype Age I)Isolated from *Arthrobacter* species G**A^CCGGT
TGGCC^A****E235
E236****100 u.a.
500 u.a.****Concentration:** 3 000 –5000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer O

Reaction condition: 1× SE Buffer O

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

250 mM NaCl; 0.1 mM EDTA;

100 µg/ml BSA; 7 mM 2-mercaptoethanol;

50% glycerol. Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:****Yes** (65° C for 20 minutes)**Ligation/recutting assay:** After 5-fold overdigestion with AsiG I, >90% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of AsiG I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|-----|--------|-------|
| Activity in SE Buffers (% of max) | 10-25 | 25-50 | 100 | 75-100 | 10-25 |

AspA2 I (prototype Avr II)Isolated from *Arthrobacter* species A2**C^CTAGG
GGATC^C****E245
E246****500 u.a.
2500 u.a.****Concentration:** 10 000 - 20 000 units/ml**Assayed on** λ DNA (Hind III–digest)**Reagents Supplied with Enzyme:**

10×SE Buffer W, BSA

Reaction conditions: 1× SE Buffer W+BSA

Incubate at 37°C.

*To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.***Storage buffer:** 10 mM Tris-HCl (pH 7.5);

100 mM NaCl; 0.1 mM EDTA;

100 µg/ml BSA; 7 mM 2-mercaptoethanol;

50% glycerol. Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:****No** (65° C for 20 minutes)**Yes** (80° C for 20 minutes)**Ligation/recutting assay:** After 20-fold overdigestion with AspA2 I, > 90% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of the DNA with 20 units of AspA2 I for 16 hours.**Do not use** BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|--------|-----|--------|
| Activity in SE Buffers (% of max) | 10-25 | 50-75 | 75-100 | 100 | 75-100 |

Qualified for human genome studies: http://science.sibenzyme.com/article8_article_31_1.phtml

New product



Mammalian Genome Qualified



New package

AspLE I (prototype Hha I)

Isolated from *Arthrobacter* species LE3860

GCG[^]C
C[^]GCG

E221
E222

500 u.a.
2500 u.a.

Concentration: 10 000 – 30 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer O

Reaction conditions: 1 \times SE Buffer O

Incubate at 37°C.

Blocked by CG methylation

5'-G(5mC)GC-3'/3'-CG(5mC)G-5' or

5'-G(5mC)GC-3'/3'-CGCG-5'

Not blocked by methylation

5'-GCG(5mC)-3'/3'-(5mC)GCG-5' or

5'-GCG(5mC)-3'/3'-CGCG-5'.

Cut hemimethylated site:

5'- G(5mC)GC-3' / 3'-CGCG-5'

Storage buffer: 10 mM Tris-HCl (pH 7.5);
100 mM NaCl; 0.1 mM EDTA;
100 μ g/ml BSA; 7 mM 2-mercaptoethanol;
50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

No (65° C for 20 minutes)

Ligation/recutting assay: After 20-fold
overdigestion with AspLE I,
>90% of the DNA fragments can be
ligated and recut.

Overdigestion assay: No nonspecific
activity was detected after incubation of
1 μ g of λ DNA with 20 units of AspLE
I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|--------|-----|-------|-------|
| Activity in SE Buffers (% of max) | 10-25 | 75-100 | 100 | 50-75 | 25-50 |

AspS9 I (prototype Sau96 I)

Isolated from *Arthrobacter* species S9



G[^]GNCC
CCNG[^]G

E117
E118

1000 u.a.
5000 u.a.

Concentration: 10 000 - 30 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer W

Reaction condition: 1 \times SE Buffer W

Incubate at 37°C.

Blocked by overlapping *dcm*-methylation

(C^mCWGG): GGNCCWGG

Storage buffer: 10 mM Tris-HCl (pH 7.5); 50
mM KCl; 0.1 mM EDTA;
200 μ g/ml BSA; 7 mM 2-mercaptoethanol;
50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 20-fold
overdigestion with AspS9 I,
>90% of the DNA fragments can be
ligated and recut.

Overdigestion assay: No nonspecific
activity was detected after incubation of
1 μ g of λ DNA with 30 units of AspS9 I
for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|--------|-----|-------|
| Activity in SE Buffers (% of max) | 50-75 | 50-75 | 75-100 | 100 | 50-75 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml

AsuC2 I (prototype Cau II)

Isolated from *Actinobacillus suis* CA

CC[^]SGG
GGG[^]CC

E257
E258

2000 u.a.
10000 u.a.

Concentration: 20 000 - 50 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer Y

Reaction conditions: 1 \times SE Buffer Y

Incubate at 37°C.

Blocked by overlapping *dcm*-methylation

(C^mCWGG): GGNCCWGG

Storage buffer: 10 mM Tris-HCl (pH 7.5);
250 mM NaCl; 0.1 mM EDTA;
7 mM 2-mercaptoethanol; 50% glycerol.
Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold
overdigestion with AsuC2 I, ~20% of the
DNA fragments can be ligated. In the
presence of 10% PEG ligation is better.

Overdigestion assay: No nonspecific
activity was detected after incubation of
1 μ g of λ DNA with 50 units of AsuC2 I
for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 75-50 | 50-75 | 10-25 | 25-50 | 100 |

AsuHP I (prototype Hph I)

Isolated from *Actinobacillus suis* HP



GGTGA(N)₈[^]
CCACT(N)₇[^]

E231
E232

200 u.a.
1000 u.a.

Concentration: 2 000 – 5 000 units/ml

Assayed on λ DNA (dam-)

Reagents Supplied with Enzyme:

10 \times SE Buffer O

Reaction conditions: 1 \times SE Buffer O

Incubate at 37°C.

Blocked by overlapping *dam*-methylation

(G^mATC): GGTGATC.

Enzyme may cleave at N₉/N₈ depending on
the sequence between the recognition and
cleave sites.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
250 mM NaCl; 0.1 mM EDTA;
7 mM 2-mercaptoethanol; 200 μ g/ml BSA;
50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 2-fold
overdigestion with AsuHP I, about 30%
of the DNA fragments can be ligated and
recut.

Overdigestion assay: No nonspecific
activity was detected after incubation of
1 μ g of λ DNA with 10 units of AsuHP I
for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|-----|--------|-------|
| Activity in SE Buffers (% of max) | 10-25 | 50-75 | 100 | 75-100 | 25-50 |

Qualified for human genome studies: http://science.sibenzyme.com/article8_article_31_1.phtml

AsuNH I (prototype Nhe I)

Isolated from *Actinobacillus suis* NH

G[^]CTAGC
CGATC[^]G

E063
E064

1000 u.a.
5000 u.a.

Concentration: 10 000 - 20 000 units/ml

Assayed on λ DNA (HindIII-digest)

Reagents Supplied with Enzyme:

10 \times SE Buffer Y, BSA

Reaction conditions: 1 \times SE Buffer Y+BSA

Incubate at 37°C.

*To obtain 100% activity, BSA should be
added to the 1 \times reaction mix to a final
concentration of 100 μ g/ml.*

Storage buffer: 10 mM Tris-HCl (pH 7.5);
250 mM NaCl; 0.1 mM EDTA;
7 mM 2-mercaptoethanol; 100 μ g/ml BSA;
50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold
overdigestion with AsuNH I,
> 90% of the DNA fragments can be
ligated and recut.

Overdigestion assay: No nonspecific
activity was detected after incubation of
1 μ g of the DNA with 20 units of
AsuNH I for 16 hours.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|-------|------|------|-----|
| Activity in SE Buffers (% of max) | 75-100 | 50-75 | 0-10 | 0-10 | 100 |



BamH I (prototype BamH I)

Isolated from an *E.coli* strain that carries the cloned BamHI gene from *Bacillus amyloliquefaciens* H

G[^]GATCC **E021** **4000 u.a.**
CCTAG[^]G **E022** **20000 u.a.**

For high concentration

E021X **4000 u.a.**
E022X **20000 u.a.**

Concentration: 20 000 and 50 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer G, BSA

Reaction conditions: 1 \times SE Buffer G+BSA
Incubate at 37°C.

Not blocked by overlapping

dam-methylation (G^mATC): **GGATCC**.

To obtain 100% activity, BSA should be added to the 1 \times reaction mix to a final concentration of 100 μ g/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.6); 50 mM NaCl; 0,1 mM EDTA; 100 μ g/ml BSA, 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 50-fold overdigestion with BamH I, > 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 40 units of BamH I for 16 hours.

Star activity: High enzyme concentration may result in star activity.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-----|--------|--------|-------|
| Activity in SE Buffers (% of max) | 25-50 | 100 | 75-100 | 75-100 | 25-50 |

Bar I (prototype Bar I)

Isolated from *Bacillus sphaericus*

^(N)₇GAAG(N)₆TAC(N)₁₂^

^(N)₁₂CTTC(N)₆ATG(N)₇^

E547 **100 u.a.**

E548 **500 u.a.**

Concentration: 500 - 2 000 units/ml

Assayed on T7 DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer 2K

Reaction conditions: 1 \times SE Buffer 2K

Incubate at 37°C.

Storage buffer: 20 mM KH₂PO₄ (pH 7.4); 100 mM KCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 μ g/ml BSA, 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 2-fold overdigestion with Bar I, about 90% of T7 DNA fragments can be ligated and 95% of these can be recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of T7 DNA with 2 units of Bar I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|---|------|-------|-------|-------|
| Activity in SE Buffers (% of max) | 0 | 0-10 | 25-50 | 50-75 | 10-25 |

Bbv12 I (prototype HgiA I)

Isolated from *Bacillus brevis* 12

GWGCW[^]C

C[^]WCGWG

E023 **200 u.a.**

E024 **1000 u.a.**

Concentration: 1 000 - 5 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer O

Reaction conditions: 1 \times SE Buffer O

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 μ g/ml BSA, 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (80°C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with Bbv12 I, >90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 5 units of Bbv12 I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|------|-------|-----|--------|-------|
| Activity in SE Buffers (% of max) | 0-10 | 10-25 | 100 | 75-100 | 10-25 |

Bgl I (prototype Bgl I)

Isolated from *Bacillus globigii*

GCCNNNN[^]NGGC

CGGN[^]NNNNCCG

E025 **500 u.a.**

E026 **2500 u.a.**

Concentration: 5 000 - 10 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer 2W

Reaction conditions: 1 \times SE Buffer 2W

Incubate at 37°C.

Not blocked by *dcm*-methylation

(C^mCWGG): **GCCWGGNNGGC**.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 200 mM NaCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 μ g/ml BSA; 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Bgl I, > 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 10 units of Bgl I for 16 hours.

Star activity: High enzyme concentration may result in star activity.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|------|--------|-------|
| Activity in SE Buffers (% of max) | 50-75 | 50-75 | 0-10 | 75-100 | 25-50 |

Bgl II (prototype Bgl II)

Isolated from an *E.coli* strain that carries the cloned Bgl II gene from *Bacillus globigii*

A[^]GATCT

TCTAG[^]A

E027 **1000 u.a.**

E028 **5000 u.a.**

Concentration: 10 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer O

Reaction conditions: 1 \times SE Buffer O

Incubate at 37°C.

Not blocked by overlapping

dam-methylation (G^mATC): **AGATCT**.

Storage buffer: 10mM Tris-HCl(pH 7.5); 50 mM NaCl; 0,1 mM EDTA; 1 mM DTT; 200 μ g/ml BSA; 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

No (65°C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Bgl II, > 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 20 units of Bgl II for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|------|-------|-----|-------|-------|
| Activity in SE Buffers (% of max) | 0-10 | 10-25 | 100 | 25-50 | 10-25 |



New product



Mammalian Genome Qualified



New package

Bis I (prototype Bis I)Isolated from *Bacillus subtilis* T30G(5mC)[^]NGC
CGN[^](5mC)GE485
E48640 u.a.
200 u.a.

See page 52 for more information about this enzyme.

Bls I (prototype Bls I)Isolated from *Bacillus simplex* 23

DNA sequence with at least two 5mC:

PuPyN[^]PuPy
PyPu[^]NPyPuE533
E534100 u.a.
500 u.a.

See page 53 for more information about this enzyme.

Bme18 I (prototype Ava II)Isolated from *Bacillus megaterium* 18G[^]GWCC
CCWG[^]GE029
E0301000 u.a.
5000 u.a.**Concentration:** 5 000 – 20 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer O

Reaction conditions: 1× SE Buffer O

Incubate at 37°C.

Cleaved of DNA is impaired byoverlapping dcm-methylation (C^mCWGG):**GGWCCWGG.****Storage buffer:** 10mM Tris-HCl(pH 7.5);
100 mM NaCl; 0.1 mM EDTA;7 mM 2-mercaptoethanol; 100 µg/ml BSA;
50% glycerol. Store at -20°C.**Diluent:** SE Buffer A**Heat inactivation:**

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold
overdigestion with Bme18 I, > 90% of
the DNA fragments can be ligated and
recut.**Overdigestion assay:** No nonspecific
activity was detected after incubation of
1 µg of λ DNA with 20 units of Bme18
I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|-----|--------|-------|
| Activity in SE Buffers (% of max) | 10-25 | 25-50 | 100 | 75-100 | 10-25 |

Bmt I (prototype Nhe I)Isolated from an *E.coli* strain that carries the
cloned Bmt I gene from *Bacillus megaterium* S2GCTAG[^]C
C[^]GATCGE457
E4581000 u.a.
5000 u.a.**Concentration:** 20 000 units/ml**Assayed on** λ DNA (HindIII-digest)**Reagents Supplied with Enzyme:**

10×SE Buffer W

Reaction conditions: 1×SE Buffer W

Incubate at 37°C.

BmtI is a neoschizomer of NheI.**Storage buffer:** 10mM Tris-HCl(pH 7.5);
250 mM NaCl; 0,1 mM EDTA;200 µg/ml BSA; 7 mM 2-mercaptoethanol;
50% glycerol. Store at -20°C.**Diluent:** SE Buffer A**Heat inactivation:**

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 20-fold
overdigestion with Bmt I about 95% of
the DNA fragments can be ligated and
recut.**Overdigestion assay:** No nonspecific
activity was detected after incubation of
1 µg of λ DNA with 20 units of Bmt I
for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|-------|-----|--------|
| Activity in SE Buffers (% of max) | 10-25 | 50-75 | 50-75 | 100 | 75-100 |

Bmu I (prototype Bfi I)Isolated from *Bacillus megaterium* S2ACTGGG(N)₅[^]
TGACCC(N)₄[^]E487
E48850 u.a.
250 u.a.**Concentration:** 500 - 1000 units/ml**Assayed on** λ DNA (HindIII-digest)**Reagents Supplied with Enzyme:**

10×SE Buffer Y

Reaction conditions: 1×SE Buffer Y

Incubate at 37°C.

Enzyme is active in presence of EDTA.

Storage buffer: 10mM Tris-HCl(pH 7.5);
250 mM NaCl; 0,1 mM EDTA;7 mM 2-mercaptoethanol; 50% glycerol.
Store at -20°C.**Diluent:** SE Buffer A**Heat inactivation:**

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 2-fold
overdigestion with Bmu I about 75% of
the DNA fragments can be ligated and
95% of these can be recut.**Overnight digest** with BmuI is not
recommended. A 50 µl reaction
containing 1 µg of λ DNA and 0.5 units
of enzyme incubated for 4 hours resulted
in the same pattern of DNA bands as a
reaction incubated for 1 hour.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|--------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 25-50 | 10-25 | 100 |

Bpm I (prototype Gsu I)

Isolated from *Bacillus pumilus*

Concentration: 200 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer W, BSA

Reaction conditions: 1 \times SE Buffer W+BSA

Incubate at 37°C.

To obtain 100% activity, BSA should be

added to the 1 \times reaction mix to a final

concentration of 100 μ g/ml.

CTGGAG(N)₁₆[^]

GACCTC(N)₁₄[^]

E467

E468

50 u.a.

250 u.a.

Storage buffer: 10mM Tris-HCl(pH 7.5); 50

mM KCl; 0,1 mM EDTA;

10 mM 2-mercaptoethanol;

200 μ g/ml BSA; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 2-fold overdigestion with Bpm I about 95% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 0,4 units of Bpm I for 16 hours.

Do not use BSA for long incubation.

| | | | | | |
|-----------------------------------|-------|-------|--------|-----|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 25-50 | 50-75 | 75-100 | 100 | 50-75 |

Bpu10 I (prototype Bpu10 I)

Isolated from an *Escherichia coli* strain that carries plasmids pBpu10IA and pBpu10IB



CC[^]TNAGC

GGANT[^]CG

E149

E150

200 u.a.

1000 u.a.

Concentration: 5 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer K

Reaction conditions: 1 \times SE Buffer K

Incubate at 37°C.

Storage buffer: 10mM Tris-HCl(pH 7.5);

50 mM NaCl; 0.1 mM EDTA; 1 mM DTT;

50% glycerol.

Store at -20°C.

Diluent: SE Buffer A

Heat Inactivation:

No (65° C for 20 minutes)

Yes (80° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with Bpu10 I, 80% of the DNA fragments can be ligated. Of these, 90% can be recut. In the presence of 10% PEG ligation is better.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of DNA with 5 units of enzyme for 16 hour.

Note: High enzyme concentration may result in star activity or incomplete DNA cleavage. We recommend increasing the incubation time instead of using an excess of Bpu10 I.

| | | | | | |
|-----------------------------------|-------|-------|-------|-------|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 10-25 | 25-50 | 50-75 | 50-75 | 25-50 |

Qualified for human genome studies: http://science.sibenzyme.com/article8_article_31_1.phtml

Bpu14 I (prototype Asu II)

Isolated from *Bacillus pumilus* 14

TT[^]CGAA

AAGC[^]TT

E033

E034

1000 u.a.

5000 u.a.

Concentration: 10 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer G

Reaction conditions: 1 \times SE Buffer G

Incubate at 37°C.

Storage buffer: 10mM Tris-HCl(pH 7.5);

100 mM NaCl; 0,1 mM EDTA;

7 mM 2-mercaptoethanol; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Bpu14 I, 95% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 20 units of Bpu14 I for 16 hours.

| | | | | | |
|-----------------------------------|-------|-----|-------|-------|--------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 50-75 | 100 | 25-50 | 25-50 | 75-100 |

Bsa29 I (prototype Cla I)

Isolated from *Bacillus stearothermophilus* 29

AT[^]CGAT

TAGC[^]TA

E205

E206

1000 u.a.

5000 u.a.

Concentration: 20 000 units/ml

Assayed on λ DNA (dam⁻)

Reagents Supplied with Enzyme:

10 \times SE Buffer G, BSA

Reaction conditions: 1 \times SE Buffer G+BSA

Incubate at 37°C.

Blocked by overlapping *dam*-methylation

(G^mATC): GATCGATC.

Blocked by CG methylation.

To obtain 100% activity, BSA should be

added to the 1 \times reaction mix to a final

concentration of 100 μ g/ml.

Storage buffer: 10mM Tris-HCl(pH 7.5);

250 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol;

200 μ g/ml BSA, 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 20-fold overdigestion with Bsa29 I,

> 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 40 units of Bsa29 I for 16 hour.

Do not use BSA for long incubation.

| | | | | | |
|-----------------------------------|-------|-----|-------|-------|--------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 25-50 | 100 | 50-75 | 50-75 | 75-100 |



New product



Mammalian Genome Qualified



New package

Bsc4 I (prototype BsiY I)Isolated from *Bacillus schlegelii* 4CCNNNN[^]NNGG
GGNN[^]NNNNCCE219
E220500 u.a.
2500 u.a.**Concentration:** 10 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer W, BSA

Reaction conditions: 1×SE Buffer W+BSA

Incubate at 55°C.

*To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.***Storage buffer:** 10mM Tris-HCl(pH 7.5); 50 mM KCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol. Store at -20°C.**Diluent:** SE Buffer A**Heat inactivation:**

Yes (80° C for 20 minutes)

Ligation/recutting assay: After 20-fold overdigestion with Bsc4 I, > 90% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 30 units of Bsc4 I for 16 hours.**Do not use** BSA for long incubation.

| | | | | | |
|-----------------------------------|--------|--------|-------|-----|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 50-75 | 100 | 25-50 |

Bse1 I (prototype Bsr I)Isolated from *Bacillus stearothermophilus* 1ACTGGN[^]
TGAC[^]CNE035
E0361000 u.a.
5000 u.a.**Concentration:** 10 000 – 20 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer Y

Reaction conditions: 1× SE Buffer Y

Incubate at 65°C.

Storage buffer: 10mM Tris-HCl(pH 7.5); 100 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 100 µg/ml BSA; 50% glycerol. Store at -20°C.**Diluent:** SE Buffer A**Heat inactivation:**

No (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Bse1 I, 95% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 40 units of Bse1 I for 16 hours.

| | | | | | |
|-----------------------------------|--------|--------|-------|-------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 25-50 | 10-25 | 100 |

Bse118 I (prototype Cfr10 I)Isolated from *Bacillus stearothermophilus* 118R[^]CCGGY
YGGCC[^]RE039
E040200 u.a.
1000 u.a.**Concentration:** 2 000-5 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer O

Reaction conditions: 1× SE Buffer O

Incubate at 65°C.

Storage buffer: 10 mM KH₂PO₄ (pH 7.5); 100 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 100 µg/ml BSA; 50% glycerol. Store at -20°C.**Diluent:** SE Buffer A**Heat inactivation:**

No (65° C for 20 minutes)

Ligation/recutting assay: After 2-fold overdigestion with Bse118 I, 95% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 2 units of Bse118 I for 16 hours.

| | | | | | |
|-----------------------------------|------|-------|-----|--------|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 0-10 | 50-75 | 100 | 75-100 | 25-50 |

Bse21 I (prototype Sau I)Isolated from *Bacillus species* 21CC[^]TNAGG
GGANT[^]CCE037
E038500 u.a.
2500 u.a.**Concentration:** 10 000 – 30 000 units/ml**Assayed on** λ DNA (Hind III-digest)**Reagents Supplied with Enzyme:**

10×SE Buffer Y

Reaction conditions: 1×SE Buffer Y

Incubate at 37°C.

Storage buffer: 10 mM KH₂PO₄ (pH 7.4); 50 mM KCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA, 50% glycerol. Store at -20°C.**Diluent:** SE Buffer A**Heat inactivation:**

No (65° C for 20 minutes)

Ligation/recutting assay: After 2-fold overdigestion with Bse21 I, about 50% of the DNA fragments can be ligated (by using of high concentration T4 DNA Ligase and 10% PEG). Of these, >90% can be recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 30 units of Bse21 I for 16 hours.

| | | | | | |
|-----------------------------------|-------|-------|-------|-------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 50-75 | 50-75 | 10-25 | 25-50 | 100 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml**Bse3D I (prototype BsrD I)**Isolated from *Bacillus stearothermophilus* 3DGCAATGNN[^]
CGTTAC[^]NNE253
E254200 u.a.
1000 u.a.**Concentration:** 5 000 - 10 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer G

Reaction conditions: 1× SE Buffer G

Incubate at 60°C.

Storage buffer: 10mM Tris-HCl(pH 7.5); 200 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 50% glycerol. Store at -20°C.**Diluent:** SE Buffer A**Heat inactivation:**

No (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Bse3D I, > 90% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of Bse3D I for 16 hours.

| | | | | | |
|-----------------------------------|-------|-----|-------|-------|--------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 10-25 | 100 | 25-50 | 50-75 | 75-100 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml

Bse8 I (prototype BsaB I)Isolated from *Bacillus* species 8GATNN[^]NNATC
CTANN[^]NNTAGE147
E1481000 u.a.
5000 u.a.

Concentration: 5 000 units/ml
Assayed on λ DNA
Reagents Supplied with Enzyme:
 10×SE Buffer G
Reaction conditions: 1×SE Buffer G
 Incubate at 60°C.

Storage buffer: 10mM Tris-HCl(pH 7.5);
 100 mM NaCl; 0,1 mM EDTA;
 7 mM 2-mercaptoethanol;
 100 µg/ml BSA; and 50% glycerol.
 Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
 No (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with Bse8 I, 80% of the DNA fragments can be ligated and recut.
Overdigestion assay: Long incubation is not recommended owing to occurrence of star activity.
Star activity is observed at a greater than 5-fold overdigestion of 1 µg substrate with Bse8 I.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-----|--------|--------|-------|
| Activity in SE Buffers (% of max) | 25-50 | 100 | 75-100 | 75-100 | 50-75 |

BseP I (prototype BseP I)Isolated from *Bacillus stearothermophilus* PG[^]CGCGC
CGCGC[^]GE181
E182200 u.a.
1000 u.a.

Concentration: 5 000 units/ml
Assayed on λ DNA
Reagents Supplied with Enzyme:
 10×SE Buffer G
Reaction conditions: 1× SE Buffer G
 Incubate at 50°C.
Blocked by CG methylation.

Storage buffer: 10mM Tris-HCl(pH 7.5); 50 mM KCl; 0,1 mM EDTA;
 10 mM 2-mercaptoethanol;
 200 µg/ml BSA; 50% glycerol.
 Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
 Yes (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with BseP I, 90% of the DNA fragments can be ligated and recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of BseP I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-----|--------|-------|-------|
| Activity in SE Buffers (% of max) | 50-75 | 100 | 75-100 | 50-75 | 50-75 |

BseX3 I (prototype Xma III)Isolated from *Bacillus stearothermophilus* X3C[^]GGCCG
GCCGG[^]CE263
E264200 u.a.
1000 u.a.

Concentration: 5 000 units/ml
Assayed on λ DNA
Reagents Supplied with Enzyme:
 10×SE Buffer O
Reaction conditions: 1× SE Buffer O
 Incubate at 50°C.
Blocked by CG methylation.

Storage buffer: 10mM Tris-HCl(pH 7.5);
 100 mM NaCl; 0,1 mM EDTA;
 7 mM 2-mercaptoethanol;
 200 µg/ml BSA; 50% glycerol.
 Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
 No (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with BseX3 I, > 90% of the DNA fragments can be ligated. Of these, 80% can be recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of BseX3 I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|-----|-------|-------|
| Activity in SE Buffers (% of max) | 10-25 | 25-50 | 100 | 50-75 | 10-25 |

BslF I (prototype Fin I)Isolated from *Bacillus stearothermophilus* FIGGGAC(N)₁₀[^]
CCCTG(N)₁₄[^]E479
E480100 u.a.
500 u.a.

Concentration: 1 000 units/ml
Assayed on λ DNA
Reagents Supplied with Enzyme:
 10×SE Buffer Y, BSA
Reaction conditions: 1×SE Buffer Y+BSA
 Incubate at 37°C.
 BslF I also cleaves the sequence
 GGGAC(11/15).
 To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10mM Tris-HCl(pH 7.5);
 250 mM NaCl; 0,1 mM EDTA;
 7 mM 2-mercaptoethanol; 50% glycerol.
 Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
 No (65° C for 20 minutes)
 Yes (80° C for 20 minutes)

Ligation/recutting assay: After 3-fold overdigestion with BslFI I, 90% of the DNA fragments can be ligated. Of these, 95% can be recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 2 units of BslF I for 16 hours.
Star activity: High enzyme concentration may result in star activity.
Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 25-50 | 25-50 | 10-25 | 25-50 | 100 |

Bso31 I (prototype Eco31 I)Isolated from *Bacillus stearothermophilus* 31GGTCTC(N)₁[^]
CCAGAG(N)₅[^]E285
E286200 u.a.
1000 u.a.

Concentration: 5 000 – 10 000 units/ml
Assayed on T7 DNA
Reagents Supplied with Enzyme:
 10×SE Buffer O, BSA
Reaction conditions: 1×SE Buffer O+BSA
 Incubate at 55°C.
Not blocked by methylation GGTCT^mC.
 To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10mM Tris-HCl(pH 7.5);
 100 mM NaCl; 0,1 mM EDTA;
 100 µg/ml BSA;
 7 mM 2-mercaptoethanol; 50% glycerol.
 Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
 No (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with Bso31 I, 90% of the DNA fragments can be ligated. Of these, 80% can be recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of T7 DNA with 5 units of Bso31 I for 16 hours.
Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|--------|-----|--------|-------|
| Activity in SE Buffers (% of max) | 25-50 | 75-100 | 100 | 75-100 | 25-50 |



New product



Mammalian Genome Qualified



New package

Bsp13 I (prototype BspM II)

Isolated from *Bacillus* species 13

T[^]CCGGA
AGGCC[^]T

E183
E184

1000 u.a.
5000 u.a.

Concentration: 10 000 – 20 000 units/ml
Assayed on λ DNA (dam⁻)
Reagents Supplied with Enzyme:
10×SE Buffer 2K
Reaction conditions: 1× SE Buffer 2K
Incubate at 50°C.
Blocked by overlapping *dam*-methylation
(G^mATC):
TCCGGATC and GATCCGGA.

Storage buffer: 10mM Tris-HCl(pH 7.5);
200 mM KCl; 0,1 mM EDTA;
7 mM 2-mercaptoethanol;
200 µg/ml BSA; 50% glycerol.
Store at -20°C
Diluent: SE Buffer A
Heat inactivation:
Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Bsp13 I, > 90% of the DNA fragments can be ligated and recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of Bsp13 I for 16 hours.

| | | | | | |
|-----------------------------------|-------|-------|--------|-------|------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 25-50 | 50-75 | 75-100 | 50-75 | 0-10 |

Bsp1720 I (prototype Esp I)

Isolated from *Bacillus* species 1720

GC[^]TNAGC
CGANT[^]CG

E185
E186

500 u.a.
2500 u.a.

Concentration: 10 000 units/ml
Assayed on λ DNA
Reagents Supplied with Enzyme:
10×SE Buffer G
Reaction mixture: 1× SE Buffer G
Incubate at 37°C.

Storage buffer: 10mM Tris-HCl(pH 7.5);
250 mM NaCl; 0,1 mM EDTA;
7 mM 2-mercaptoethanol;
50% glycerol. Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
No (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Bsp1720I, about 80% of the DNA fragments can be ligated. Of these, 95 % can be recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of Bsp1720 I for 16 hours.

| | | | | | |
|-----------------------------------|-------|-----|-------|-------|--------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 50-75 | 100 | 50-75 | 50-75 | 75-100 |

Bsp19 I (prototype Nco I)

Isolated from *Bacillus* species 19

C[^]CATGG
GGTAC[^]C

E047
E048

1000 u.a.
5000 u.a.

Concentration: 10 000 – 20 000 units/ml
Assayed on λ DNA
Reagents Supplied with Enzyme:
10×SE Buffer 2W, BSA
Reaction conditions:
1×SE Buffer 2W+BSA
Incubate at 37°C.
Bsp19I cuts hemimethylated site
5⁻-(5mC)CATGG-3[^]/3[^]-GGTACC-5[^]
and **doesn't cut methylated sites**
5⁻-(5mC)CATGG-3[^]/3[^]-GGTAC(5mC)-5[^] and
5⁻-(4mC)CATGG-3[^]/3[^]-GGTAC(4mC)-5[^].
To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10mM Tris-HCl(pH 7.5); 50 mM KCl; 0,1 mM EDTA;
7 mM 2-mercaptoethanol;
200 µg/ml BSA; 50% glycerol.
Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Bsp19 I, > 90% of the DNA fragments can be ligated and recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of Bsp19 I for 16 hours.
Do not use BSA for long incubation.

| | | | | | |
|-----------------------------------|------|-------|-------|--------|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 0-10 | 10-25 | 50-75 | 75-100 | 10-25 |

BspAC I (prototype Aci I)

Isolated from *Bacillus* species AC



C[^]CGC
GGC[^]G

E501
E502

200 u.a.
1000 u.a.

Concentration: 2 000 – 5 000 units/ml
Assayed on λ DNA
Reagents Supplied with Enzyme:
10×SE Buffer O, BSA
Reaction conditions: 1×SE Buffer O+BSA
Incubate at 37°C.
Blocked by CG methylation.
Note: BspACI has a non-palindromic recognition site.
To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM KH₂PO₄ (pH 7.2);
100 mM NaCl; 0,1 mM EDTA;
7 mM 2-mercaptoethanol;
200 µg/ml BSA; 50% glycerol.
Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
Yes (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with BspAC I, > 95% of λ DNA fragments can be ligated with T4 DNA Ligase at 16°C and 50% of these can be recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of BspAC I at 37°C for 16 hours.
Do not use BSA for long incubation.

| | | | | | |
|-----------------------------------|-------|-------|-----|--------|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 10-25 | 25-50 | 100 | 75-100 | 10-25 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml

BspFN I (prototype FnuD II)

Isolated from *Bacillus* species FN

CG[^]CG
GC[^]GC

E557
E558

500 u.a.
2500 u.a.

Concentration: 5 000 units/ml
Assayed on λ DNA
Reagents Supplied with Enzyme:
10×SE Buffer Y
Reaction conditions: 1×SE Buffer Y
Incubate at 37°C.
Blocked by CG methylation.

Storage buffer: 20 mM Tris-HCl (pH 7.5);
300 mM NaCl; 0,1 mM EDTA;
7 mM 2-mercaptoethanol; 200 µg/ml BSA;
50% glycerol.
Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
Yes (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with BspFN I, > 95% of λ DNA fragments can be ligated with T4 DNA Ligase and recut
Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of BspFN I at 37°C for 16 hours.

| | | | | | |
|-----------------------------------|-------|--------|--------|-------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 50-75 | 75-100 | 75-100 | 50-75 | 100 |

BssEC I (prototype Sec I)Isolated from *Bacillus stearothermophilus* ECC[^]CNNGG
GGNC[^]CE273
E274200 u.a.
1000 u.a.**Concentration:** 10 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer Y

Reaction conditions: 1×SE Buffer Y

Incubate at 60°C.

Storage buffer: 10mM Tris-HCl(pH 7.5); 50 mM KCl; 0,1 mM EDTA;

200 µg/ml BSA;

7 mM 2-mercaptoethanol; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:**

No (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with BssEC I 90% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of BssEC I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|-------|--------|-----|
| Activity in SE Buffers (% of max) | 50-75 | 50-75 | 50-75 | 75-100 | 100 |

Qualified for human genome studies: http://science.sibenzyme.com/article8_article_31_1.phtml**BssNA I (prototype Sna I)**Isolated from *Bacillus stearothermophilus* NAGTA[^]TAC
CAT[^]ATGE261
E2621000 u.a.
5000 u.a.**Concentration:** 10 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer W, BSA

Reaction conditions:

1× SE Buffer W+BSA

Incubate at 37°C.

*To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.***Storage buffer:** 10mM Tris-HCl(pH 7.5); 100 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol;

100 µg/ml BSA; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:**

No (65° C, 80° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with BssNA I, > 90% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of BssNA I for 16 hours**Star activity:** High enzyme concentration results in star activity.**Do not use BSA for long incubation.**

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|--------|-----|--------|
| Activity in SE Buffers (% of max) | 50-75 | 50-75 | 75-100 | 100 | 75-100 |

BssT1 I (prototype Sty I)Isolated from *Bacillus stearothermophilus* T1C[^]CWWGG
GGW[^]C[^]CE207
E2081000 u.a.
5000 u.a.**Concentration:** 20 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer 2K

Reaction conditions: 1× SE Buffer 2K

Incubate at 60°C.

Storage buffer: 10mM Tris-HCl(pH 7.5); 50 mM KCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol;

200 µg/ml BSA; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:**

Yes (80° C for 20 minutes)

Ligation/recutting assay: After 20-fold overdigestion with BssT1 I, > 90% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 40 units of BssT1 I for 16 hours.**Star activity:** High enzyme concentration results in star activity.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|-------|--------|-------|
| Activity in SE Buffers (% of max) | 10-25 | 25-50 | 25-50 | 75-100 | 10-25 |

Bst2B I (prototype Bsi I)Isolated from *Bacillus stearothermophilus* 2BC[^]TCGTG
GAGCA[^]CE043
E044200 u.a.
1000 u.a.**Concentration:** 5 000 - 10 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer Y, BSA

Reaction conditions: 1×SE Buffer Y+BSA

Incubate at 60°C.

*To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.***Storage buffer:** 10mM Tris-HCl(pH 7.5); 200 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:**

No (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Bst2B I, 95% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of Bst2B I for 16 hours.**Do not use BSA for long incubation.**

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|-------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 75-100 | 25-50 | 10-25 | 25-50 | 100 |

Bst2U I (prototype BstN I)Isolated from *Bacillus stearothermophilus* 2UCC[^]WGG
GGW[^]CCE051
E0521000 u.a.
5000 u.a.**Concentration:** 10 000 – 20 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer G, BSA

Reaction conditions: 1× SE Buffer G+BSA

Incubate at 60°C.

Not blocked by overlapping*dcm*-methylation (C^mCWGG): **CCWGG**.*To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.***Storage buffer:** 10mM Tris-HCl(pH 7.5); 200 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol;

200 µg/ml BSA, 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:**

No (80° C for 20 minutes)

Ligation/recutting assay: After 2-fold overdigestion with Bst2U I, none of the DNA fragments can be ligated.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of Bst2U I for 16 hours.**Do not use BSA for long incubation.**

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|-----|-------|-------|-------|
| Activity in SE Buffers (% of max) | 75-100 | 100 | 50-75 | 50-75 | 10-25 |

Qualified for human genome studies: http://science.sibenzyme.com/article8_article_31_1.phtml

New product



Mammalian Genome Qualified



New package

Bst4C I (prototype Tsp4C I)Isolated from *Bacillus stearothermophilus* 4CACN[^]GT
TG[^]NCAE265
E266500 u.a.
2500 u.a.**Concentration:** 10 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer Y

Reaction conditions: 1×SE Buffer Y

Incubate at 65°C.

Storage buffer 10 mM Tris-HCl (pH 7.5);

200 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol;

200 µg/ml BSA; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:**

No (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Bst4C I 50% of the DNA fragments can be ligated and recut. In the presence of 10% PEG ligation is better.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of Bst4C I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|--------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 10-25 | 25-50 | 100 |

Bst6 I (prototype Ksp632 I)Isolated from *Bacillus stearothermophilus* 6CTCTTC(N)₁[^]
GAGAAG(N)₄[^]E239
E240200 u.a.
1000 u.a.**Concentration:** 1 000 – 5 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer Y, BSA

Reaction conditions:

1× SE Buffer Y+BSA

Incubate at 65°C.

*To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.***Storage buffer:** 10 mM Tris-HCl (pH 7.5);

50 mM KCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol;

200 µg/ml BSA, 50% glycerol.

Store at -20°C (*see note).

Diluent: SE Buffer A**Heat inactivation:**

No (65° C for 20 minutes)

Yes (80° C for 20 minutes)

Ligation/recutting assay: After 2-fold overdigestion with Bst6 I, 80% of the DNA fragments can be ligated and 80% of those can be recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 5 units of Bst6I for 16 hours.**Do not use** BSA for long incubation.***Note: For long term storage (more than 30 days), store at -70°C .**

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|--------|-------|--------|-----|
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 50-75 | 75-100 | 100 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml**BstAC I (prototype Acy I)**Isolated from *Bacillus stearothermophilus* ACGR[^]CGYC
CYGC[^]RGE093
E094500 u.a.
2500 u.a.**Concentration:** 10 000 – 20 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer W

Reaction conditions: 1×SE Buffer W

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

200 mM KCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol;

200 µg/ml BSA; and 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:**

Yes (80° C for 20 minutes)

Ligation/recutting assay: After 40-fold overdigestion with BstAC I, 95% of λ DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 40 units of BstACI for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|--------|-------|-----|--------|
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 50-75 | 100 | 75-100 |

BstAF I (prototype Afl II)Isolated from *Bacillus stearothermophilus* AFC[^]TTAAG
GAATT[^]CE135
E1361000 u.a.
5000 u.a.**Concentration:** 20 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer W , BSA

Reaction conditions:

1×SE Buffer W + BSA

Incubate at 55°C.

*To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.***Storage buffer:** 10 mM Tris-HCl (pH 7.5);

200 mM KCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 µg/ml BSA;

50% glycerol. Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:**

Yes (80° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with BstAF I, about 40% of the DNA fragments can be ligated and 95% of those can be recut. In the presence of 10% PEG ligation is better.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 40 units of BstAF I for 16 hours.**Do not use** BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|--------|-----|-------|
| Activity in SE Buffers (% of max) | 10-25 | 25-50 | 75-100 | 100 | 25-50 |

BstAP I (prototype ApaB I)Isolated from *Bacillus stearothermophilus* APGCANNNN[^]NTGC
CGTN[^]NNNNACGE259
E260200 u.a.
1000 u.a.**Concentration:** 5 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer W

Reaction conditions: 1×SE Buffer W

Incubate at 60°C.

BstAPI is a neoschizomer of ApaBI.**Storage buffer:** 10 mM Tris-HCl (pH 7.5);
50 mM KCl; 0,1 mM EDTA; 1 mM DTT;
200 µg/ml BSA; and 50% glycerol.
Store at -20°C.**Diluent:** SE Buffer A**Heat inactivation:****No** (65° C for 20 minutes)**Yes** (80° C for 20 minutes)**Ligation/recutting assay:** After 5-fold
overdigestion with BstAP I, > 90% of
the DNA fragments can be ligated and
recut.**Overdigestion assay:** No nonspecific
activity was detected after incubation of
1 µg of λ DNA with 5 units of BstAPI
for 16 hours.**Star activity:** High enzyme
concentration (>5 units for 16 hours) on
1µg of DNA may result in star activity.

| | | | | | |
|-----------------------------------|-------|-------|--------|-----|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 25-50 | 25-50 | 75-100 | 100 | 25-50 |

BstAU I (prototype Bsp1407 I)Isolated from *Bacillus stearothermophilus* AUT[^]GTACA
ACATG[^]TE267
E2681000 u.a.
5000 u.a.**Concentration:** 20 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer W

Reaction conditions: 1×SE Buffer W

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
100 mM NaCl; 0,1 mM EDTA;
7 mM 2-mercaptoethanol;
100 µg/ml BSA; 50% glycerol.
Store at -20°C.**Diluent:** SE Buffer A**Heat inactivation:****No** (65° C for 20 minutes)**Yes** (80° C for 20 minutes)**Ligation/recutting assay:** After 10-fold
overdigestion with BstAU I about 90%
of λ DNA fragments can be ligated and
recut.**Overdigestion assay:** No nonspecific
activity was detected after incubation of
1 µg of λ DNA with 20 units of BstAU I
for 16 hours.

| | | | | | |
|-----------------------------------|-------|-------|-------|-----|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 10-25 | 50-75 | 25-50 | 100 | 25-50 |

BstBA I (prototype BsaA I)Isolated from *Bacillus stearothermophilus* BAYAC[^]GTR
RTG[^]CAYE237
E238500 u.a.
2500 u.a.**Concentration:** 5 000 - 20 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer W, BSA

Reaction Conditions:

1×SE Buffer W+BSA

Incubate at 65°C.

Blocked by CG methylation.*To obtain 100% activity, BSA should be
added to the 1× reaction mix to a final
concentration of 100 µg/ml.***Storage buffer:** 10 mM Tris-HCl (pH 7.5);
50 mM KCl; 0.1 mM EDTA;
10 mM 2-mercaptoethanol;
200 µg/ml BSA; 50% glycerol.
Store at -20°C.**Diluent:** SE Buffer A**Heat inactivation:****Yes** (80° C for 20 minutes)**Ligation/recutting assay:** After 10-fold
overdigestion with BstBA I, > 90% of
the DNA fragments can be ligated and
recut.**Overdigestion assay:** No nonspecific
activity was detected after incubation of
1 µg of λ DNA with 20 units of BstBA I
for 16 hours.**Do not use** BSA for long incubation.

| | | | | | |
|-----------------------------------|-------|-------|--------|-----|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 25-50 | 25-50 | 75-100 | 100 | 25-50 |

BstC8 I (prototype Cac8 I)Isolated from *Bacillus stearothermophilus* BAGCN[^]NGC
CGN[^]NCGE305
E306500 u.a.
2500 u.a.**Concentration:** 10 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer Y

Reaction conditions: 1×SE Buffer Y

Incubate at 55°C.

Incubation at 37°C results in 50% activity.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
50 mM KCl; 0.1 mM EDTA;
10 mM 2-mercaptoethanol;
200 µg/ml BSA; and 50% glycerol.
Store at -20°C.**Diluent:** SE Buffer A**Heat inactivation:****No** (65° C for 20 minutes)**Ligation/recutting assay:** After 10-fold
overdigestion with BstC8 I, > 95% of the
DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific
activity was detected after incubation of
1 µg of λ DNA with 20 units of BstC8 I
for 16 hours.

| | | | | | |
|-----------------------------------|-------|-------|-------|--------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 10-25 | 25-50 | 50-75 | 75-100 | 100 |



New product



Mammalian Genome Qualified



New package

BstDE I (prototype Dde I)C[^]TNAG
GANT[^]CE227
E228500 u.a.
2500 u.a.Isolated from *Bacillus stearothermophilus* DE**Concentration:** 10 000 – 20 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer G

Reaction conditions: 1× SE Buffer G

Incubate at 60°C.

Incubation at 37°C results in 10% activity.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

100 mM NaCl; 0,1 mM EDTA;

7 mM 2-mercaptoethanol; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:**

No (65° C for 20 minutes)

Yes (80° C for 20 minutes)

Ligation/recutting assay: After 30-fold

overdigestion with BstDE I, 90% of the

DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific

activity was detected after incubation of

1 µg of λ DNA with 60 units of BstDE I

for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|-----|-------|-------|-------|
| Activity in SE Buffers (% of max) | 75-100 | 100 | 25-50 | 50-75 | 10-25 |

Qualified for human genome studies: http://science.sibenzyme.com/article8_article_31_1.phtml**BstDS I (prototype Dsa I)**C[^]CRYGG
GGYRC[^]CE083
E0841000 u.a.
5000 u.a.Isolated from *Bacillus stearothermophilus* DS**Concentration:** 10 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer Y

Reaction conditions: 1×SE Buffer Y

Incubate at 65°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

100 mM KCl; 0,1 mM EDTA;

7 mM 2-mercaptoethanol;

200 µg/ml BSA; and 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:**

No (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold

overdigestion with BstDS I, 95% of the

DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific

activity was detected after incubation of

1 µg of λ DNA with 20 units of BstDS I

for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|------|--------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 0-10 | 75-100 | 50-75 | 25-50 | 100 |

BstEN I (prototype EcoN I)CCTNN[^]NNNAGG
GGANN[^]NNTCCE103
E104200 e.a.
1000 e.a.Isolated from *Bacillus stearothermophilus* EN**Concentration:** 2 000 – 5 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer Y

Reaction conditions: 1× SE Buffer Y

Incubate at 65°C.

Not blocked by overlapping*dcm*-methylation (C^mCWGG):CCWGGNNNAGG orCCTNNNCCAGG.**Storage buffer:** 10 mM Tris-HCl (pH 7.5); 50

mM KCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol;

200 µg/ml BSA; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:**

No (65° C for 20 minutes)

Yes (80° C for 20 minutes)

Ligation/recutting assay: After 3-fold

overdigestion with BstEN I, about 60%

of the DNA fragments can be ligated. Of

these, 90 % can be recut.

Overdigestion assay: No nonspecific

activity was detected after incubation of

1 µg of λ DNA with 5 units of BstEN I

for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 50-75 | 50-75 | 25-50 | 25-50 | 100 |

BstF5 I (prototype Fok I)GGATGNN[^]
CCTAC[^]NNE031
E032500 u.a.
2500 u.a.Isolated from *Bacillus stearothermophilus* F5**Concentration:** 10 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer Y

Reaction conditions: 1×SE Buffer Y

Incubate at 65°C.

BstF5I is a neoschizomer of FokI.**Storage buffer:** 10 mM Tris-HCl (pH 7.5);

50 mM KCl; 0,1 mM EDTA;

7 mM 2-mercaptoethanol;

200 µg/ml BSA; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:**

No (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold

overdigestion with BstF5 I,

> 90% of the λ DNA fragments can

be ligated. Of these, 95 % can be recut.

Overdigestion assay: No nonspecific

activity was detected after incubation of

1 µg of λ DNA with 20 units of BstF5 I

for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|-------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 75-100 | 50-75 | 25-50 | 50-75 | 100 |

BstFN I (prototype FnuD II)CG[^]CG
GC[^]GCE283
E284300 u.a.
1500 u.a.Isolated from *Bacillus stearothermophilus* FN**Concentration:** 2 000 - 10 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer Y

Reaction conditions: 1×SE Buffer Y

Incubate at 60°C.

Blocked by CG methylation.**Storage buffer:** 20 mM Tris-HCl (pH 7.5);

300 mM NaCl; 0,1 mM EDTA;

10 mM MgCl₂; 7 mM 2-mercaptoethanol;

200 µg/ml BSA; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:**

No (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold

overdigestion with BstFN I, 95% of the

DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific

activity was detected after incubation of

1 µg of λ DNA with 10 units of BstFN I

for 16 hours.


| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|-------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 75-100 | 50-75 | 25-50 | 25-50 | 100 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml

| BstH2 I (prototype Hae II) | RGCGC[^]Y Y[^]CGCGR | E171 E172 | 500 u.a. 2500 u.a. | | |
|---|---|---|-------------------------------|-------|-----|
| Isolated from <i>Bacillus stearothermophilus</i> H2 | | | | | |
| Concentration: 10 000 – 30 000 units/ml Assayed on λ DNA Reagents Supplied with Enzyme: 10×SE Buffer Y, BSA Reaction conditions: 1×SE Buffer Y+BSA Incubate at 65°C. <i>To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.</i> | Storage buffer: 10 mM KH ₂ PO ₄ (pH 7.2); 100 mM NaCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: No (65° C for 20 minutes) | Ligation/recutting assay: After 20-fold overdigestion with BstH2 I, > 90% of the DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of BstH2 I for 16 hours. Star activity: High enzyme concentration results in star activity. Do not use BSA for long incubation. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 50-75 | 50-75 | 0-10 | 10-25 | 100 |

| BstHH I (prototype Hha I) | GCG[^]C C[^]GCG | E143 E144 | 2000 u.a. 10000 u.a. | | |
|---|---|---|---------------------------------|-------|-----|
| Isolated from <i>Bacillus stearothermophilus</i> HH | | | | | |
| Concentration: 50 000 units/ml Assayed on λ DNA Reagents Supplied with Enzyme: 10×SE Buffer Y, BSA Reaction conditions: 1×SE Buffer Y+BSA Incubate at 50°C. Blocked by CG methylation 5'-G(5mC)GC-3'/3'-CG(5mC)G-5' or 5'-G(5mC)GC-3'/3'-CGCG-5' Not blocked by methylation 5'-GCG(5mC)-3'/3'-(5mC)GCG-5' or 5'-GCG(5mC)-3'/3'-CGCG-5'. <i>To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.</i> | Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA, 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: No (65°C, 80° C for 20 minutes) | Ligation/recutting assay: After 40-fold overdigestion with BstHH I, > 90% of the DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 100 units of BstHH I for 16 hours. Do not use BSA for long incubation. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 50-75 | 25-50 | 50-75 | 100 |

| BstKT I (prototype Mbo I) | GAT[^]C C[^]TAG | E151 E152 | 200 u.a. 1000 u.a. | | |
|---|---|--|-------------------------------|-----|-------|
| Isolated from <i>Bacillus stearothermophilus</i> KT | | | | | |
| Concentration: 2 000 – 5 000 units/ml Assayed on λ DNA (dam ⁻) Reagents Supplied with Enzyme: 10×SE Buffer W Reaction conditions: 1 × SE Buffer W Incubate at 37°C. Blocked by overlapping <i>dam</i> -methylation (G ^m ATC): GATC . Not blocked by CG methylation. Cut hemimethylated site: 5'- G ^m ATC-3' / 5'-GATC-3' BstKTI is a neoschizomer of MboI. | Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 100 µg/ml BSA; 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat Inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 5-fold overdigestion with BstKT I, > 90% of the DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 5 units of BstKT I for 16 hours. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 25-50 | 50-75 | 75-100 | 100 | 50-75 |

| BstMA I (prototype BsmA I) |  GTCTC(N)₁[^] CAGAG(N)₅[^] | E291 E292 | 2000 u.a. 10 000 u.a. | | |
|--|---|--|----------------------------------|-----|--------|
| Isolated from <i>Bacillus stearothermophilus</i> MA | | | | | |
| Concentration: 30 000 – 100 000 units/ml Assayed on λ DNA Reagents Supplied with Enzyme: 10×SE Buffer W, BSA Reaction conditions: 1×SE Buffer W+BSA Incubate at 55°C. <i>To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.</i> | Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 100 µg/ml BSA, 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 50-fold overdigestion with BstMA I, >90% of the DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 100 units of BstMA I for 16 hours. Do not use BSA for long incubation. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 25-50 | 50-75 | 50-75 | 100 | 75-100 |

Qualified for human genome studies: http://science.sibenzyme.com/article8_article_31_1.phtml

| BstMB I (prototype Mbo I) | [^]GATC CTAG[^] | E119 E120 | 200 u.a. 1000 u.a. | | | |
|---|--|--|-------------------------------|-----|--------|-------|
| Isolated from <i>Bacillus stearothermophilus</i> MB | | | | | | |
| Concentration: 5 000 - 10 000 units/ml | Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0.1 mM EDTA; | Ligation/recutting assay: After 10-fold overdigestion with BstMB I, > 95% of the DNA fragments can be ligated and recut. | | | | |
| Assayed on λ DNA (dam ⁻) | 7 mM 2-mercaptoethanol; 200 μ g/ml BSA, 50% glycerol. Store at -20°C. | Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 20 units of BstMB I for 16 hours. | | | | |
| Reagents Supplied with Enzyme: 10 \times SE Buffer O | Diluent: SE Buffer A | | | | | |
| Reaction conditions: 1 \times SE Buffer O Incubate at 65°C. | Heat inactivation: No (65° C for 20 minutes) Yes (80° C for 20 minutes) | | | | | |
| Blocked by overlapping <i>dam</i> -methylation (G ^m ATC): GATC . | | | | | | |
| Not blocked by CG methylation. | | | | | | |
| Not cut hemimethylated site: 5'-G ^m ATC-3' / 3'-CTAG-5' | | | | | | |
| | SE Buffers | B | G | O | W | Y |
| | Activity in SE Buffers (% of max) | 10-25 | 25-50 | 100 | 75-100 | 10-25 |

| BstMC I (prototype Mcr I) | CGRY[^]CG GC[^]YRGC | E071 E072 | 500 u.a. 2500 u.a. | | | |
|--|---|--|-------------------------------|-------|-------|-------|
| Isolated from <i>Bacillus stearothermophilus</i> MC | | | | | | |
| Concentration: 5 000 units/ml | Storage buffer: 10 mM Tris-HCl (pH 7.5); 200 mM KCl; 0.1 mM EDTA; | Ligation/recutting assay: After 5-fold overdigestion with BstMC I, > 90% of the DNA fragments can be ligated and recut. | | | | |
| Assayed on λ DNA | 7 mM 2-mercaptoethanol; 200 μ g/ml BSA; 50% glycerol. Store at -20°C. | Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 10 units of BstMC I for 16 hours. | | | | |
| Reagents Supplied with: 10 \times SE Buffer B, BSA | Diluent: SE Buffer A | Do not use BSA for long incubation. | | | | |
| Reaction conditions: 1 \times SE Buffer B+BSA Incubate at 50°C. | Heat inactivation: No (65° C for 20 minutes) | | | | | |
| <i>To obtain 100% activity, BSA should be added to the 1\times reaction mix to a final concentration of 100 μg/ml.</i> | | | | | | |
| | SE Buffers | B | G | O | W | Y |
| | Activity in SE Buffers (% of max) | 100 | 75-100 | 10-25 | 10-25 | 50-75 |

| BstMW I (prototype Mwo I) | GCNNNN[^]NNGC CGNN[^]NNNNCG | E459 E460 | 500 u.a. 2500 u.a. | | | |
|--|--|--|-------------------------------|-------|-------|-----|
| Isolated from <i>Bacillus stearothermophilus</i> MW | | | | | | |
| Concentration: 2 000 – 10 000 units/ml | Storage buffer: 10 mM Tris-HCl(pH 7.5); 50 mM KCl; 0.1 mM EDTA; | Ligation/recutting assay: After 5-fold overdigestion with BstMW I, > 95% of the DNA fragments can be ligated and recut. | | | | |
| Assayed on λ DNA | 7 mM 2-mercaptoethanol; 200 μ g/ml BSA, and 50% glycerol. | Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 10 units of BstMW I for 16 hours. | | | | |
| Reagents Supplied with Enzyme: 10 \times SE Buffer Y | Store at -20°C (*see note). | *Note: For long term storage (more than 7 days), store at -70°C. | | | | |
| Reaction conditions: 1 \times SE Buffer Y Incubate at 55°C. | Diluent: SE Buffer A | | | | | |
| Incubation at 37° results in 20% activity. | Heat inactivation: No (65° C for 20 minutes) Yes (80° C for 20 minutes) | | | | | |
| | SE Buffers | B | G | O | W | Y |
| | Activity in SE Buffers (% of max) | 10-25 | 25-50 | 25-50 | 50-75 | 100 |

| BstNS I (prototype Nsp I) | RCATG[^]Y Y[^]GTACR | E251 E252 | 200 u.a. 1000 u.a. | | | |
|--|---|--|-------------------------------|-------|-------|--------|
| Isolated from <i>Bacillus stearothermophilus</i> NS | | | | | | |
| Concentration: 10 000 units/ml | Storage buffer: 10 mM Tris-HCl (pH 7.5); 250 mM NaCl; 0.1 mM EDTA; | Ligation/recutting assay: After 10-fold overdigestion with BstNS I, > 95% of the DNA fragments can be ligated and recut. | | | | |
| Assayed on λ DNA | 7 mM 2-mercaptoethanol; 100 μ g/ml BSA; 50% glycerol. Store at -20°C. | Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 20 units of BstNS I for 16 hours. | | | | |
| Reagents Supplied with Enzyme: 10 \times SE Buffer B, BSA | Diluent: SE Buffer A | Do not use BSA for long incubation. | | | | |
| Reaction conditions: 1 \times SE Buffer B+BSA Incubate at 37°C. | Heat inactivation : Yes (65° C for 20 minutes) | | | | | |
| <i>To obtain 100% activity, BSA should be added to the 1\times reaction mix to a final concentration of 100 μg/ml.</i> | | | | | | |
| | SE Buffers | B | G | O | W | Y |
| | Activity in SE Buffers (% of max) | 100 | 50-75 | 10-25 | 10-25 | 75-100 |

| BstPA I (prototype PshA I) | GACNN[^]NNGTC CTGNN[^]NNCAG | E299 E300 | 1000 u.a. 5000 u.a. | | | |
|--|---|--|--------------------------------|-------|-------|-----|
| Isolated from <i>Bacillus stearothermophilus</i> PA | | | | | | |
| Concentration: 10 000 –20 000 units/ml | Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0.1 mM EDTA; | Ligation/recutting assay: After 5-fold overdigestion with BstPA I, < 5% of the DNA fragments can be ligated. | | | | |
| Assayed on λ DNA | 7 mM 2-mercaptoethanol; 200 μ g/ml BSA, 50% glycerol. Store at -20°C. | Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 10 units of BstPA I for 16 hours at 25°C. | | | | |
| Reagents Supplied with Enzyme: 10 \times SE Buffer Y | Diluent: SE Buffer A | Star activity: Incubation at 65°C for 16 hours results in star activity. | | | | |
| Reaction conditions: 1 \times SE Buffer Y Incubate at 65°C. | Heat inactivation : No (65° C, 80° C for 20 minutes) | | | | | |
| | SE Buffers | B | G | O | W | Y |
| | Activity in SE Buffers (% of max) | 50-75 | 25-50 | 50-75 | 50-75 | 100 |

BstSC I (prototype ScrF I)Isolated from *Bacillus stearothermophilus* SC**^CCNGG
GGNCC^****E307
E308****100 u.a.
500 u.a.****Concentration:** 2 000 – 5 000 units/ml**Assayed on** λ DNA (dcm-)**Reagents Supplied with Enzyme:**10 \times SE Buffer Y**Reaction conditions:** 1 \times SE Buffer Y

Incubate at 55°C.

Incubation at 37°C results in 10% activity.

Blocked by overlapping *dcm*-methylation(C^mCWGG): **CCWGG**.**BstSCI is a neoschizomer of ScrFI.****Qualified for human genome studies:** http://science.sibenzyme.com/article8_article_31_1.phtml**Storage buffer:** 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0,1 mM EDTA;7 mM 2-mercaptoethanol; 200 μ g/ml BSA; and 50% glycerol. Store at -20°C.**Diluent:** SE Buffer A**Heat Inactivation:****No** (65° C for 20 minutes)**Yes** (80° C for 20 minutes)**Ligation/recutting assay:** After 5-fold overdigestion with BstSC I, > 95% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 3 units of BstSC I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 50-75 | 50-75 | 50-75 | 50-75 | 100 |

BstSF I (prototype Sfe I)Isolated from *Bacillus stearothermophilus* SF**C^TRYAG
GAYRT^C****E197
E198****200 u.a.
1000 u.a.****Concentration:** 2 000 - 5 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**10 \times SE Buffer O, BSA**Reaction conditions:** 1 \times SE Buffer O+BSA

Incubate at 60°C.

*To obtain 100% activity, BSA should be added to the 1 \times reaction mix to a final concentration of 100 μ g/ml.***Storage buffer:** 10 mM Tris-HCl (pH 7.5); 100 mM KCl; 0,1 mM EDTA;7 mM 2-mercaptoethanol; 200 μ g/ml BSA; and 50% glycerol. Store at -20°C.**Diluent:** SE Buffer A**Heat Inactivation:****No**(65°C, 80°C for 20 minutes)**Ligation/recutting assay:** After 3-fold overdigestion with BstSF I, > 95% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 5 units of BstSF I for 16 hours.**Do not use BSA for long incubation.**

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|-------|-----|-------|-------|
| Activity in SE Buffers (% of max) | 75-100 | 25-50 | 100 | 50-75 | 50-75 |

BstSL I (prototype BseS I)Isolated from *Bacillus stearothermophilus* S**GKGCM^C
C^MCGKG****E561
E562****500 u.a.
2500 u.a.****Concentration:** 10 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**10 \times SE Buffer G, BSA**Reaction conditions:** 1 \times SE Buffer G + BSA

Incubate at 55°C.

Not blocked by overlapping *dcm*-methylation (C^mCWGG):**GKGCCWGG****Blocked** by **GKG^mCMC** methylation.*To obtain 100% activity, BSA should be added to the 1 \times reaction mix to a final concentration of 100 μ g/ml.***Storage buffer:** 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0.1 mM EDTA;7 mM 2-mercaptoethanol; 200 μ g/ml BSA, 50% glycerol. Store at -20°C.**Diluent:** SE Buffer A**Heat inactivation:****Yes** (65° C for 20 minutes)**Ligation/recutting assay:** After 5-fold overdigestion with BstSL I, ~80% of the DNA fragments can be ligated. Of these, 95% can be recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 10 units of BstSL I for 16 hours.**Do not use BSA for long incubation.**

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-----|-------|--------|--------|
| Activity in SE Buffers (% of max) | 50-75 | 100 | 50-75 | 75-100 | 75-100 |

BstSN I (prototype SnaB I)Isolated from *Bacillus stearothermophilus* SN**TAC^GTA
ATG^CAT****E065
E066****200 u.a.
1000 u.a.****Concentration:** 5 000 – 10 000 units/ml**Assayed on** T7 DNA**Reagents Supplied with Enzyme:**10 \times SE Buffer B**Reaction conditions:** 1 \times SE Buffer B

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 250 mM NaCl; 0.1 mM EDTA;7 mM 2-mercaptoethanol; 100 μ g/ml BSA, 50% glycerol. Store at -20°C.**Diluent:** SE Buffer A**Heat inactivation:****No** (65° C for 20 minutes)**Ligation/recutting assay:** After 5-fold overdigestion with BstSN I, ~70% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 μ g of T7 DNA with 5 units of BstSN I for 16 hours.**Star activity:** High enzyme

concentration results in star activity.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-----|-------|------|-------|-------|
| Activity in SE Buffers (% of max) | 100 | 50-75 | 0-10 | 10-25 | 50-75 |

BstV1 I (prototype Bbv I)Isolated from *Bacillus stearothermophilus* V1**GCAGC(N)₈[^]
CGTCG(N)₁₂[^]****E303
E304****100 u.a.
500 u.a.****Concentration:** 1 000 – 2 000 units/ml**Assayed on** pBR322 DNA**Reagents Supplied with Enzyme:**10 \times SE Buffer G**Reaction conditions:** 1 \times SE Buffer G

Incubate at 55°C.

Incubation at 37° results in 10% activity.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 200 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 50% glycerol. Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:****No** (65° C for 20 minutes)**Yes** (80° C for 20 minutes)**Ligation/recutting assay:** After 3-fold overdigestion with BstV1 I, > 90% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 μ g of the DNA with 2 units of enzyme for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|-----|--------|--------|--------|
| Activity in SE Buffers (% of max) | 75-100 | 100 | 75-100 | 75-100 | 75-100 |



New product



Mammalian Genome Qualified



New package

BstV2 I (prototype Bbv II)Isolated from *Bacillus stearothermophilus* V2GAAGAC(N)₂[^]
CTTCTG(N)₆[^]E297
E298200 u.a.
1000 u.a.**Concentration:** 5 000 – 15 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Buffer:**

10×SE Buffer Y, BSA

Reaction conditions: 1×SE Buffer Y+BSA

Incubate at 55°C.

*To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.***Storage buffer:** 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 µg/ml BSA, 50% glycerol. Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:**

Yes (65°C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with BstV2 I, > 90% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 5 units of BstV2 I for 16 hours.**Star activity:** High enzyme concentration results in star activity.**Do not use BSA for long incubation.**

| | | | | | |
|-----------------------------------|--------|--------|-------|-------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 25-50 | 25-50 | 100 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml**BstX I (prototype BstX I)**Isolated from *Bacillus stearothermophilus* XCCANNNNN[^]NTGG
GGTN[^]NNNNNACCE465
E466200 u.a.
1000 u.a.**Concentration:** 5 000 -15 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Buffer:**

10×SE Buffer O

Reaction conditions: 1×SE Buffer O

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol. Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:**

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with BstX I, > 95% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of BstX I for 16 hours.

| | | | | | |
|-----------------------------------|-------|-------|-----|--------|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 10-25 | 10-25 | 100 | 75-100 | 25-50 |

BstX2 I (prototype Xho II)Isolated from an *E.coli* strain that carries the cloned BstX2I gene from *Bacillus stearothermophilus* X2R[^]GATCY
YCTAG[^]RE229
E230500 u.a.
2500 u.a.**Concentration:** 10 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer G

Reaction conditions: 1×SE Buffer G

Incubate at 60°C.

Not blocked by overlapping *dam*-methylation (G^mATC): **RGATCY**.**Storage buffer:** 10 mM Tris-HCl (pH 7.6); 100 mM NaCl; 0.1 mM EDTA; 1 mM DTT;

200 µg/ml BSA, 50% glycerol. Store at -20°C.

Diluent: SE Buffer A**Heat Inactivation:**

No (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with BstX2 I, > 95% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 5 units of BstX2 I for 16 hours.

| | | | | | |
|-----------------------------------|--------|-----|------|-------|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 100 | 0-10 | 10-25 | 25-50 |

Bsu I (prototype BciVI)Isolated from *Bacillus sphaericus*GTATCC(N)₆
CATAGG(N)₅E581
E582200 u.a.
1000 u.a.**Concentration:** 2 000 – 5 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer Y

Reaction conditions: 1×SE Buffer Y

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 µg/ml BSA; 0,05% Triton X-100, 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A**Heat Inactivation:**

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with BsuI 10% of the DNA fragments can be ligated with T4 DNA Ligase and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 6 units of BsuI for 16 hours at 37°C.

| | | | | | |
|-----------------------------------|--------|-------|-------|-------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 50-75 | 10-25 | 25-50 | 100 |

BsuR I (prototype Hae III)Isolated from *Bacillus subtilis* RGG[^]CC
CC[^]GGE053
E0541000 u.a.
5000 u.a..**Concentration:** 20 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer G

Reaction conditions: 1×SE Buffer G

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0,1 mM EDTA;

7 mM 2-mercaptoethanol; 200 µg/ml BSA, 50% glycerol. Store at -20°C.

Diluent: SE Buffer A**Heat Inactivation:**

No (65° C for 20 minutes)

Ligation/recutting assay: After 20-fold overdigestion with BsuR I, > 90% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of BsuR I for 16 hours.

| | | | | | |
|-----------------------------------|--------|-----|-------|-------|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 100 | 25-50 | 50-75 | 50-75 |

Btr I (prototype Btr I)

Isolated from *Bacillus stearothermophilus* SE-U62

CAC[^]GTC
GTG[^]CAG

E277
E278

100 u.a.
500 u.a.

Concentration: 2 000 – 5 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer O, BSA

Reaction conditions: 1×SE Buffer O+BSA

Incubate at 60°C.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; and 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

No (65° C for 20 minutes)

Yes (80° C for 20 minutes)

Ligation/recutting assay: After 3-fold overdigestion with Btr I, 80% of the DNA fragments can be ligated. Of these, 90 % can be recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 5 units of Btr I for 16 hours.

Star activity: High enzyme concentration results in star activity.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|--------|-----|--------|--------|
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 100 | 75-100 | 75-100 |

Cci I (prototype BspH I)

Isolated from *Curtobacterium citreum*

T[^]CATGA
AGTAC[^]T

E565
E566

1000 u.a.
5000 u.a.

Concentration: 20000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer W, BSA

Reaction conditions:

1×SE Buffer W + BSA

Incubate at 55°C.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (80° C for 20 minutes)

Ligation/recutting assay: After 20-fold overdigestion with Cci I, > 90% of the DNA fragments can be ligated and recut with.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 40 units of Cci I for 16 hours.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|------|-------|-------|-----|--------|
| Activity in SE Buffers (% of max) | 0-10 | 10-25 | 25-50 | 100 | 75-100 |

CciN I (prototype Not I)

Isolated from *Curtobacterium citreum* N

GC[^]GGCCGC
CGCCGG[^]CG

E203
E204

200 u.a.
1000 u.a.

Concentration: 2000 –5000 units/ml

Assayed on Ad2 DNA

Reagents Supplied with Enzyme:

10×SE Buffer Y

Reaction conditions: 1×SE Buffer Y

Incubate at 37°C.

Blocked by CG methylation.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with CciN I, > 95% of the DNA fragments can be ligated and recut with.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of Ad2 DNA with 10 units of CciN I for 16 hours.

Star activity: High enzyme concentration may result in star activity.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|--------|--------|-----|
| Activity in SE Buffers (% of max) | 25-50 | 50-75 | 75-100 | 75-100 | 100 |

Dra I (prototype Aha III)

Isolated from *Deinococcus radiophilus*

TTT[^]AAA
AAA[^]TTT

E055
E056

1000 u.a.
5000 u.a..

Concentration: 10 000 - 30 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer G, BSA

Reaction conditions: 1×SE Buffer G+BSA

Incubate at 37°C.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 250 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 100 µg/ml BSA, 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with DraI ~70% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 30 units of Dra I for 16 hours.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|-----|-------|--------|--------|
| Activity in SE Buffers (% of max) | 75-100 | 100 | 25-50 | 75-100 | 75-100 |



New product



Mammalian Genome Qualified



New package

Dra III (prototype Dra III)

Isolated from an *E.coli* strain that carries the cloned DraIII gene from *Deinococcus radiophilus*

CACNNN[^]GTG
GTG[^]NNNCAC

E309
E310

500 u.a.
2500 u.a.

Concentration: 5 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer 2K, BSA

Reaction conditions: 1×SE Buffer 2K+BSA

Incubate at 37°C.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 300 mM NaCl; 0.1 mM EDTA;

1 mM DTT; 200 µg/ml BSA, 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with Dra III, about 70% of the DNA fragments can be ligated and recut. In the presence of 10% PEG ligation is better.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of DraIII for 16 hours.

Star activity: High concentration of enzyme may result in star activity.

Do not use BSA for long incubation.

| | | | | | |
|-----------------------------------|-------|-------|--------|--------|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 25-50 | 50-75 | 75-100 | 75-100 | 50-75 |

Dri I (prototype Eam1105 I)

Isolated from *Deinococcus radiophilus* EA

GACNNN[^]NNGTC
CTGNN[^]NNNCAG

E193
E194

200 u.a.
1000 u.a.

Concentration: 5 000 -10 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer Y

Reaction conditions: 1×SE Buffer Y

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 250 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with DriI, about 5% of the DNA fragments can be ligated. In the presence of 10% PEG ligation is better.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of Dri I for 16 hours.

| | | | | | |
|-----------------------------------|--------|--------|-------|-------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 10-25 | 10-25 | 100 |

DseD I (prototype Drd I)

Isolated from *Deinococcus* species D2

GACNNNN[^]NNGTC
CTGNN[^]NNNNCAG

E241
E242

500 u.a.
2500 u.a.

Concentration: 10 000 - 30 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer Y, BSA

Reaction conditions: 1×SE Buffer Y+BSA

Incubate at 37°C.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0,1 mM EDTA;

7 mM 2-mercaptoethanol; 200 µg/ml BSA, 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

No (65°C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with DseD I, > 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of DseD I for 16 hours.

Do not use BSA for long incubation.

| | | | | | |
|-----------------------------------|--------|--------|-------|-------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 25-50 | 50-75 | 100 |

EcoICR I (prototype Sac I)

Isolated from *Escherichia coli* ICR

GAG[^]CTC
CTC[^]GAG

E469
E470

200 u.a.
1000 u.a.

Concentration: 2 000 - 10 000 units/ml

Assayed on λ DNA (Hind III-digest)

Reagents Supplied with Enzyme:

10×SE Buffer G, BSA

Reaction conditions: 1×SE Buffer G+BSA

Incubate at 37°C.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

EcoICRI is a neoschizomer of SacI.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol, 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A

Heat Inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with EcoICR I, > 95% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of EcoICR I for 16 hours.

Do not use BSA for long incubation.

| | | | | | |
|-----------------------------------|--------|-----|------|------|--------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 100 | 0-10 | 0-10 | 75-100 |

EcoR I (prototype EcoR I)

Isolated from an *E.coli* strain that carries the cloned EcoR I gene from *Escherichia coli*



G[^]AATTC
CTTAA[^]G

E057 5000 u.a.
E058 25000 u.a.

For high concentration

E057X 5000 u.a.
E058X 25000 u.a.

Concentration: 20 000 and 50 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

5×SE Buffer EcoR I, BSA

Reaction conditions:

1×SE Buffer EcoRI+BSA

Incubate at 37°C.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 200 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 µg/ml BSA, 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat Inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 50-fold overdigestion with EcoR I, 95% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 40 units of EcoR I for 16 hours.

Star activity: High concentration of enzyme results in star activity.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|--------|--------|--------|-------|
| Activity in SE Buffers (% of max) | 50-75 | 75-100 | 75-100 | 75-100 | 50-75 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml

EcoR V (prototype EcoR V)

Isolated from an *E.coli* strain that carries the cloned EcoR V gene from *Escherichia coli*

GAT[^]ATC
CTA[^]TAG

E059 2000 u.a.
E060 10000 u.a.

Concentration: 20 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer W, BSA

Reaction conditions: 1×SE Buffer W+BSA

Incubate at 37°C.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM NaCl; 0,1 mM EDTA; 1 mM DTT;

200 µg/ml BSA; 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

No (65° C for 20 minutes)

Ligation/recutting assay: After 20-fold overdigestion with EcoR V, 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of EcoR V for 16 hours.

Star activity: High enzyme concentration results in star activity.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|------|-------|-------|-----|-------|
| Activity in SE Buffers (% of max) | 0-10 | 25-50 | 50-75 | 100 | 25-50 |

Ege I (prototype Nar I)

Isolated from *Enterobacter gergoviae*

GGC[^]GCC
CCG[^]CGG

E243 200 u.a.
E244 1000 u.a.

Concentration: 5 000 – 10 000 units/ml

Assayed on λ DNA (HindIII-digest)

Reagents Supplied with Enzyme:

10×SE Buffer B, BSA

Reaction conditions: 1×SE Buffer B+BSA

Incubate at 37°C.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

EgeI is a neoschizomer of NarI.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 250 mM NaCl; 0,1 mM EDTA;

7 mM 2-mercaptoethanol; 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with Ege I, 70% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of the DNA with 10 units of Ege I for 16 hours.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-----|--------|-------|-------|--------|
| Activity in SE Buffers (% of max) | 100 | 75-100 | 10-25 | 50-75 | 75-100 |

Erh I (prototype Sty I)

Isolated from *Erwinia rhapontici*

C[^]CWWGG
GGWWC[^]C

E061 1000 u.a.
E062 5000 u.a.

Concentration: 10 000 - 20 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer 2W, BSA

Reaction conditions:

1 ×SE Buffer 2W+BSA

Incubate at 37°C.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0,1 mM EDTA;

7 mM 2-mercaptoethanol; 200 µg/ml BSA, 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 20-fold overdigestion with Erh I, 95% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 50 units of ErhI for 16 hours.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|-------|--------|-------|
| Activity in SE Buffers (% of max) | 10-25 | 25-50 | 25-50 | 75-100 | 10-25 |



New product



Mammalian Genome Qualified



New package

Fae I (prototype Nla III)

Isolated from *Flavobacterium aquatile* N3

Concentration: 500 - 2 000 units/ml

Assayed on pUC19 DNA

Reagents Supplied with Enzyme:

10×SE Buffer FaeI, BSA

Reaction conditions:

1 ×SE Buffer FaeI+BSA

Incubate at 37°C.

Blocked by C^mATG methylation.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 2-fold overdigestion with Fae I, >90% of the DNA fragments can be ligated with T4 DNA Ligase and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of pUC19 DNA with 1 units of Fae I for 16 hours.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|-------|-------|--------|
| Activity in SE Buffers (% of max) | 25-50 | 50-75 | 10-25 | 10-25 | 75-100 |

Fai I (prototype Fai I)

Isolated from *Flavobacterium aquatile* B15

FaiI cleaves 4 expected recognition sites as well as several other sites with a weaker activity

Concentration: 2 000 units/ml

One unit is defined as the amount of enzyme required to cleave 1 pmol of the double-stranded oligonucleotide with the following structure

5'- CGAGTTCA[^]TAGCTGGGCCCAAC -3'

3'- GCTCAAGT[^]ATCGACCCGGGTG -5'

in 1 hour at 50°C in a total reaction volume of 20 µl.

Reagents Supplied with Enzyme:

10×SE Buffer B

Reaction conditions: 1 ×SE Buffer B

Incubate at 50°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM KCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol. Store at -20°C.

Heat inactivation:

Yes (80° C for 20 minutes)

Ligation/recutting assay: After 3-fold overdigestion with Fai I, about 90% of the pUC19 DNA fragments can be ligated with DNA Ligase and recut.

Note! In the case of long incubation with Fai I DNA can be digested to small oligos.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-----|-------|-------|-------|-------|
| Activity in SE Buffers (% of max) | 100 | 50-75 | 10-25 | 25-10 | 25-50 |

Fal I (prototype Fal I)

Isolated from *Flavobacterium aquatile* Ob10

Concentration: 1 000 - 3 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer W, SAM

Reaction conditions:

1 ×SE Buffer W+SAM

Incubate at 37°C.

To obtain 100% activity, SAM should be added to a final concentration 0.01 mM.

[^](N)₈AAGN₅CTT(N)₁₃[^]

[^](N)₁₃TTCN₅GAA(N)₈[^]

E153

100 u.a.

E154

500 u.a.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 250 mM NaCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 3-fold overdigestion with Fal I, 20% of the DNA fragments can be ligated. Of these, 80 % can be recut. In the presence of 10% PEG ligation is better.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 5 units of Fal I for 16 hours.

Star activity: High enzyme concentration results in star activity.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|------|-------|--------|-----|-------|
| Activity in SE Buffers (% of max) | 0-10 | 25-50 | 75-100 | 100 | 50-75 |

Fat I (prototype Nla III)

Isolated from an *E.coli* strain that carries the cloned Fat I gene from *Flavobacterium aquatile* NL3



[^]CATG
GTAC[^]

E155

100 u.a.

E156

500 u.a.

Concentration: 2 000 - 5 000 units/ml

Assayed on DNA pUC19

Reagents Supplied with Enzyme:

10×SE Buffer G

Reaction conditions: 1×SE Buffer G

Incubate at 55°C.

Blocked by ^mCATG methylation.

FatI is a neoschizomer of NlaIII.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0,1 mM EDTA; 1 mM DTT; 200 µg/ml BSA; 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat Inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 2-fold overdigestion with Fat I, > 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of the DNA with 3 units of Fat I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-----|-------|-------|-------|
| Activity in SE Buffers (% of max) | 10-25 | 100 | 25-50 | 10-25 | 50-75 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml

Fau I (prototype Fau I)

Isolated from an *E.coli* strain that carries the cloned Fau I gene from *Flavobacterium aquatili*

CCCGC(N)₄[^]

E209

100 u.a.

GGGCG(N)₆[^]

E210

500 u.a.

Concentration: 2 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer B

Reaction conditions: 1×SE Buffer B

Incubate at 55°C.

Blocked by CG methylation.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM NaCl; 0.1 mM EDTA; 1 mM DTT; 200 µg/ml BSA and 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A

Heat Inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 2-fold overdigestion with Fau I, > 90% of the DNA fragments can be ligated. Of these, 95 % can be recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 4 units of Fau I for 16 hours.

| | | | | | |
|-----------------------------------|-----|-------|------|------|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 100 | 25-50 | 0-10 | 0-10 | 50-75 |

FauND I (prototype Nde I)

Isolated from an *E.coli* strain that carries the cloned FauND I gene from *Flavobacterium aquatili* ND



CA[^]TATG

E009

1000 u.a.

GTAT[^]AC

E010

5000 u.a.

Concentration: 10 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer Y, BSA

Reaction conditions: 1× SE Buffer Y+BSA

Incubate at 37°C.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0,1 mM EDTA; 1mM DTT; 200 µg/ml BSA, 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with FauND I, 80% of the DNA fragments can be ligated and recut. In the presence of 10% PEG ligation is better.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of FauND I for 16 hours.

Note: *Sensitive to impurities present in some DNA preparations. For example, DNA purified by standard miniprep procedures is cleaved at lower rate.*

Do not use BSA for long incubation.

| | | | | | |
|-----------------------------------|-------|--------|-------|-------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 50-75 | 75-100 | 10-25 | 50-75 | 100 |

Qualified for human genome studies: http://science.sibenzyme.com/article8_article_31_1.phtml

Fbl I (prototype Acc I)

Isolated from an *E.coli* strain that carries the cloned Fbl I gene from *Flavobacterium balustinum*

GT[^]MKAC

E271

100 u.a.

CAKM[^]TG

E272

500 u.a.

Concentration: 1 000 – 2 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer Y

Reaction conditions: 1×SE Buffer Y

Incubate at 55°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 200 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA, 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat Inactivation:

No (65° C for 20 minutes)

Ligation/recutting assay: After 2-fold overdigestion with Fbl I, 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 2 units of Fbl I for 16 hours.

| | | | | | |
|-----------------------------------|-------|--------|------|-------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 50-75 | 75-100 | 0-10 | 50-75 | 100 |

Fok I (prototype Fok I)

Isolated from an *E.coli* strain that carries the cloned Fok I gene from *Flavobacterium okeanoikoites*

GGATG(N)₉[^]

E247

100 u.a.

CCTAC(N)₁₃[^]

E248

500 u.a.

Concentration: 1 000 – 2 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer Y

Reaction conditions: SE Buffer Y

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM NaCl; 0,1 mM EDTA; 1 mM DTT ; 200 µg/ml BSA, 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 2-fold overdigestion with Fok I, > 95% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 1 units of Fok I for 16 hours.

Note: *Overdigestions of > 5 units of Fok I per 1 µg of DNA and incubation times > 2 hours are not recommended.*

| | | | | | |
|-----------------------------------|-------|-------|-------|-------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 50-75 | 50-75 | 25-50 | 25-50 | 100 |



New product



Mammalian Genome Qualified



New package

FriO I (prototype Ban II) **GRGCY[^]C**
C[^]YCGRG **E157**
E158 **1000 u.a.**
5000 u.a.

Isolated from *Flavobacterium rigense* O

Concentration: 10 000 - 40 000 units/ml **Storage buffer:** 10 mM Tris-HCl (pH 7.5); 200 mM KCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA, 50% glycerol. Store at -20°C.

Assayed on λ DNA **Ligation/recutting assay:** After 20-fold overdigestion with FriO I, 90% of the DNA fragments can be ligated and recut.

Reagents Supplied with Enzyme: 10×SE Buffer Y, BSA **Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 40 units of FriOI for 16 hours.

Reaction conditions: 1×SE Buffer Y+BSA Incubate at 37°C. **Diluent:** SE Buffer A **Heat inactivation:** Yes (65° C for 20 minutes)

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|--------|-------|------|-----|
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 10-25 | 0-10 | 100 |

Fsp4H I (prototype Fnu4H I) **GC[^]NGC**
CGN[^]CG **E095**
E096 **200 u.a.**
1000 u.a.

Isolated from an *E.coli* strain that carries the cloned Fsp4H I gene from *Flavobacterium* species 4H

Concentration: 3 000 – 5 000 units/ml **Storage buffer:** 10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0.1 mM EDTA; 1 mM DTT; 200 µg/ml BSA; 50% glycerol. Store at -20°C.

Assayed on λ DNA **Ligation/recutting assay:** After 5-fold overdigestion with Fsp4H I, about 5% of the DNA fragments can be ligated and recut.

Reagents Supplied with Enzyme: 10×SE Buffer Y **Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of Fsp4H I for 16 hours.

Reaction conditions: 1× SE Buffer Y Incubate at 37°C. **Diluent:** SE Buffer A **Heat inactivation:** Yes (65° C for 20 minutes)

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|--------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 50-75 | 75-100 | 10-25 | 25-50 | 100 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml

Gla I (prototype Glai) **Pu(5mC)[^]GPy**
PyG[^](5mC)Pu **E493**
E494 **100 u.a.**
500 u.a.

Isolated from *Glacial ice bacterium* GL 29

See page 54 for more information about this enzyme.

Glu I (prototype GluI) **G(5mC)[^]NG(5mC)**
(5mC)GN[^](5mC)G **E519**
E520 **100 u.a.**
500 u.a.

Isolated from *Glacial ice bacterium* GL 24

See page 55 for more information about this enzyme.

Gsa I (prototype BseY I) **CCCAG[^]C**
G[^]GGTCC **E563**
E564 **1000 u.a.**
5000 u.a.

Isolated from *Geobacillus stearothermophilus* Y

Concentration: 10 000-20 000 units/ml **Storage buffer:** 10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol. Store at -20°C.

Assayed on λ DNA **Ligation/recutting assay:** After 20-fold overdigestion with Gsa I, > 90% of the DNA fragments can be ligated and recut.

Reagents Supplied with Enzyme: 10×SE Buffer W, BSA **Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of Gsa I for 16 hours.

Reaction conditions: 1×SE Buffer W+BSA Incubate at 70°C. **Diluent:** SE Buffer A **Heat inactivation:** No (80° C for 20 minutes)

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|--------|-----|--------|
| Activity in SE Buffers (% of max) | 10-25 | 25-50 | 75-100 | 100 | 75-100 |

Hae III (prototype Hae III) **GG[^]CC**
CC[^]GG **E067**
E068 **2000 u.a.**
10000 u.a.

Isolated from an *E.coli* strain that carries the cloned Hae III gene from *Haemophilus aegyptius*

Concentration: 10 000 and 50 000 units/ml **Storage buffer:** 10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0.1 mM EDTA; 1 mM DTT; 200 µg/ml BSA; 50% glycerol. Store at -20°C.

Assayed on λ DNA **Ligation/recutting assay:** After 20-fold overdigestion with Hae III, > 90% of the DNA fragments can be ligated and recut.

Reagents Supplied with Enzyme: 10×SE Buffer G **Overdigestion assay:** No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of Hae III for 16 hours.

Reaction conditions: 1× SE Buffer G Incubate at 37°C. **Diluent:** SE Buffer A **Heat inactivation:** No (65° C for 20 minutes)
Yes (80° C for 20 minutes)

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|-----|-------|-------|-------|
| Activity in SE Buffers (% of max) | 75-100 | 100 | 25-50 | 50-75 | 50-75 |

Qualified for human genome studies: http://science.sibenzyme.com/article8_article_27_1.phtml

Hga I (prototype Hga I)

Isolated from *E.coli* strain that carries the cloned HgaI gene from *Haemophilus gallinarum*

GACGC(N)₅[^]
CTGCG(N)₁₀[^]

E461
E462

50 u.a.
250 u.a.

Concentration: 1 000 units/ml

Assayed on pBR322 DNA

Reagents Supplied with Enzyme:

10×SE Buffer B

Reaction conditions: 1× SE Buffer B

Incubate at 37°C.

Blocked by CG methylation.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
50 mM NaCl; 0.1 mM EDTA; 1 mM DTT;
200 µg/ml BSA; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 3-fold overdigestion with Hga I, > 90% of the DNA fragments can be ligated and recut.

Note: Incubation with >2 units of HgaI per 1 µg of DNA and digestion > 1 hour is not recommended.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-----|--------|-------|-------|-------|
| Activity in SE Buffers (% of max) | 100 | 75-100 | 10-25 | 25-50 | 50-75 |

Hind II (prototype Hind II)

Isolated from an *E.coli* strain that carries the cloned Hind II gene from *Haemophilus influenzae* Rd

GTY[^]RAC
CAR[^]YTG

E201
E202

1000 u.a.
5000 u.a.

Concentration: 10 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer G, BSA

Reaction conditions: 1×SE Buffer G+BSA

Incubate at 37°C.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
200 mM NaCl; 0.1 mM EDTA; 1 mM DTT;
200 µg/ml BSA; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Hind II, > 60% of the DNA fragments can be ligated and recut. In the presence of 10% PEG ligation is better.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of Hind II for 16 hours.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|-----|-------|-------|--------|
| Activity in SE Buffers (% of max) | 75-100 | 100 | 25-50 | 25-50 | 75-100 |

Hind III (prototype Hind III)

Isolated from an *E.coli* strain that carries the cloned Hind III gene from *Haemophilus influenzae* Rd



A[^]AGCTT
TTCGA[^]A

E073
E074

5000 u.a.
25000 u.a.

For high concentration

E073X 5000 u.a.

E074X 25000 u.a.

Concentration: 20 000 and 100 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer W, BSA

Reaction conditions:

1 ×SE Buffer W+BSA

Incubate at 37°C.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
250 mM NaCl; 0.1 mM EDTA;
7 mM 2-mercaptoethanol; 200 µg/ml BSA,
50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

No (65° C for 20 minutes)

Yes (80° C for 20 minutes)

Ligation/recutting assay: After 50-fold overdigestion with Hind III, 95% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 40 units of Hind III for 16 hours.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|------|-----|------|
| Activity in SE Buffers (% of max) | 10-25 | 25-50 | 0-10 | 100 | 0-10 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml

Hinf I (prototype Hinf I)

Isolated from an *E.coli* strain that carries the cloned Hinf I gene from *Haemophilus influenzae*



G[^]ANTC
CTNA[^]G

E075
E076

2000 u.a.
10000 u.a.

For high concentration

E076X 10000 u.a.

Concentration: 20 000 and 40 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer O

Reaction conditions: 1× SE Buffer O

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
50 mM NaCl; 0.1 mM EDTA; 1 mM DTT;
200 µg/ml BSA; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

No (65° C for 20 minutes)

Yes (80° C for 20 minutes)

Ligation/recutting assay: After 20-fold overdigestion with Hinf I, > 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 40 units of Hinf I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|--------|-----|--------|--------|
| Activity in SE Buffers (% of max) | 25-50 | 75-100 | 100 | 75-100 | 75-100 |

Qualified for human genome studies: http://science.sibenzyme.com/article8_article_31_1.phtml



New product



Mammalian Genome Qualified



New package

Hpa I (prototype Hpa I)

Isolated from an *E.coli* strain that carries the cloned Hpa I gene from *Haemophilus parainfluenzae*

GTT[^]AAC
CAA[^]TTG

E077
E078

500 u.a.
2500 u.a.

Concentration: 5 000 units/ml
Assayed on λ DNA
Reagents Supplied with Enzyme:
10 \times SE Buffer Y
Reaction conditions: 1 \times SE Buffer Y
Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
50 mM NaCl; 0.1 mM EDTA; 1 mM DTT;
200 μ g/ml BSA, 50% glycerol. Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
Yes (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with Hpa I, 60% of the DNA fragments can be ligated and recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 5 units of Hpa I for 16 hours.
Star activity: High enzyme concentration results in star activity.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|------|-------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 0-10 | 50-75 | 10-25 | 25-50 | 100 |

Hpa II (prototype Hpa II)

Isolated from an *E.coli* strain that carries the cloned Hpa II gene from *Haemophilus parainfluenzae*



C[^]CGG
GGC[^]C

E161
E162

500 u.a.
2500 u.a.

Concentration: 10 000 units/ml
Assayed on λ DNA
Reagents Supplied with Enzyme:
10 \times SE Buffer B
Reaction conditions: 1 \times SE Buffer B
Incubate at 37°C.
Blocked by CG methylation.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
50 mM NaCl; 0.1 mM EDTA; 1 mM DTT;
100 μ g/ml BSA; 50% glycerol.
Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
No (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Hpa II, > 95% of the DNA fragments can be ligated and recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 10 units of Hpa II for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-----|-------|-------|-------|-------|
| Activity in SE Buffers (% of max) | 100 | 50-75 | 10-25 | 25-50 | 50-75 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml

HpySE526 I (prototype Mae II)

Isolated from an *E.coli* strain that carries the cloned HpySE526 I gene from *Helicobacter pylori* SE526.



A[^]CGT
TGC[^]A

E583
E584

200 u.a.
1000 u.a.

Concentration: 5 000 units/ml
Assayed on pUC19 DNA
Reagents Supplied with Enzyme:
10 \times SE Buffer Y
Reaction Conditions: 1 \times SE Buffer Y
Incubate at 37°C.
Blocked by CG methylation

Storage buffer: 10 mM Tris-HCl (pH 7.6);
100 mM NaCl; 0.1 mM EDTA;
200 μ g/ml BSA; 50% glycerol.
Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
Yes (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with HpySE526 I about 95% of the DNA fragments can be ligated with T4 DNA Ligase and recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 5 units of HpySE526 I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|--------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 10-25 | 25-50 | 100 |

HspA I (prototype Hha I)

Isolated from an *E.coli* strain that carries the cloned HspA I gene from *Haemophilus* species A1



G[^]CGC
CGC[^]G

E069
E070

1000 u.a.
5000 u.a.

Concentration: 20 000 units/ml
Assayed on λ DNA
Reagents Supplied with Enzyme:
10 \times SE Buffer Y
Reaction Conditions: 1 \times SE Buffer Y
Incubate at 37°C.
Blocked by CG methylation
5'-G(5mC)GC-3'/3'-CG(5mC)G-5'.
Not blocked by methylation
5'-GCG(5mC)-3'/3'-CGCG-5' or
5'-GCG(5mC)-3'/3'-(5mC)GCG-5'.
HspAI is a neoschizomer of HhaI.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
50 mM NaCl; 0.1 mM EDTA;
7 mM 2-mercaptoethanol; 200 μ g/ml BSA;
50% glycerol. Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
No (65° C for 20 minutes)

Ligation/recutting assay: After 20-fold overdigestion with HspA I, > 90% of the DNA fragments can be ligated and recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 20 units of HspA I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 50-75 | 50-75 | 25-50 | 25-50 | 100 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml

Kpn I (prototype Kpn I)

Isolated from an *E.coli* strain that carries the cloned Kpn I gene from *Klebsiella pneumonia*



GGTAC[^]C
C[^]CATGG

E079
E080

2000 u.a.
10000 u.a.

For high concentration

E079X 2000 u.a.
E080X 10000 u.a.

Concentration: 20 000 and 40 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer B, BSA

Reaction conditions: 1 \times SE Buffer B+BSA

Incubate at 37°C.

Not blocked by overlapping *dcm*-methylation (C^mCWGG): **GGTACCWGG**.

To obtain 100% activity, BSA should be added to the 1 \times reaction mix to a final concentration of 100 μ g/ml.

Qualified for human genome studies: http://science.sibenzyme.com/article8_article_31_1.phtml

Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM NaCl; 0,1 mM EDTA; 1 mM DTT ; 200 μ g/ml BSA; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

No (65° C for 20 minutes)

Ligation/recutting assay: After 20-fold overdigestion with Kpn I, > 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 20 units of Kpn I for 16 hours.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-----|-------|-------|-------|--------|
| Activity in SE Buffers (% of max) | 100 | 25-50 | 25-50 | 25-50 | 75-100 |

Kro I (prototype Kro I)

Isolated from *Kocurea rosea* 307

G[^]C(5mC)GGC
CGG(5mC)C[^]G

E541
E542

50 u.a.
250 u.a.

See page 56 for more information about this enzyme.

Ksp22 I (prototype Bcl I)

Isolated from *Kurthia* species 22



T[^]GATCA
ACTAG[^]T

E081
E082

1000 u.a.
5000 u.a.

Concentration: 10 000 - 30 000 units/ml

Assayed on λ DNA (dam-)

Reagents Supplied with Enzyme:

10 \times SE Buffer 2K, BSA

Reaction conditions: 1 \times SE Buffer 2K+BSA

Incubate at 37°C.

Blocked by overlapping *dam*-methylation (G^mATC): **TGATCA**.

To obtain 100% activity, BSA should be added to the 1 \times reaction mix to a final concentration of 100 μ g/ml.

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml

Storage buffer: 10 mM Tris-HCl (pH 7.5); 250 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 200 μ g/ml BSA, 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Ksp22 I, >90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 20 units of Ksp22 I for 16 hours.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-----|-------|-------|-------|
| Activity in SE Buffers (% of max) | 50-75 | 100 | 50-75 | 50-75 | 25-50 |

Kzo9 I (prototype Mbo I)

Isolated from *Kurthia zopfii* 9



[^]GATC
CTAG[^]

E187
E188

200 u.a.
1000 u.a.

Concentration: 1 000 – 5 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer G

Reaction conditions: 1 \times SE Buffer G

Incubate at 37°C.

Not blocked by overlapping *dam*-methylation (G^mATC): **GATC**.

Blocked by CG methylation.

Cleaved of DNA is impaired by overlapping CG methylation: **GAT^mCG**.

Qualified for human genome studies: http://science.sibenzyme.com/article8_article_27_1.phtml

Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 μ g/ml BSA; and 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat Inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with Kzo9 I, > 95% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 2 units of Kzo9I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-----|-------|-------|-------|
| Activity in SE Buffers (% of max) | 50-75 | 100 | 50-75 | 50-75 | 50-75 |

Lmn I (prototype Lmn I)

Isolated from *Lysinibacillus manganicus* An22



GCTCCN[^]
CGAG[^]GN

E593
E594

50 e.a.
250 e.a.

Concentration: 500 – 1 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer B

Reaction conditions: 1 \times SE Buffer B

Incubate at 37°C.

Storage buffer: 20 mM Tris-HCl (pH 7.5); 250 mM NaCl; 100 μ g/ml BSA; 7 mM 2-mercaptoethanol; and 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat Inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 3-fold overdigestion with Lmn I 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 2 units of Lmn I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-----|--------|-------|-------|--------|
| Activity in SE Buffers (% of max) | 100 | 75-100 | 50-75 | 50-75 | 75-100 |



New product



Mammalian Genome Qualified



New package

Mab I (prototype SexA I)

Isolated from *Microbacterium arborescens* SE

A[^]CCWGGT
TGGWCC[^]A

E121
E122

200 u.a.
1000 u.a.

Concentration: 1 000 – 5 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer W, BSA

Reaction conditions: 1×SE Buffer W+BSA

Incubate at 37°C.

Not blocked by overlapping

dcm-methylation

(C^mC(A/T)GG): **ACCWGGT**.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

100 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 µg/ml BSA,

50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold

overdigestion with Mab I, 90% of the

DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific

activity was detected after incubation of

1 µg of λ DNA with 10 units of Mab I

for 16 hours.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|--------|-----|-------|
| Activity in SE Buffers (% of max) | 25-50 | 50-75 | 75-100 | 100 | 50-75 |

Mal I (prototype Dpn I)

Isolated from *Marinococcus albus* I

G(mA)[^]TC
CT[^](mA)G

E489
E490

50 u.a.
250 u.a.

See page 56 for more information about this enzyme.

Mbo II (prototype Mbo II)

Isolated from an *E.coli* strain that carries the cloned Mbo II gene from *Moraxella bovis*

GAAGA(N)₈[^]
CTTCT(N)₇[^]

E471
E472

200 u.a.
1000 u.a.

Concentration: 5 000 units/ml

Assayed on λ DNA (*dam*⁻)

Reagents Supplied with Enzyme:

10×SE Buffer Y

Reaction conditions: 1×SE Buffer Y

Incubate at 37°C.

Blocked by overlapping *dam*-methylation

(G^mATC): **GAAGATC**.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

50 mM NaCl; 0.1 mM EDTA; 1 mM DTT;

200 µg/ml BSA; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold

overdigestion with Mbo II, 60% of the

DNA fragments can be ligated and recut.

In presence of 10% PEG ligation is

better.

Overdigestion assay: No nonspecific

activity was detected after incubation of

1 µg of DNA with 5 units of enzyme for

16 hours at 37°C.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|--------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 25-50 | 50-75 | 100 |

Mfe I (prototype Mfe I)

Isolated from an *E.coli* strain that carries the cloned MfeI gene from *Mycoplasma fermentans*

C[^]AATTG
GTAA[^]C

E295
E296

1000 u.a.
5000 u.a.

Concentration: 20 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer B, BSA

Reaction conditions: 1×SE Buffer B+BSA

Incubate at 37°C.

To obtain 100% activity, BSA should be

added to the 1× reaction mix to a final

concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.6);

50 mM NaCl; 0.1 mM EDTA; 1 mM DTT;

200 µg/ml BSA, 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

No (65° C for 20 minutes)

Ligation/recutting assay: After 20-fold

overdigestion with Mfe I, about 90% of

the DNA fragments can be ligated and

recut.

Overdigestion assay: No nonspecific

activity was detected after incubation of

1 µg of λ DNA with 20 units of Mfe I

for 16 hours.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-----|--------|-------|-------|--------|
| Activity in SE Buffers (% of max) | 100 | 75-100 | 10-25 | 25-50 | 75-100 |

Mhl I (prototype Sdu I)

Isolated from *Micrococcus halobius* SD

GDGCH[^]C
C[^]HCGDG

E049
E050

500 u.a.
2500 u.a.

Concentration: 5 000 – 10 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer W

Reaction conditions: 1 × SE Buffer W

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

100 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

No (65° C for 20 minutes)

Yes (80° C for 20 minutes)

Ligation/recutting assay: After 10-fold

overdigestion with Mhl I, >90% of the

DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific

activity was detected after incubation of

1 µg of λ DNA with 10 units of Mhl I

for 16 hours.

Star activity: High enzyme

concentration results in star activity.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|--------|-----|-------|
| Activity in SE Buffers (% of max) | 10-25 | 25-50 | 75-100 | 100 | 10-25 |

Mlu I (prototype Mlu I)

Isolated from *Micrococcus luteus*

A[^]CGCGT
TGCGC[^]A

E085
E086

1000 u.a.
5000 u.a.

Concentration: 10 000 – 30 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer O

Reaction conditions: 1× SE Buffer O

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

100 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 µg/ml BSA;

50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold

overdigestion with Mlu I, 95% of the

DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific

activity was detected after incubation of

1 µg of λ DNA with 20 units of Mlu I

for 16 hours.

| | | | | | |
|-----------------------------------|------|-------|-----|-------|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 0-10 | 10-25 | 100 | 25-50 | 10-25 |

Mly113 I (prototype Nar I)

Isolated from *Micrococcus lylae* 113

Concentration: 3 000 – 5 000 units/ml
Assayed on T7 DNA
Reagents Supplied with Enzyme:
 10×SE Buffer B
Reaction conditions: 1× SE Buffer B
 Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
 50 mM KCl; 0.1 mM EDTA;
 7 mM 2-mercaptoethanol; 200 µg/ml BSA;
 50% glycerol. Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
 Yes (65° C for 20 minutes)

Ligation/recutting assay: After 3-fold overdigestion with Mly113 I, > 80% of the DNA fragments can be ligated and recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of T7 DNA with 10 units of Mly113 I for 16 hours.
Star activity: High enzyme concentration may result in star activity.

| | | | | | |
|-----------------------------------|-----|-------|-------|-------|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 100 | 25-50 | 10-25 | 10-25 | 50-75 |

Mnl I (prototype Mnl I)

Isolated from an *E.coli* strain that carries the cloned MnlI gene from *Moraxella nonliquefaciens*

Concentration: 10 000 units/ml
Assayed on λ DNA
Reagents Supplied with Enzyme:
 10×SE Buffer G, BSA
Reaction conditions: 1×SE Buffer G+BSA
 Incubate at 37°C.
Blocked by overlapping CG methylation:
CCT^mCG.
To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
 200 mM NaCl; 0,1 mM EDTA; 1 mM DTT;
 200 µg/ml BSA, 50% glycerol.
 Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
 Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Mnl I, about 50% of the DNA fragments can be ligated and recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of Mnl I for 16 hours.
Do not use BSA for long incubation.

| | | | | | |
|-----------------------------------|--------|-----|-------|-------|--------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 100 | 25-50 | 25-50 | 75-100 |

Mox20 I (prototype Bal I)

Isolated from *Microbacterium oxydans*

Concentration: 10 000 – 30 000 units/ml
Assayed on λ DNA (dcm-)
Reagents Supplied with Enzyme:
 10×SE Buffer O
Reaction conditions: 1×SE Buffer O
 Incubate at 37°C.
Blocked by overlapping *dcm*-methylation
 (C^mCWGG) : **TGGCCAGG.**

Storage buffer: 10 mM Tris-HCl (pH 7.5);
 200 mM KCl; 0,1 mM EDTA;
 7 mM 2-mercaptoethanol; 200 µg/ml BSA;
 50% glycerol. Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
 No (80° C for 20 minutes)

Ligation/recutting assay: After 20-fold overdigestion with Mox20 I, 80% of the DNA fragments can be ligated and recut. In the presence of 10% PEG ligation is better.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 40 units of Mox20 I for 16 hours..

| | | | | | |
|-----------------------------------|-------|-------|-----|--------|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 10-25 | 25-50 | 100 | 75-100 | 25-50 |

Qualified for human genome studies: http://science.sibenzyme.com/article8_article_31_1.phtml

MroN I (prototype Nae I)

Isolated from *Micrococcus roseus* NO

Concentration: 2 000 –10 000 units/ml
Assayed on Adenovirus-2 DNA
Reagents Supplied with Enzyme:
 10×SE Buffer B
Reaction conditions: 1×SE Buffer B
 Incubate at 37°C.
MroNI is a neoschizomer of NaeI.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
 250 mM NaCl; 0.1 mM EDTA;
 7 mM 2-mercaptoethanol; 200 µg/ml BSA;
 50% glycerol. Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
 No (65°C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with MroN I, > 90% of the Ad-2 DNA fragments can be ligated and recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of Ad-2 DNA with 5 units of MroN I for 16 hours.

| | | | | | |
|-----------------------------------|-----|-------|-------|------|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 100 | 50-75 | 10-25 | 0-10 | 10-25 |

MroX I (prototype Xmn I)

Isolated from *Micrococcus roseus* X


Concentration: 5 000 - 15 000 units/ml
Assayed on λ DNA
Reagents Supplied with Enzyme:
 10×SE Buffer W
Reaction mixture: 1×SE Buffer W
 Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
 200 mM NaCl; 0,1 mM EDTA;
 7 mM 2-mercaptoethanol; 200 µg/ml BSA;
 50% glycerol. Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
 Yes (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with MroX I, 50% of the DNA fragments can be ligated and recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of MroX I for 16 hours.

| | | | | | |
|-----------------------------------|-------|-------|-------|-----|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 50-75 | 50-75 | 50-75 | 100 | 25-50 |

Qualified for human genome studies: http://science.sibenzyme.com/article8_article_31_1.phtml

| | | | | | |
|--|---|--|--------------|------------------------|--------|
| Msp I (prototype Hpa II) Isolated from <i>Moraxella</i> species |  | C [^] CGG GGC [^] C | E091 E092 | 1000 u.a. 5000 u.a. | |
| Concentration: 10 000 - 20 000 units/ml Assayed on λ DNA Reagents Supplied with Enzyme: 10×SE Buffer B Reaction conditions: 1×SE Buffer B Incubate at 37°C. | Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 10-fold overdigestion with Msp I, 95% of the DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of DNA with 20 units of Msp I for 16 hours. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 100 | 75-100 | 50-75 | 75-100 | 75-100 |

Qualified for human genome studies: http://science.sibenzyme.com/article8_article_27_1.phtml

| | | | | |
|---|------------------------|--|--|--|
| Msp20 I (prototype Bal I) Isolated from <i>Micrococcus</i> species 20 | Replaced by Mox20 I | TGG [^] CCA ACC [^] GGT | | |
|---|------------------------|--|--|--|

| | | | | | |
|---|--|--|-----------------------|-------|-----|
| MspA1 I (prototype NspB II) Isolated from <i>Moraxella</i> species A1 | CMG [^] CKG GKC [^] GMC | E191 E192 | 500 u.a. 2500 u.a. | | |
| Concentration: 5 000 – 10 000 units/ml Assayed on λ DNA Reagents Supplied with Enzyme: 10×SE Buffer Y, BSA Reaction conditions: 1× SE BufferY + BSA Incubate at 37°C. <i>To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.</i> | Storage buffer: 20 mM Tris-HCl (pH 7.6); 300 mM NaCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 10 mM MgCl ₂ , 200 µg/ml BSA, 50% glycerol. Store at -20°C (*see note). Diluent: SE Buffer A Heat inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 10-fold overdigestion with MspA1 I, 90% of the DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of MspA1 I for 16 hours. Do not use BSA for long incubation *Note: Store For long term storage (more than 30 days), store at -70°C . | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 10-25 | 75-100 | 10-25 | 25-50 | 100 |

| | | | | | |
|---|--|--|------------------------|-------|--------|
| MspR9 I (prototype ScrF I) Isolated from <i>Moraxella</i> species R9 | CC [^] NGG GGN [^] CC | E175 E176 | 1000 u.a. 5000 u.a. | | |
| Concentration: 10 000 – 20 000 units/ml Assayed on λ DNA (dcm-) Reagents Supplied with Enzyme: 10×SE Buffer O Reaction conditions: 1×SE Buffer O Incubate at 37°C. Blocked by overlapping dcm-methylation (C^mCWGG): <u>CCWGG</u> . | Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 100 µg/ml BSA; 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat Inactivation: No (65° C for 20 minutes) | Ligation/recutting assay: After 2-fold overdigestion with MspR9 I none of the DNA fragments can be ligated. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 30 units of MspR9I for 16 hours. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 50-75 | 50-75 | 100 | 50-75 | 75-100 |

| | | | |
|---|--|--------------|-----------------------|
| Mte I (prototype Mte I) Isolated from <i>Microbacterium testaceum</i> 17B | G(5mC)G(5mC) [^] NG(5mC)G(5mC) (5mC)G(5mC)GN [^] (5mC)G(5mC)G | E553 E554 | 500 u.a. 2500 u.a. |
|---|--|--------------|-----------------------|

See page 57 for more information about this enzyme.

| | | | | | |
|--|---|--|-----------------------|-----|-------|
| Nru I (prototype Nru I) Isolated from <i>Nocardia rubra</i> | TCG [^] CGA AGC [^] GCT | E099 E100 | 500 u.a. 2500 u.a. | | |
| Concentration: 5 000 - 10 000 units/ml Assayed on λ DNA (dam-) Reagents Supplied with Enzyme: 10×SE Buffer W Reaction conditions: 1× SE Buffer W Incubate at 37°C. Blocked by overlapping dam-methylation (G^mATC): <u>TCGCGATC</u> . | Storage buffer: 10 mM Tris-HCl (pH 7.5); 250 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: No (65° C for 20 minutes) | Ligation/recutting assay: After 5-fold overdigestion with Nru I, 20% of the DNA fragments can be ligated. Of these, 90% can be recut. In the presence of 10% PEG ligation is better. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of Nru I for 16 hours. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 0-10 | 10-25 | 75-100 | 100 | 10-25 |

PalA I (prototype Asc I)Isolated from *Pseudomonas alcaligenes* BS17**GG[^]CGCGCC
CCGCGC[^]GG****E483
E484****100 u.a.
500 u.a.****Concentration:** 500 – 2 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**10 \times SE Buffer Y**Reaction condition:** 1 \times SE Buffer Y

Incubate at 37°C.

Blocked by CpG methylation.**Storage buffer:** 10 mM Tris-HCl (pH 7.5);

100 mM KCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 μ g/ml BSA,

50% glycerol. Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:****Yes** (65° C for 20 minutes)**Ligation/recutting assay:** After 2-fold overdigestion with PalAI,

> 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 2 units of PalAI for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|---|---|-----|
| Activity in SE Buffers (% of max) | 25-50 | 10-25 | 0 | 0 | 100 |

Pce I (prototype Stu I)Isolated from *Planococcus citreus* 55**AGG[^]CCT
TCC[^]GGA****E105
E106****1000 u.a.
5000 u.a.****Concentration:** 10 000 – 20 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**10 \times SE Buffer Y**Reaction Condition:** 1 \times SE Buffer Y

Incubate at 50°C.

Incubation at 37°C results in 50% activity.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

100 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 100 μ g/ml BSA;

50% glycerol. Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:****No** (65° C for 20 minutes)**Yes** (80° C for 20 minutes)**Ligation/recutting assay:** After 10-fold overdigestion with Pce I 70% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 20 units of Pce I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|--------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 50-75 | 25-50 | 100 |

Pci I (prototype BspLU11 I)Isolated from an *E.coli* strain that carries the cloned Pci I gene from *Planococcus citreus* SE-F45**A[^]CATGT
TGTAC[^]A****E275
E276****300 u.a.
1500 u.a.****Concentration:** 10 000 units/ml**Assayed on** T7 DNA**Reagents Supplied with Enzyme:**10 \times SE Buffer O**Reaction conditions:** 1 \times SE Buffer O

Incubate at 37°C.

Blocked by ^mACATGT methylation.**Not blocked by** AC^mATGT methylation.**Storage buffer:** 10 mM Tris-HCl (pH 7.5);

250 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; and 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A**Heat Inactivation:****Yes** (65° C for 20 minutes)**Ligation/recutting assay:** After 10-fold overdigestion with Pci I, 90% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 μ g of T7 DNA with 20 units of Pci I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|--------|-----|--------|-------|
| Activity in SE Buffers (% of max) | 50-75 | 75-100 | 100 | 75-100 | 50-75 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml**PciS I (prototype Sap I)**Isolated from *Planococcus citreus* S**GCTCTTC(N)₁[^]
CGAGAAG(N)₄[^]****E497
E498****50 u.a.
250 u.a.****Concentration:** 500 - 2 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**10 \times SE Buffer B**Reaction conditions:** 1 \times SE Buffer B

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

50 mM KCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 μ g/ml BSA;

50% glycerol. Store at -20°C.

Diluent: SE Buffer A**Heat Inactivation:****Yes** (65° C for 20 minutes)**Ligation/recutting assay:** After 3-fold overdigestion with PciS I, 90% of the DNA fragments can be ligated with T4 DNA Ligase at 16°C and 95% of these can be recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 1 units of PciS I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-----|-------|------|------|--------|
| Activity in SE Buffers (% of max) | 100 | 50-75 | 0-10 | 0-10 | 75-100 |

Pcs I (prototype Pcs I)Isolated from *Paracoccus carotinifaciens* 3K**(5mC)GNNNNN[^]NN(5mC)G
G(5mC)NN[^]NNNNNG(5mC)****E505
E506****50 u.a.
250 u.a.****See page 57 for more information about this enzyme.****Pct I (prototype Bsm I)**Isolated from *Planococcus citreus* SM**GAATGCN[^]
CTTAC[^]GN****E045
E046****1000 u.a.
5000 u.a.****Concentration:** 10 000 - 40 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**10 \times SE Buffer O**Reaction conditions:** 1 \times SE Buffer O

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

250 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 100 μ g/ml BSA;

50% glycerol. Store at -20°C.

Diluent: SE Buffer A**Heat Inactivation:****Yes** (65° C for 20 minutes)**Ligation/recutting assay:** After 20-fold overdigestion with Pct I, 90% of the DNA fragments can be ligated and recut.**Overdigestion assay:** No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 20units of Pct I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|-----|--------|-------|
| Activity in SE Buffers (% of max) | 25-50 | 50-75 | 100 | 75-100 | 10-25 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml

New product



Mammalian Genome Qualified



New package

Pkr I (prototype Pkr I)Isolated from *Planomicrobium koreense* 78k

DNA sequence with at least

three 5mC:

G(5mC)N^G(5mC)

(5mC)G^N(5mC)G

E579

E580

50 u.a.

250 u.a.

See page 58 for more information about this enzyme.

Ple19 I (prototype Pvu I)Isolated from *Pseudomonas lemoignei* 19

CGAT^CG

GC^TAGC

E195

E196

100 u.a.

500 u.a.

Concentration: 2 000 - 5 000 units/ml**Assayed on** λ DNA (Hind III-digest)**Reagents Supplied with Enzyme:**

10×SE Buffer Y

Reaction conditions: 1×SE Buffer Y

Incubate at 37°C.

Not blocked by overlapping*dam*-methylation (G^mATC): **CGATCG**.**Storage buffer:** 10 mM Tris-HCl (pH 7.5);

100 mM KCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 μ g/ml BSA;

50% glycerol. Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:****Yes** (65° C for 20 minutes)**Ligation/recutting assay:** After 3-fold

overdigestion with Ple19 I, 90% of the

DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific

activity was detected after incubation of

1 μ g of λ DNA with 5 units of Ple19 I

for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|--------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 25-50 | 25-50 | 100 |

Pps I (prototype Ple I)Isolated from *Pseudomonas pseudoalcaligenes*GAGTC(N)₄^CTCAG(N)₅^

E269

E270

25 u.a.

125 u.a.

Concentration: 500 - 1 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer Y, BSA

Reaction conditions: 1×SE Buffer Y+BSA

Incubate at 37°C.

*To obtain 100% activity, BSA should be**added to the 1× reaction mix to a final**concentration of 100 μ g/ml.***Storage buffer:** 10 mM Tris-HCl (pH 7.5);

50 mM KCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 μ g/ml BSA ,

50% glycerol. Store at -20°C.

Diluent: SE Buffer A**Heat Inactivation:****Yes** (65° C for 20 minutes)**Ligation/recutting assay:** After 2-fold

overdigestion with Pps I, 20% of the

DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific

activity was detected after incubation of

1 μ g of λ DNA with 2 unit of Pps I for

16 hours.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|------|-------|-----|
| Activity in SE Buffers (% of max) | 50-75 | 10-25 | 0-10 | 25-50 | 100 |

Psi I (prototype Psi I)Isolated from *Pseudomonas* species SE-G49

TTA^TAA

AAT^ATT

E279

E280

200 u.a.

1000 u.a.

Concentration: 5 000 - 10 000 units/ml**Assayed on** λ DNA**Reagents Supplied with Enzyme:**

10×SE Buffer B

Reaction conditions: 1×SE Buffer B

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

200 mM KCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 μ g/ml BSA;

50% glycerol. Store at -20°C.

Diluent: SE Buffer A**Heat Inactivation:****Yes** (65° C for 20 minutes)**Ligation/recutting assay:** After 5-fold

overdigestion with Psi I, about 50% of

the DNA fragments can be ligated. Of

these, 95% can be recut.

Overdigestion assay: No nonspecific

activity was detected after incubation of

1 μ g of λ DNA with 20 units of Psi I for

16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-----|-------|-------|-------|--------|
| Activity in SE Buffers (% of max) | 100 | 25-50 | 10-25 | 25-50 | 75-100 |

Psp124B I (prototype Sac I)Isolated from *Pseudomonas* species 124B

GAGCT^C

C^TCGAG

E107

E108

1000 u.a.

5000 u.a.

Concentration: 10 000 - 30 000 units/ml**Assayed on** λ DNA (Hind III-digest)**Reagents Supplied with Enzyme:**

10×SE Buffer G

Reaction mixture: 1×SE Buffer G

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

250 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 100 μ g/ml BSA;

50% glycerol. Store at -20°C.

Diluent: SE Buffer A**Heat inactivation:****No** (65° C for 20 minutes)**Yes** (80° C for 20 minutes)**Ligation/recutting assay:** After 20-fold

overdigestion with Psp124B I, > 90% of

the DNA fragments can be ligated and

recut.

Overdigestion assay: No nonspecific

activity was detected after incubation of

1 μ g of λ DNA with 20 units of Psp124B

I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|-----|-------|------|--------|
| Activity in SE Buffers (% of max) | 75-100 | 100 | 10-25 | 0-10 | 75-100 |

Psp6 I (prototype EcoR II)Isolated from *Pseudomonas* species 6

^CCWGG

GGWCC^

E453

E454

100 u.a.

500 u.a.

Concentration: 1 000 - 3 000 units/ml**Assayed on** λ DNA (*dcm*-)**Reagents Supplied with Enzyme:**

10×SE Buffer B

Reaction conditions: 1×SE Buffer B

Incubate at 55°C.

Blocked by overlapping *dcm*-methylation(C^mCWGG): **CCWGG**.**Storage buffer:** 10 mM Tris-HCl (pH 7.5);

100 mM KCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 μ g/ml BSA;

50% glycerol. Store at -20°C.

Diluent: SE Buffer A**Heat Inactivation:****No** (65° C for 20 minutes)**Yes** (80° C for 20 minutes)**Ligation/recutting assay:** After 3-fold

overdigestion with Psp6 I, 95% of the

DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific

activity was detected after incubation of

1 μ g of λ DNA with 6 units of Psp6 I for

16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-----|-------|-------|-------|--------|
| Activity in SE Buffers (% of max) | 100 | 50-75 | 10-25 | 25-50 | 75-100 |



PspC I (prototype PmaC I)

Isolated from *Pseudomonas* species C

Concentration: 10 000 - 30 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer B, BSA

Reaction conditions: 1 \times SE Buffer B+BSA

Incubate at 37°C.

To obtain 100% activity, BSA should be

added to the 1 \times reaction mix to a final

concentration of 100 μ g/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

50 mM KCl; 0.1 mM EDTA;

10 mM 2-mercaptoethanol; 200 μ g/ml BSA;

50% glycerol. Store at -20°C (*see note).

Diluent: SE Buffer A

Heat Inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold

overdigestion with PspC I, > 90% of the

DNA fragments can be ligated and recut.

In the presence of 10% PEG ligation is

better.

Overdigestion assay: No nonspecific

activity was detected after incubation of

1 μ g of λ DNA with 40 units of PspC I

for 16 hours.

Do not use BSA for long incubation.

***Note: For long term storage (more**

than 30 days), store at -70°C .

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-----|-------|---|---|-------|
| Activity in SE Buffers (% of max) | 100 | 50-75 | 0 | 0 | 50-75 |

PspE I (prototype BstE II)

Isolated from *Pseudomonas* species E

Concentration: 5 000 - 10 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer B

Reaction conditions: 1 \times SE Buffer B

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

50 mM KCl; 0,1 mM EDTA;

7 mM 2-mercaptoethanol; 200 μ g/ml BSA;

50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold

overdigestion with PspE I

>90% of the DNA fragments can be

ligated and recut.

Overdigestion assay: No nonspecific

activity was detected after incubation of

1 μ g of λ DNA with 10 units of PspE I

for 16 hours.

Star activity: High enzyme

concentration may result in star activity.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-----|-------|-------|-------|-------|
| Activity in SE Buffers (% of max) | 100 | 50-75 | 25-50 | 50-75 | 50-75 |

PspL I (prototype Spl I)

Isolated from *Pseudomonas* species L

Concentration: 2 000 - 5 000 units/ml

Assayed on λ DNA (HindIII-digest)

Reagents Supplied with Enzyme:

10 \times SE Buffer Y, BSA

Reaction conditions: 1 \times SE Buffer Y+BSA

Incubate at 37°C.

To obtain 100% activity, BSA should be

added to the 1 \times reaction mix to a final

concentration of 100 μ g/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

100 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 100 μ g/ml BSA;

50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 2-fold

overdigestion with PspL I, 95% of the

DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific

activity was detected after incubation of

1 μ g of λ DNA with 2 units of PspL I for

16 hours.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|--------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 25-50 | 10-25 | 100 |

PspN4 I (prototype Nla IV)

Isolated from *Pseudomonas* species N4

Concentration: 10 000 - 20 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer Y

Reaction conditions: 1 \times SE Buffer Y

Incubate at 37°C.

Blocked by methylation

5'-GGNN(5mC)C-3'/3'-C(5mC)NNGG-5' or.

5'-GGNN(5mC)C-3'/3'-CCNNGG-5'

Storage buffer: 10 mM Tris-HCl (pH 7.5);

50 mM KCl; 0,1 mM EDTA;

7 mM 2-mercaptoethanol; 200 μ g/ml BSA;

50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold

overdigestion with PspN4 I, 95% of the

DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific

activity was detected after incubation of

1 μ g of λ DNA with 30 units of PspN4 I

for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 10-25 | 10-25 | 10-25 | 25-50 | 100 |

PspOM I (prototype Apa I)

Isolated from an *E.coli* strain that carries the

cloned PspOM I gene from *Pseudomonas* species

OM2164

Concentration: 10 000 units/ml

Assayed on λ DNA (Bam HI-digest)

Reagents Supplied with Enzyme:

10 \times SE Buffer Y

Reaction conditions: 1 \times SE Buffer Y

Incubate at 37°C

PspOMI is a neoschizomer of ApaI.

Storage buffer: 20 mM Tris-HCl (pH 7.5);

50 mM NaCl; 0,1 mM EDTA; 1 mM DTT;

200 μ g/ml BSA; 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat Inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold

overdigestion with PspOM I, 95% of the

DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific

activity was detected after incubation of

1 μ g of λ DNA with 20 units of PspOM

I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|-------|------|------|-----|
| Activity in SE Buffers (% of max) | 75-100 | 10-25 | 0-10 | 0-10 | 100 |



New product



Mammalian Genome Qualified



New package

PspPP I (prototype PpuM I)

Isolated from *Pseudomonas* species PP

Concentration: 2 000 - 5 000 units/ml
Assayed on λ DNA (dcm-,HindIII-digest)
Reagents Supplied with Enzyme:
10 \times SE Buffer Y, BSA
Reaction conditions: 1 \times SE Buffer Y+BSA
Incubate at 37°C.
Blocked by overlapping *dcm*-methylation
(C^mCWGG): **RGWCCTGG**.
To obtain 100% activity, BSA should be added to the 1 \times reaction mix to a final concentration of 100 μ g/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
200 mM NaCl; 0.1 mM EDTA;
7 mM 2-mercaptoethanol; 200 μ g/ml BSA;
50% glycerol. Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
Yes (65° C for 20 minutes)

Ligation/recutting assay: After 2-fold overdigestion with PspPP I, >70% of the DNA fragments can be ligated and 80% of these can be recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 5 units of PspPP I for 16 hours.
Do not use BSA for long incubation.

| | | | | | |
|-----------------------------------|-------|-------|---|------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 10-25 | 10-25 | 0 | 0-10 | 100 |

Psp X I (prototype PspX I)

Isolated from an *E.coli* strain that carries the cloned PspX I gene from *Pseudomonas* species A1-1

Concentration: 10 000 and 50 000 units/ml
Assayed on λ DNA (HindIII-digest)
Reagents Supplied with Enzyme:
10 \times SE Buffer Y, BSA
Reaction conditions: 1 \times SE Buffer Y+BSA.
Incubate at 37°C.
To obtain 100% activity, BSA should be added to the 1 \times reaction mix to a final concentration of 100 μ g/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
200 mM KCl; 0.1 mM EDTA;
7 mM 2-mercaptoethanol, 200 μ g/ml BSA;
50% glycerol. Store at -20°C.
Diluent: SE Buffer A
Heat Inactivation:
No (65° C for 20 minutes)
Yes (80° C for 20 minutes)

For high concentration
E478X **1000 u.a.**
Ligation/recutting assay: After 20-fold overdigestion with PspX I, > 90% of the DNA fragments can be ligated and recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 20 units of PspXI for 16 hours.
Do not use BSA for long incubation.

| | | | | | |
|-----------------------------------|-------|-------|-------|--------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 50-75 | 50-75 | 25-50 | 75-100 | 100 |

PsrI (prototype Psr I)

Isolated from *Pseudomonas stutzeri* N2

Concentration: 1 000 – 3 000units/ml
Assayed on T7 DNA
Reagents Supplied with Enzyme:
10 \times SE Buffer Y, BSA
Reaction mixture: 1 \times SE Buffer Y+BSA
Incubate at 30°C.
Incubation at 37°C results in 20% activity.
To obtain 100% activity, BSA should be added to the 1 \times reaction mix to a final concentration of 100 μ g/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
50 mM KCl; 0,1 mM EDTA;
7 mM 2-mercaptoethanol; 200 μ g/ml BSA,
50% glycerol. Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
Yes (65° C for 20 minutes)

Ligation/recutting assay: After 2-fold overdigestion with Psr I, > 70% of the DNA fragments can be ligated. Of these, 80% of these can be recut. In the presence of 10% PEG ligation is better.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of T7 DNA with 2 units of Psr I for 16 hours.
Star activity: High enzyme concentration may result in star activity.
Do not use BSA for long incubation.

| | | | | | |
|-----------------------------------|-------|-------|---|------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 10-25 | 10-25 | 0 | 0-10 | 100 |

Pst I (prototype Pst I)

Isolated from an *E.coli* strain that carries the cloned PstI gene from *Providencia stuartii*


Concentration: 20 000 and 50 000 units/ml
Assayed on λ DNA
Reagents Supplied with Enzyme:
10 \times SE Buffer O, BSA
Reaction conditions: 1 \times SE Buffer O+BSA
Incubate at 37°C.
To obtain 100% activity, BSA should be added to the 1 \times reaction mix to a final concentration of 100 μ g/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
200 mM NaCl; 0,1 mM EDTA;
7 mM 2-mercaptoethanol; 200 μ g/ml BSA;
50% glycerol.
Store at -20°C.
Diluent: SE Buffer A
Heat inactivation:
No (65° C for 20 minutes)

For high concentration
E109X **4000 u.a.**
E110X **20000 u.a.**
Ligation/recutting assay: After 20-fold overdigestion with Pst I, > 90% of the DNA fragments can be ligated and recut.
Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 40 units of Pst I for 16 hours.
Star activity: High enzyme concentration may result in star activity.
Do not use BSA for long incubation.

| | | | | | |
|-----------------------------------|-------|-------|-----|-------|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 10-25 | 25-50 | 100 | 25-50 | 25-50 |


| | | | | | |
|---|---|---|-------------------------------|-------|-----|
| PstN I (prototype AlwN I) Isolated from <i>Bacillus sphaericus</i> | CAGNNN[^]CTG GTC[^]NNNGAC | E571 E572 | 500 u.a. 2500 u.a. | | |
| Concentration: 5 000 – 10 000 units/ml Assayed on λ DNA Reagents Supplied with Enzyme: 10×SE Buffer Y Reaction conditions: 1×SE Buffer Y Incubate at 37°C. | Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat Inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 10-fold overdigestion with PstNI > 95% of Lambda DNA fragments can be ligated with T4 DNA Ligase and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of PstNI for 16 hours at 37°C. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 25-50 | 50-75 | 10-25 | 25-50 | 100 |

| | | | | | |
|---|---|--|---------------------------------|-------|-------|
| Pvu II (prototype Pvu II) Isolated from an <i>E.coli</i> strain that carries the cloned Pvu II gene from <i>Proteus vulgaris</i> |  CAG[^]CTG GTC[^]GAC | E111 E112 | 2000 u.a. 10000 u.a. | | |
| Concentration: 10 000 units/ml Assayed on λ DNA Reagents Supplied with enzyme: 10×SE Buffer G, BSA Reaction conditions: 1× SE Buffer G+BSA Incubate at 37°C. <i>To obtain 100% activity, BSA should be added to the 1×reaction mix to a final concentration of 100 µg/ml.</i> | Storage buffer: 10 mM Tris-HCl (pH 7.5); 300 mM NaCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 100 µg/ml BSA; 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: No (65° C for 20 minutes) | Ligation/recutting assay: After 10-fold overdigestion with Pvu II, 70% of the DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of Pvu II for 16 hours. Star activity: High enzyme concentration may result in star activity. Do not use BSA for long incubation. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 25-50 | 100 | 25-50 | 25-50 | 25-50 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml

| | | | | | |
|---|---|---|-------------------------------|-------|-----|
| Rga I (prototype Sgf I) Isolated from <i>Rhizobium galegae</i> | GCGAT[^]CGC CGC[^]TAGCG | E491 E492 | 200 u.a. 1000 u.a. | | |
| Concentration: 5 000 - 10 000 units/ml Assayed on Ad2 DNA Reagents Supplied with Enzyme: 10×SE Buffer Y Reaction conditions: 1×SE Buffer Y Incubate at 55°C. Not blocked by overlapping <i>dam</i> -methylation (G ^m ATC): GCGATCGC . | Storage buffer: 10 mM Tris-HCl (pH 7.5); 250 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat Inactivation: No (65° C for 20 minutes) Yes (80° C for 20 minutes) | Ligation/recutting assay: After 5-fold overdigestion with Rga I, > 90% of the DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of Ad2 DNA with 10 units of Rga I for 16 hours. Star activity: High enzyme concentration may result in star activity. | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 50-75 | 10-25 | 25-50 | 100 |

| | | | | | |
|---|---|---|------------------------------|-------|-------|
| Rig I (prototype Fse I) Isolated from <i>Rhizobium yangligense</i> | GGCCGG[^]CC CC[^]GGCCGG | E529 E530 | 100 u.a. 500 u.a. | | |
| Concentration: 1 000 - 5 000 units/ml Assayed on Adenovirus-2 DNA Reagents Supplied with Enzyme: 10×SE Buffer RigI, BSA Reaction conditions: 1×SE Buffer RigI+BSA Incubate at 37°C. Blocked by mCG or GmC methylation: 5'-GGC(m5C)GGCC-3'/3'-CCGG(m5C)CGG-5' or 5'-GG(m5C)CGG(m5C)C-3'/ 3'-C(m5C)GGC(m5C)GG-5' <i>To obtain 100% activity, BSA should be added to the 1×reaction mix to a final concentration of 100 µg/ml.</i> | Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM KCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA, 50% glycerol. Store at -20°C(*see note) Diluent: SE Buffer A Heat Inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 5-fold overdigestion with Rig I, > 95% of Ad2 DNA fragments can be ligated with T4 DNA Ligase and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of Ad2 DNA with 6 units of Rig I for 16 hours. Do not use BSA for long incubation. *Note: For long term storage (more than 7 days), store at -70°C . | | | |
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 50-75 | 0-10 | 10-25 | 50-75 |

| | | | | |
|---|--|---|----------------------|--------------------------------|
| Rsa I (prototype Rsa I) Isolated from <i>Rhodospseudomonas sphaeroides</i> |  | GT[^]AC CA[^]TG | E113 E114 | 1000 u.a. 5000 u.a. |
| Concentration: 10 000 – 30 000 units/ml Assayed on λ DNA Reagents Supplied with Enzyme: 10×SE Buffer B Reaction conditions: 1× SE Buffer B Incubate at 37°C. | Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 100 µg/ml BSA, 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: No (80° C for 20 minutes) | Ligation/recutting assay: After 20-fold overdigestion with Rsa I, > 90% of the DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 40 units of Rsa I for 16 hours. | | |

| | | | | | |
|-----------------------------------|-----|-------|------|-------|--------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 100 | 50-75 | 0-10 | 50-75 | 75-100 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml

| | | | |
|---|---|---|-------------------------------|
| RsaN I (prototype Rsa I) Isolated from <i>Rhodospseudomonas sphaeroides</i> N | G[^]TAC CAT[^]G | E555 E556 | 200 u.a. 1000 u.a. |
| Concentration: 5 000 units/ml Assayed on λ DNA Reagents Supplied with Enzyme: 10×SE Buffer B Reaction conditions: 1× SE Buffer B Incubate at 37°C. | Storage buffer: 10 mM Tris-HCl (pH 7.5); 250 mM NaCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 100 µg/ml BSA, 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: Yes (80° C for 20 minutes) | Ligation/recutting assay: After 5-fold overdigestion with RsaN I, > 90% of λ DNA fragments can be ligated with T4 DNA Ligase and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of RsaN I for 16 hours. | |

| | | | | | |
|-----------------------------------|-----|--------|-------|-------|--------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 100 | 75-100 | 50-75 | 50-75 | 75-100 |

| | | | |
|--|---|---|--------------------------------|
| Rsr2 I (prototype Rsr II) Isolated from <i>Rhodobacter sphaeroides</i> 12 | CG[^]GWCCG GCCWG[^]GC | E281 E282 | 1000 u.a. 5000 u.a. |
| Concentration: 10 000 - 30 000 units/ml Assayed on λ DNA Reagents Supplied with Enzyme: 10×SE Buffer Y, BSA Reaction conditions: 1×SE Buffer Y+BSA Incubate at 37°C. <i>To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.</i> | Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 100 µg/ml BSA, 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat Inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 20-fold overdigestion with Rsr2 I, > 90% of the DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 60 units of Rsr2 I for 16 hours. Do not use BSA for long incubation. | |

| | | | | | |
|-----------------------------------|-------|--------|------|-------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 50-75 | 75-100 | 0-10 | 10-25 | 100 |

| | | | |
|--|--|--|---------------------------------|
| Sal I (prototype Sal I) Isolated from an <i>E.coli</i> strain that carries the cloned Sal I gene from <i>Streptomyces albus</i> | G[^]TCGAC CAGCT[^]G | E115 E116 | 2000 u.a. 10000 u.a. |
| Concentration: 10 000 units/ml Assayed on λ DNA (Hind III-digest) Reagents Supplied with Enzyme: 10×SE Buffer O Reaction conditions: 1× SE Buffer O Incubate at 37°C. | Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM NaCl; 0,1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat inactivation: Yes (65° C for 20 minutes) | Ligation/recutting assay: After 10-fold overdigestion with Sal I, 95% of the DNA fragments can be ligated and recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of Sal I for 16 hours. Star activity: High enzyme concentration may result in star activity. | |

| | | | | | |
|-----------------------------------|------|-------|-----|-------|------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 0-10 | 10-25 | 100 | 25-50 | 0-10 |

| | | | |
|--|---|---|-------------------------------|
| Sbf I (prototype Sse8387 I) Isolated from <i>Streptomyces</i> species Bf61 | CCTGCA[^]GG GG[^]ACGTCC | E101 E102 | 200 u.a. 1000 u.a. |
| Concentration: 5 000 units/ml Assayed on λ DNA Reagents Supplied with Enzyme: 10×SE Buffer Y Reaction conditions: 1×SE Buffer Y Incubate at 37°C. | Storage buffer: 10 mM Tris-HCl (pH 7.6); 50 mM NaCl; 0.1 mM EDTA; 1 mM DTT; 200 µg/ml BSA; 50% glycerol. Store at -20°C. Diluent: SE Buffer A Heat Inactivation: No (65° C for 20 minutes) Yes (80° C for 20 minutes) | Ligation/recutting assay: After 5-fold overdigestion with Sbf I, > 90% of the DNA fragments can be ligated. Of these, 90% can be recut. Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of Sbf I for 16 hours. Star activity: High enzyme concentration may result in star activity. | |

| | | | | | |
|-----------------------------------|--------|-------|------|------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 50-75 | 0-10 | 0-10 | 100 |

Set I (prototype Set I)

Isolated from an *E.coli* strain that carries the cloned Set I gene from *Streptomyces werraensis*
37

ASST[^] E537 200 u.a.
^TSSA E538 1000 u.a.

For high concentration
E538X 1000 u.a.

Set I is a restriction endonuclease with 4 expected recognition sites as well as relaxed non-cognate sites (star sites).

Concentration: 5 000 and 20 000 units/ml
One unit is defined as the amount of enzyme required to cleave 1 pmol of the double-stranded oligonucleotide with the following structure

5'-CGAGTTATAGCTGGGCCCAAC-3'

3'-GCTCAAATATCGACCCGGGTTG-5'

in 1 hour at 50°C in a total reaction volume of 20 µl.

Reagents Supplied with Enzyme:

10×SE Buffer Y

Reaction conditions: 1×SE Buffer Y

Incubate at 55°C.

Storage buffer: 10 mM Tris-HCl (pH 7.6); 100 mM NaCl; 0.1 mM EDTA; 1mM DTT, 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat Inactivation:

No (65° C for 20 minutes)

Yes (80° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with Set I, 50% of the pBR322 DNA fragments can be ligated with T4 DNA Ligase and recut.

Note! In the case of long incubation with Set I DNA can be digested to small oligos.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|--------|--------|-----|
| Activity in SE Buffers (% of max) | 25-50 | 25-50 | 75-100 | 75-100 | 100 |

SfaN I (prototype SfaN I)

Isolated from an *E.coli* strain that carries the cloned SfaN I gene from *Streptococcus faecalis* N



GCATC(N)₅[^] E165 500 u.a.
CGTAG(N)₉[^] E166 2500 u.a.

Concentration: 10 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer O

Reaction conditions: 1× SE Buffer O

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 300 mM NaCl; 0,1 mM EDTA; 1 mM DTT; 200 µg/ml BSA; 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

No (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with SfaN I, > 95% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of SfaN I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|-----|--------|------|
| Activity in SE Buffers (% of max) | 10-25 | 25-50 | 100 | 75-100 | 0-10 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml

Sfi I (prototype Sfi I)

Isolated from an *E.coli* strain that carries the cloned Sfi I gene from *Streptomyces fimbriatus*

GGCCNNNN[^]NGGCC E123 1000 u.a.
CCGGN[^]NNNNCCGG E124 5000 u.a.

For high concentration
E123X 1000 u.a.
E124X 5000 u.a.

Concentration: 10 000 and 40 000 units/ml
Assayed on T7 DNA

Reagents Supplied with Enzyme:

10×SE Buffer G, BSA

Reaction conditions: 1×SE Buffer G+BSA

Incubate at 50°C.

Blocked by overlapping *dcm*-methylation (C^mCWGG): GGCCWGGNNGGCC.

Not blocked by overlapping

dcm-methylation

(C^mCWGG): GGCCNNNNNGGCCWGG.

To obtain 100% activity, BSA should be

added to the 1× reaction mix to a final

concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 200 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA, 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Sfi I, about 70% of the DNA fragments can be ligated and recut. In the presence of 10% PEG ligation is better.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of T7 DNA with 20 units of Sfi I for 16 hours.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|-----|-------|-------|-------|
| Activity in SE Buffers (% of max) | 75-100 | 100 | 25-50 | 25-50 | 25-50 |

Sfr274 I (prototype Xho I)

Isolated from *Streptomyces fradiae* 274

C[^]TCGAG E125 2000 u.a.
GAGCT[^]C E126 10000 u.a.

Concentration: 10 000 – 30 000 units/ml

Assayed on λ DNA (Hind III-digest)

Reagents Supplied with Enzyme:

10×SE Buffer B

Reaction conditions: 1×SE Buffer B

Incubate at 50°C.

Note: At 37°C activity is about 70% from maximum.

Blocked by CTCG^mAG methylation.

Not blocked by CT^mCGAG methylation.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 20-fold overdigestion with Sfr274 I, 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of Sfr274 I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-----|--------|-------|-------|--------|
| Activity in SE Buffers (% of max) | 100 | 75-100 | 50-75 | 50-75 | 75-100 |



New product



Mammalian Genome Qualified



New package

Sfr303 I (prototype Sac II)

Isolated from *Streptomyces fradiae* 303

CCGC[^]GG

GG[^]CGCC

E127

E128

1000 u.a.

5000 u.a.

Concentration: 5 000 – 20 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer B

Reaction conditions: 1×SE Buffer B

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

100 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 µg/ml BSA;

50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Sfr303 I,

> 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of Sfr303 I for 16 hours.

| | | | | | |
|-----------------------------------|-----|-------|-------|-------|--------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 100 | 50-75 | 10-25 | 10-25 | 75-100 |

Sma I (prototype Sma I)

Isolated from an *E.coli* strain that carries the cloned Sma I gene from *Serratia marcescens*

CCC[^]GGG

GGG[^]CCC

E177

E178

2000 u.a.

10000 u.a.

Concentration: 20 000 units/ml

Assayed on λ DNA (Hind III-digest)

Reagents Supplied with Enzyme:

10×SE Buffer Y

Reaction conditions: 1× SE Buffer Y

Incubate at 25°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

50 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 µg/ml BSA,

50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 20-fold overdigestion with Sma I,

> 90% of the DNA fragments can be ligated by high concentration T4 DNA Ligase and 10% PEG, and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 40 units of Sma I for 16 hours.

| | | | | | |
|-----------------------------------|------|------|------|------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 0-10 | 0-10 | 0-10 | 0-10 | 100 |

Smi I (prototype Swa I)

Isolated from *Streptococcus milleri* S

ATTT[^]AAAT

TAAA[^]TTTA

E225

E226

1000 u.a.

5000 u.a.

Concentration: 10 000 - 30 000 units/ml

Assayed on T7 DNA (Ssp I-digest)

Reagents Supplied with Enzyme:

10×SE Buffer O, BSA

Reaction conditions: 1×SE Buffer O+BSA

Incubate at 37°C.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

250 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 100 µg/ml BSA;

50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with Smi I, about 80% of the DNA fragments can be ligated and recut. In the presence of 10% PEG ligation is better.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of T7 DNA with 30 units of Smi I for 16 hours.

Do not use BSA for long incubation.

| | | | | | |
|-----------------------------------|-------|-------|-----|--------|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 25-50 | 25-50 | 100 | 75-100 | 25-50 |

SmiM I (prototype Msl I)

Isolated from *Spingobacterium mizutae* M

CAYNN[^]NNRTG

GTRNN[^]NNYAC

E293

E294

500 u.a.

2500 u.a.

Concentration: 5 000 – 10 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer W

Reaction conditions: 1×SE Buffer W

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

250 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 µg/ml BSA;

and 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat Inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with SmiM I, 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of DNA with 20 units of SmiM I for 16 hours.

| | | | | | |
|-----------------------------------|-------|-------|--------|-----|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 10-25 | 10-25 | 75-100 | 100 | 10-25 |

Sph I (prototype Sph I)

Isolated from an *E.coli* strain that carries the cloned SphI gene from *Streptomyces phaeochromogenes*

GCATG[^]C

C[^]GTACG

E129

E130

500 u.a.

2500 u.a.

Concentration: 5 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer G, BSA

Reaction conditions: 1×SE Buffer G+BSA

Incubate at 37°C.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

100 mM NaCl; 0.1 mM EDTA; 1 mM DTT;

100 µg/ml BSA; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Sph I, 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of Sph I for 16 hours.

Do not use BSA for long incubation.

| | | | | | |
|-----------------------------------|-------|-----|--------|--------|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 25-50 | 100 | 75-100 | 75-100 | 50-75 |

Sse9 I (prototype Tsp509 I)

Isolated from an *E.coli* strain that carries the cloned Sse9I gene from *Sporosarcina* species 9



[^]AATT
TTAA[^]

E217
E218

500 u.a.
2500 u.a.

Concentration: 5 000 units/ml

Assayed on pBR322 DNA

Reagents Supplied with Enzyme:

10×SE Buffer B, BSA

Reaction conditions: 1×SE Buffer B+BSA

Incubate at 55°C.

Incubation at 37°C results in 75% activity.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM NaCl; 0.1 mM EDTA; 1 mM DTT; 200 µg/ml BSA; 50% glycerol. Store at -20° C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with Sse9 I, > 95% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 5 units of Sse9 I for 16 hours.

Do not use BSA for long incubation.

| | | | | | |
|-----------------------------------|-----|--------|-------|-------|--------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 100 | 75-100 | 50-75 | 50-75 | 75-100 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml

Ssp I (prototype Ssp I)

Isolated from an *E.coli* strain that carries the cloned SspI gene from *Sphaerotilus* species

AAT[^]ATT
TTA[^]TAA

E041
E042

500 u.a.
2500 u.a.

Concentration: 10 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer K, BSA

Reaction conditions: 1×SE Buffer K+BSA

Incubate at 37°C.

Blocked by A^mATATT methylation.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 200 mM NaCl; 0,1 mM EDTA; 1 mM DTT; 200 µg/ml BSA; 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with SspI 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of Ssp I for 16 hours.

Star activity: High enzyme concentration may result in star activity.

Do not use BSA for long incubation.

| | | | | | |
|-----------------------------------|--------|-------|-------|-------|--------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 50-75 | 25-50 | 50-75 | 75-100 |

SspM I (prototype Mae I)

Isolated from *Sporosarcina* species M



C[^]TAG
GAT[^]C

E591
E592

100 u.a.
500 u.a.

Concentration: 1 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer Y

Reaction conditions: 1×SE Buffer Y

Incubate at 55°C.

At 37°C activity is 75% from maximum.

Storage buffer: 20 mM Tris-HCl (pH 7.6); 100 mM KCl; 0,1 mM EDTA; 0,01% Triton X-100; 500 µg/ml BSA; 7 mM 2-mercaptoethanol; 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

No (80° C for 20 minutes)

Ligation/recutting assay: After 3-fold overdigestion with SspMI 5% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 2 units of SspM I for 16 hours.

***Note: For long term storage (more than 30 days), store at -70°C ..**

| | | | | | |
|-----------------------------------|-------|-------|-------|-------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 50-75 | 25-50 | 10-25 | 50-75 | 100 |

Taq I (prototype Taq I)

Isolated from an *E.coli* strain that carries the cloned Taq I gene from *Thermus aquaticus*



T[^]CGA
AGC[^]T

E133
E134

2000 u.a.
10000 u.a.

Concentration: 20 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer Y, BSA

Reaction conditions: 1×SE Buffer Y+BSA

Incubate at 65°C.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 300 mM NaCl; 0,1 mM EDTA; 1 mM DTT ; 200 µg/ml BSA; 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

No (65° C for 20 minutes)

Ligation/recutting assay: After 20-fold overdigestion with Taq I, > 95% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of Taq I for 16 hours.

Do not use BSA for long incubation.

| | | | | | |
|-----------------------------------|-------|--------|--------|-------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 50-75 | 75-100 | 75-100 | 50-75 | 100 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml



New product



Mammalian Genome Qualified



New package

Tru9 I (prototype Mse I)

Isolated from an *E.coli* strain that carries the cloned Tru9 I gene from *Thermus ruber* 9



T[^]TAA
AAT[^]T

E199
E200

500 u.a.
2500 u.a.

Concentration: 20 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer W

Reaction conditions: 1 \times SE Buffer W

Incubate at 65°C.

Blocked by TTA^mA methylation.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
50 mM NaCl; 0,1 mM EDTA; 1 mM DTT;
200 μ g/ml BSA; 50% glycerol.
Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (80° C for 20 minutes)

Ligation/recutting assay: After 20-fold overdigestion with Tru9 I,
> 95% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 40 units of Tru9 I for 16 hours.

| | | | | | |
|-----------------------------------|--------|-------|-------|-----|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 25-50 | 25-50 | 100 | 50-75 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml

TseF I (prototype Tsp45 I)

Isolated from an *E.coli* strain that carries the cloned TseF I gene from *Thermus species* F35



[^]GTSAC
CASTG[^]

E589
E590

200 e.a.
1000 e.a.

Concentration: 5 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer B

Reaction conditions: 1 \times SE Buffer B

Incubate at 65°C.

Storage buffer: 10 mM Tris-HCl (pH 7.6);
50 mM NaCl; 0,1 mM EDTA; 1 mM DTT;
200 μ g/ml BSA; 50% glycerol.
Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

No (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with TseF I,
about 95% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 5 units of TseFI for 16 hours.

| | | | | | |
|-----------------------------------|-----|-------|------|-------|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 100 | 50-75 | 0-10 | 25-50 | 50-75 |

Tth111 I (prototype Tth111 I)

Isolated from an *E.coli* strain that carries the cloned Tth111 I gene from *Thermus thermophilus* 111

GACN[^]NNGTC
CTGNN[^]NCAG

E097
E098

400 u.a.
2000 u.a.

Concentration: 5 000 units/ml

Assayed on λ DNA (HindIII-digest)

Reagents Supplied with Enzyme:

10 \times SE Buffer Y

Reaction conditions: 1 \times SE Buffer Y

Incubate at 65°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
500 mM NaCl; 0,1 mM EDTA; 1 mM DTT ;
50% glycerol.
Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

No (65° C for 20 minutes)

Ligation/recutting assay: After 2-fold overdigestion with Tth111 I, about 10% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 5 units of Tth111 I for 16 hours.

Star activity: Long incubation or conditions of high enzyme concentration may result in star activity. Star activity is observed at a greater than 5-fold overdigestion of 1 μ g substrate with Tth1111 for 1 hour.

| | | | | | |
|-----------------------------------|--------|-------|-------|-------|-----|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 75-100 | 50-75 | 10-25 | 10-25 | 100 |

Vne I (prototype ApaL I)

Isolated from an *E.coli* strain that carries the cloned Vnel gene from *Vibrio nereis* 18

G[^]TGCAC
CACGT[^]G

E137
E138

1000 u.a.
5000 u.a.

Concentration: 10 000 – 20 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer O

Reaction conditions: 1 \times SE Buffer O

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0,1 mM EDTA;
7 mM 2-mercaptoethanol; 200 μ g/ml BSA;
50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Vne I, 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 20 units of Vne I for 16 hours.

| | | | | | |
|-----------------------------------|-------|-------|-----|-------|-------|
| SE Buffers | B | G | O | W | Y |
| Activity in SE Buffers (% of max) | 10-25 | 25-50 | 100 | 25-50 | 25-50 |

Vsp I (prototype Vsp I)

Isolated from an *E.coli* strain that carries the cloned Vsp I gene from *Vibrio* species 343



AT[^]TAAT
TAAT[^]TA

E139
E140

1000 u.a.
5000 u.a.

Concentration: 10 000 units/ml

Assayed on λ DNA

Reagents Supplied with Buffer:

10×SE Buffer W

Reaction conditions: 1×SE Buffer W

Incubate at 37°C.

Blocked by ATTA^mAT methylation.

Storage buffer: 10 mM Tris HCl (pH 7.6);
50 mM NaCl; 0.1 mM EDTA; 1 mM DTT;
200 µg/ml BSA; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Vsp I, 70% of the DNA fragments can be ligated. Of these, 90% can be recut. In the presence of 10% PEG ligation is better.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 20 units of Vsp I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|------|-------|-------|-----|-------|
| Activity in SE Buffers (% of max) | 0-10 | 10-25 | 50-75 | 100 | 25-50 |

Qualified for mammalian genome studies: http://science.sibenzyme.com/article8_article_28_1.phtml

Xba I (prototype Xba I)

Isolated from an *E.coli* strain that carries the cloned Xba I gene from *Xanthomonas badrii*



T[^]CTAGA
AGATC[^]T

E141
E142

2000 u.a.
10000 u.a.

For high concentration

E141X

2000 u.a.

E142X

10000 u.a.

Concentration:

20 000 and 50 000 units/ml

Assayed on λ DNA (dam⁻/ Hind III-digest)

Reagents Supplied with Enzyme:

10×SE Buffer O, BSA

Reaction conditions: 1×SE Buffer O+BSA

Incubate at 37°C.

Blocked by overlapping dam-methylation

(G^mATC): TCTAGATC.

To obtain 100% activity, BSA should be added to the 1× reaction mix to a final concentration of 100 µg/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5);
50 mM NaCl; 0.1 mM EDTA; 1 mM DTT;-
200 µg/ml BSA; 50% glycerol.

Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 20-fold overdigestion with Xba I, > 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of DNA with 40 units of Xba I for 16 hours.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|--------|-----|-------|--------|
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 100 | 50-75 | 75-100 |

Qualified for human genome studies: http://science.sibenzyme.com/article8_article_31_1.phtml

Xma I (prototype Xma I)

Isolated from an *E.coli* strain that carries the cloned Xma I gene from *Xanthomonas malvacearum*

C[^]CCGGG
GGGCC[^]C

E233
E234

300 u.a.
1500 u.a.

Concentration: 3000 units/ml

Assayed on Adenovirus-2 DNA

Reagents Supplied with Enzyme:

10×SE Buffer Y

Reaction conditions: 1×SE Buffer Y

Incubate at 37°C.

Storage buffer: 20 mM Tris-HCl (pH 7.6);
50 mM NaCl; 0.1 mM EDTA; 1 mM DTT;
200 µg/ml BSA; 50% glycerol. Store at -20°
C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 3-fold overdigestion with XmaI, 95% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 3 units of Xma I for 16 hours.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|-------|---|------|-----|
| Activity in SE Buffers (% of max) | 75-100 | 50-75 | 0 | 0-10 | 100 |

Zra I (prototype Aat II)

Isolated from an *E.coli* strain that carries the cloned Zra I gene from *Zoogloea ramigera*11

GAC[^]GTC
CTG[^]CAG

E463
E464

200 u.a.
1000 u.a.

Concentration: 10 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10×SE Buffer B

Reaction conditions: 1×SE Buffer B

Incubate at 37°C.

ZraI is a neoschizomer of AatII.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 50
mM NaCl; 0.1 mM EDTA; 1 mM DTT;
200 µg/ml BSA; 50% glycerol. Store at -20°
C.

Diluent: SE Buffer A

Heat Inactivation:

No (65° C for 20 minutes)

Yes (80° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Zra I, about 90% of the DNA fragments can be ligated and recut. In the presence of 10% PEG ligation is better.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 10 units of Zra I for 16 hours.

Star activity: High enzyme concentration results in star activity.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-----|-------|-------|-------|--------|
| Activity in SE Buffers (% of max) | 100 | 50-75 | 25-50 | 25-50 | 75-100 |



New product



Mammalian Genome Qualified



New package

Zrm I (prototype Sca I)

Isolated from *Zoogloea ramigera* SCA

AGT[^]ACT
TCA[^]TGA

E005
E006

1000 u.a.
5000 u.a.

Concentration: 10 000 – 20 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer Y, BSA

Reaction conditions: 1 \times SE Buffer Y+BSA

Incubate at 37°C.

To obtain 100% activity, BSA should be

added to the 1 \times reaction mix to a final

concentration of 100 μ g/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0,1 mM EDTA;

7 mM 2-mercaptoethanol; 200 μ g/ml BSA; 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Zrm I, 70% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 20 units of Zrm I for 16 hours.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|------|------|-----|
| Activity in SE Buffers (% of max) | 50-75 | 25-50 | 0-10 | 0-10 | 100 |

Zsp2 I (prototype Ava III)

Isolated from *Zoogloea* species 2

ATGCA[^]T
T[^]ACGTA

E145
E146

1000 u.a.
5000 u.a.

Concentration: 5 000 – 20 000 units/ml

Assayed on λ DNA

Reagents Supplied with Enzyme:

10 \times SE Buffer B, BSA

Reaction conditions: 1 \times SE Buffer B+BSA

Incubate at 37°C.

To obtain 100% activity, BSA should be

added to the 1 \times reaction mix to a final

concentration of 100 μ g/ml.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 200 mM KCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 μ g/ml BSA; 50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

Yes (65° C for 20 minutes)

Ligation/recutting assay: After 10-fold overdigestion with Zsp2 I, > 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 μ g of λ DNA with 20 units of Zsp2 I for 16 hours.

Do not use BSA for long incubation.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-----|-------|-------|-------|-------|
| Activity in SE Buffers (% of max) | 100 | 50-75 | 25-50 | 25-50 | 25-50 |

Nickases

Nicking Endonuclease N.Bst9 I

Isolated from *Bacillus stearothermophilus* T9

GAGTCNNNN[^]NN
CTCAGNNNNNN

E401

100 u.a.

E402

500 u.a.

Concentration: 2 000 – 5 000 units/ml

Assayed on T7 DNA

Reagents Supplied with Enzyme:

10×SE Buffer N.Bst9I

Reaction conditions: 1×SE Buffer N.Bst9 I

Incubate at 55°C.

Incubation at 37°C results in 20% activity.

Storage buffer: 10 mM Tris-HCl (pH 7.5);

50 mM KCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 µg/ml BSA;

50% glycerol. Store at -20°C.

Diluent: SE Buffer A

Heat inactivation:

No (65° C for 20 minutes)

Yes (80° C for 20 minutes)

Ligation/recutting assay: After 5-fold overdigestion with N.Bst9 I, 90% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of T7 DNA with 2 units of N.Bst9I for 16 hours.

Star activity: High enzyme concentration results in star activity.

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|--------|-----|-----|-------|
| Activity in SE Buffers (% of max) | 10-25 | 75-100 | 100 | 100 | 50-75 |

Methyl-directed DNA endonucleases

Aox I (prototype Aox I)

Isolated from *Arthrobacter oxydans* 25K



[^]PuG(5mC)Py
Py(5mC)GPu[^]

E569
E570

50 u.a.
250 u.a.

The enzyme cleaves C5-methylated DNA and doesn't cut unmodified DNA!

Concentration: 500 units/ml

Assayed on DNA pMHaeIII/DriI is a linearized plasmid pMHaeIII. pMHaeIII carries a gene of DNA-methyltransferase M.HaeIII, which methylates sites 5'-GGCC-3' producing 5'-GG(5mC)C-3'/3'-C(5mC)GG-5'.

One unit is defined as the amount of enzyme required to hydrolyze in 1 µg of linearized plasmid pMHaeIII/DriI in 1 hour at 60°C in a total reaction volume of 50 µl.

Reagents Supplied with Enzyme: 10×SE Buffer AoxI

Reaction conditions: 1× SE Buffer AoxI

Incubate at 60°C.

Storage buffer: 10 mM Tris-HCl (pH 7.4); 100 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 50% glycerol. **Store at -20°C.**

Heat inactivation:

No (80° C for 20 minutes)

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|-------|-------|-------|--------|
| Activity in SE Buffers (% of max) | 75-100 | 25-50 | 10-25 | 25-50 | 75-100 |

Overdigestion assay: No detectable degradation of 1µg of Lambda DNA was observed after incubation with 1 units of enzyme for 16 hours at 60°C in a total reaction volume of 50 µl.

Bis I (prototype Bis I)

Isolated from *Bacillus subtilis* T30

G(5mC)[^]NGC
CGN[^](5mC)G

E485
E486

40 u.a.
200 u.a.

The enzyme cleaves only C5-methylated DNA and doesn't cut unmodified DNA! [1]

Concentration: 1 000 – 2 000 units/ml

Assayed on Double-stranded oligonucleotide

5' GCTTGACTTTA G(5mC)G G C ATTGATTCTCACCACG 3'

3' CGAACATGAAAT C G C(5mC)G TAACTAAGAGTGGTGC 5'

One unit is defined as the amount of enzyme required to cleave 1 pmol of the double-stranded oligonucleotide with the following structure

5' GCTTGACTTTA G(5mC)G G CATTGATTCTCACCACG 3'

3' CGAACATGAAAT C G C(5mC)G TAACTAAGAGTGGTGC 5'

in 1 hour at 37°C in a total reaction volume of 20 µl.

Reagents Supplied with Enzyme: 10×SE Buffer BisI

Reaction conditions: 1× SE Buffer BisI

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 50% glycerol. **Store at -20°C.**

I.Chmuzh E.V., Kashirina J.G., Tomilova J.E., Mezentseva N.V., Dedkov V.S., Gonchar D.A., Abdurashitov M.A., Degtyarev S.Kh. A Novel Restriction Endonuclease BisI from *Bacillus subtilis* T30, Recognizes a Methylated DNA Sequence 5'-G(m5C)[^]NGC-3'// Biotekhnologia (Moscow), No.3, p.22-26 (2005) (In Russian). Online version in English: http://science.sibenzyme.com/article8_article_7_1.phtml

Heat inactivation:

Yes (65° C for 20 minutes)

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|-------|--------|-------|
| Activity in SE Buffers (% of max) | 10-25 | 25-50 | 50-75 | 75-100 | 50-75 |

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of pFsp4HI1 DNA (BamHI digest) with 1 unit of BisI for 16 hours at 37°C. The pFsp4HI1 plasmid carries a gene for Fsp4HI DNA-methyltransferase, which modifies DNA forming 5'-G(5mC)NGC-3'/3'-CGN(5mC)G-5'.



Bls I (prototype Bls I)Isolated from *Bacillus simplex* 23*DNA sequence with at least two 5mC:*PuPyN[^]PuPyPyPu[^]NPuPu

E533

100 u.a.

E534

500 u.a.

The enzyme cleaves C5-methylated DNA and doesn't cut unmodified DNA! [1]**Concentration:** 5000 units/ml**Assayed on** DNA pFsp4HI3/DriI is a linearized plasmid pFsp4HI3, which carries a gene of DNA-methyltransferase M.Fsp4HI and includes three canonical sites:

5'-G(5mC)NG(5mC)-3'/3'-CGN(5mC)G-5' [2].

One unit is defined as the amount of enzyme required to hydrolyze at least one of three canonical sites

5'-G(5mC)NG(5mC)-3'/3'-CGN(5mC)G-5'

in 1 µg of linearized plasmid pFsp4HI3 in 1 hour at 30°C

in a total reaction volume of 50 µl.

Reagents Supplied with Enzyme: 10×SE Buffer **W****Reaction conditions:** 1× SE Buffer **W**

Incubate at 30°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 200 mM NaCl; 0.1 mM EDTA;7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol. **Store at -20°C.****Overdigestion assay:** No detectable degradation of 1µg of Lambda DNA was observed after incubation with 5 units of enzyme for 16 hours at 30°C in a total reaction volume of 50 µl.1. Chernukhin V.A., Tomilova J.E., Chmuzh E.V., Sokolova O.O., Dedkov V.S., Degtyarev S.Kh. Bacterial strain *Bacillus simplex* - producer of BlsI site specific endonuclease. // Russian Federation patent RU 2322494 C1 (2006).

2. Chmuzh E.V., Kashirina Yu.G., Tomilova Yu.E., Chernukhin V.A., Okhapkina S.S., Gonchar D. A., Dedkov V.S., Abdurashitov M. A., Degtyarev S. Kh. Gene cloning, comparative analysis of the protein structures from Fsp4HI restriction-modification system and biochemical characterization of the recombinant DNA methyltransferase M.Fsp4HI. // Molecular Biology, V.41, No 1, p. 43-50 (2007)

Heat inactivation:**Yes** (65° C for 20 minutes)

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|-------|-----|--------|
| Activity in SE Buffers (% of max) | 10-25 | 10-25 | 50-75 | 100 | 75-100 |



New product



Mammalian Genome Qualified



New package

Gla I (prototype GlaI)

Isolated from Glacial ice bacterium GL 29

Pu(5mC)[^]GPy
PyG[^](5mC)Pu

E493
E494

100 u.a.
500 u.a.

The enzyme cleaves only C5-methylated DNA and does not cut unmodified DNA and DNA with N4-methylcytosines![1]

Concentration: 10 000 units/ml

Assayed on DNA pHspAI2/GsaI is a linearized plasmid pHspAI2, which carries a gene of DNA-methyltransferase M.HspAI (recognition sequence 5'-GCGC-3') and includes a unique GlaI recognition site:

5'-G(5mC)G(5mC)-3'/3'-(5mC)G(5mC)G-5' [2].

One unit is defined as the amount of enzyme required to hydrolyze completely a unique 5'-G(5mC)G(5mC)-3'/3'-(5mC)G(5mC)G-5' site in 1 µg of pHspAI2 plasmid DNA, which is linearized with GsaI, in 1 hour at 30°C in a total reaction volume of 50 µl.

Reagents Supplied with Enzyme:

10×SE Buffer Y, pHspAI2/GsaI DNA.

pHspAI2/GsaI DNA is a linearized plasmid pHspAI2, which includes a gene of DNA-methyltransferase M.HspAI (recognition sequence 5'-GCGC-3') and contains a unique GlaI canonical site: 5'-G(5mC)G(5mC)-3'/3'-(5mC)G(5mC)G-5'.

Reaction conditions: 1× SE Buffer Y.

Incubate at 30°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 200 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 0,05% Triton X-100; 100 µg/ml BSA, 50% glycerol.

Store at -20°C.

Overdigestion assay: No detectable degradation of 1µg of Lambda DNA was observed after incubation with 8 units of enzyme for 16 hours at 30°C in a total reaction volume of 50 µl.

1. Chernukhin V.A., Nayakshina T.N., Tomilova J.E., Mezentseva N.V., Dedkov V.S., Degtyarev S.Kh. Bacterial strain Glacial ice bacterium I - producer of Glal restriction endonuclease. // Russian Federation patent RU 2287012 C1 (2006).
2. Chernukhin V.A., Najakshina T.N., Abdurashitov M.A., Tomilova J.E., Mezentseva N.V., Dedkov V.S., Mikhnenkova N.A., Gonchar D.A., Degtyarev S. Kh A novel restriction endonuclease Glal recognizes methylated sequence 5'-G(5mC)[^]GC-3' // Biotechnologia V.4. p.31-35(2006) (In Russian). Online version in English: http://science.sibenzyme.com/article8_article_11_1.phtml
3. Tarasova G. V., Nayakshina T. N., Degtyarev S. Kh. Substrate specificity of new methyl-directed DNA endonuclease GlaI . // BMC Molecular Biology 2008, 9:7.
4. Tomilova J.E., Chernukhin V.A., Degtyarev S.Kh. Dependence of site-specific endonuclease GlaI activity on quantity and location of methylcytosines in the recognition sequence 5'-GCGC-3'. // Bulletin of biotechnology and physico-chemical biology V.2, No 1, p.30-39 (2006) (In Russian). Online version in English: http://science.sibenzyme.com/article10_article_22_1.phtml
5. Chernukhin V.A., Abdurashitov M.A., Tomilov V.N., Gonchar D.A., Degtyarev S.Kh. Comparative analysis of mouse chromosomal DNA digestion with restriction endonucleases in vitro and in silico // "Ovchinnikov bulletin of biotechnology and physical and chemical biology" V.3, No 4, pp 19-27, 2007. Online version in English: http://science.sibenzyme.com/article14_article_46_1.phtml
6. Abdurashitov M.A., Chernukhin V.A., Gonchar D.A., Degtyarev S.Kh. GlaI digestion of mouse γ-satellite DNA: study of primary structure and ACGT sites methylation.// BMC Genomics 2009, 10:322.
7. D. A. Gonchar, A. G. Akishev, S. Kh. Degtyarev BIsI- and GlaI-PCR assays – a new method of DNA methylation study // "Ovchinnikov bulletin of biotechnology and physical and chemical biology" V.6, No 1, pp 5-12, 2010
Online version in English: http://science.sibenzyme.com/article12_article_53_1.phtml

Heat inactivation:

Yes (65° C for 20 minutes)

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|--------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 25-50 | 25-50 | 100 |

The enzyme activity depending on number and position of methylated nucleotides in the recognition sequence:

| recognition sequence | G(mC)G(mC) (mC)G(mC) G | R(mC)G(mC) Y G(mC) G | G (mC)R(mC) (mC) G Y G |
|----------------------|---------------------------|-------------------------|---------------------------|
| Activity, % | 100 | > 25 | > 6 |

Glu I (prototype GluI)

Isolated from Glacial ice bacterium GL 24

G(5mC)[^]NG(5mC)
(5mC)GN[^](5mC)G

E519

50 u.a.

E520

250 u.a.

The enzyme cleaves C5-methylated DNA and does not cut unmodified DNA! [1]

Concentration: 1 000 units/ml

Assayed on DNA pFsp4HI3/DriI is a linearized plasmid pFsp4HI3, which includes a gene of DNA-methyltransferase M.Fsp4HI and contains a unique GluI site:

5'-G(5mC)NG(5mC)-3'/3'-(5mC)GN(5mC)G-5'. [2]

One unit is defined as the amount of enzyme required to hydrolyze completely a unique canonical site:

5'-G(5mC)NG(5mC)-3'/3'-(5mC)GN(5mC)G-5'

in 1 µg of linearized plasmid pFsp4HI3/DriI in 1 hour at 37°C

in a total reaction volume of 50 µl.

Reagents Supplied with Enzyme:

10×SE Buffer Y

Reaction conditions: 1× SE Buffer Y

Incubate at 37°C.

Storage buffer: 10 mM KH₂PO₄ (pH 7.45); 200 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol.

Store at -20°C.

Overdigestion assay: No detectable degradation of 1µg of Lambda DNA was observed after incubation with 1 units of enzyme for 16 hours at 37°C in a total reaction volume of 50 µl.

1. Chernukhin V.A., Chmuzh E.V., Tomilova J.E., Nayakshina T.N., Dedkov V.S., Degtyarev S.Kh. Bacterial strain Glacial ice bacterium - producer of GluI site specific endonuclease. // Russian Federation patent RU 2322492 C1 (2006).

2. Chernukhin V.A., Chmuzh E.V., Tomilova Yu.E., Nayakshina T.N., Gonchar D.A., Dedkov V.S., Degtyarev S.Kh. A novel site-specific endonuclease GluI recognizes methylated DNA sequence 5'-G(5mC)[^]NG(5mC)-3'/3'-(5mC)GN[^](5mC)G.// Bulletin of biotechnology and physico-chemical biology named by Yu.A.Ovchinnikov (Moscow), V.3, No.2, p.13-17 (2007) (In Russian).

Online version in English: http://science.sibenzyme.com/article8_article_24_1.phtml

3. Chernukhin V.A., Abdurashitov M.A., Tomilov V.N., Gonchar D.A., Degtyarev S.Kh. Comparative analysis of mouse chromosomal DNA digestion with restriction endonucleases in vitro and in silico // "Ovchinnikov bulletin of biotechnology and physical and chemical biology" V.3, No 4, pp 19-27, 2007. Online version in English: http://science.sibenzyme.com/article14_article_46_1.phtml

Heat inactivation:

Yes (80° C for 20 minutes)

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|--------|--------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 75-100 | 75-100 | 25-50 | 50-75 | 100 |



New product



Mammalian Genome Qualified



New package

Kro I (prototype Kro I)

Isolated from *Kocurea rosea* 307

G[^]C(5mC)GGC
CGG(5mC)C[^]G

E541 50 u.a.
E542 250 u.a.

The enzyme cleaves C5-methylated DNA and doesn't cut unmodified DNA! [1]

Kro I doesn't cleave DNA modified with MspI DNA-methyltransferase

Concentration: 1 000 units/ml

Assayed on DNA pMHpaII1/DriI is a linearized plasmid pMHpaII1.

pMHpaII1 carries a gene of DNA-methyltransferase M.HpaII, which methylates sites

5'-CCGG-3' producing 5'-C(5mC)GG-3'/3'-GG(5mC)C-5',

and includes three canonical sites

5'-GC(5mC)GGC-3'/3'-CGG(5mC)CG-5'.

One unit One unit is defined as the amount of enzyme required to hydrolyze completely 1 µg of linearized plasmid pMHpaII1 in 1 hour at 37°C in a total reaction volume of 50 µl.

Reagents Supplied with Enzyme:

10×SE Buffer G

Reaction conditions: 1× SE Buffer G

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 200 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 200 µg/ml BSA; 50% glycerol.

Store at -20°C.

Overdigestion assay: No detectable degradation of 1µg of Lambda DNA was observed after incubation with 1 units of enzyme for 16 hours at 30°C in a total reaction volume of 50 µl.

1. Chernukhin V.A., Zhuravleva R.O., Tarasova G.V., Boltengagen A. A., Akishev A.G., Mikhnenkova N.A., Degtyarev S.Kh. Bacterial strain *Kocuria rosea* - producer of Krol site specific endonuclease. // Russian Federation patent RU 2394099 C1 (2010).

Heat inactivation:

Yes (65° C for 20 minutes)

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-----|-------|-------|--------|
| Activity in SE Buffers (% of max) | 50-75 | 100 | 25-50 | 50-75 | 75-100 |

Mal I (prototype Dpn I)

Isolated from *Marinococcus albus* I

G(mA)[^]TC
CT[^](mA)G

E489 50 u.a.
E490 250 u.a.

The enzyme cleaves only methylated DNA

Concentration: 500 - 1 000 units/ml

Assayed on

pBR322 DNA (dam-methylated)

One unit One unit of the enzyme is the amount required to hydrolyze 1 µg of pBR322 DNA (dam-methylated) in 1 hour at 37°C in a total reaction volume of 50 µl.

Reagents Supplied with Enzyme:

10×SE Buffer Mal I

Reaction conditions: 1×SE Buffer Mal I

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 100 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 100 µg/ml BSA, 50% glycerol.

Diluent: SE Buffer A

Store at -20°C.

Heat inactivation:

Yes (65° C for 20 minutes)

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|-------|--------|-------|
| Activity in SE Buffers (% of max) | 25-50 | 25-50 | 50-75 | 75-100 | 50-75 |

Ligation/recutting assay: After 2-fold overdigestion with Mal I, ~80% of the DNA fragments can be ligated and recut.

Overdigestion assay: No nonspecific activity was detected after incubation of 1 µg of λ DNA with 2 units of Mal I for 16 hours.

Mte I (prototype Mte I)Isolated from *Microbacterium testaceum* 17B**G(5mC)G(5mC)^NG(5mC)G(5mC)
(5mC)G(5mC)GN^(5mC)G(5mC)G****E553 500 u.a.
E554 2500 u.a.****The enzyme cleaves only C5-methylated DNA
and does not cut unmodified DNA ! [1]****Concentration:** 10 000 units/ml**Assayed on** pHspAI10/DriI+M.Fsp4HI is a plasmid pHspAI10, which is linearized with DriI, and, additionally, modified with Fsp4HI DNA methyltransferase. pHspAI10 carries a gene of HspAI DNA methyltransferase, that modifies the sequence

5'-GCGC-3', producing 5'-G(5mC)GC-3'.

M.Fsp4HI modifies the sequence

5'-GCNGC-3', producing 5'-G(5mC)NGC-3'.

A substrate pHspAI10/DriI+M.Fsp4HI includes one site

5'-G(5mC)G(5mC)NG(5mC)G(5mC)-3'/3'-(5mC)G(5mC)GN(5mC)G(5mC)G-5',

which is MteI canonical site [1]. The enzyme activity depends on a number and positions of methylated nucleotides in the recognition sequence. For example, MteI cuts the recognition site

with six 5-methylcytosines, but the enzyme activity is reduced for more than one order [1].

One unit is defined as the amount of enzyme required to hydrolyze completely 1 µg of linearized plasmid pHspAI10/DriI+M.Fsp4HI in 1 hour at 55°C in a total reaction volume of 50 µl.**Reagents Supplied with Enzyme:**

10×SE Buffer W

Reaction conditions: 1 x SE-buffer W

Incubate at 55°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 200 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 200 µg/ml BSA, 50% glycerol.

Store at -20°C.**Overdigestion assay:** No detectable degradation of 1µg of Lambda DNA was observed after incubation with 10 units of enzyme for 16 hours at 55°C in a total reaction volume of 50 µl.

1. V.A. Chernukhin, E.V. Kileva, V.A. Sokolova., D.A. Gonchar, L.N. Golikova, V.S. Dedkov, N.A. Mikhnenkova, S.Kh. Degtyarev.

A new methyl-directed site-specific DNA endonuclease MteI cleaves nine nucleotides sequence 5'-G(5mC)G(5mC)^NG(5mC)GC-3'/3'-CG(5mC)GN^(5mC)G(5mC)G-5' // "Ovchinnikov bulletin of biotechnology and physical and chemical biology" V.8, No 1, pp 16-26, 2012

Heat inactivation:**No** (80° C for 20 minutes)

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|--------|--------|-----|-------|
| Activity in SE Buffers (% of max) | 25-50 | 75-100 | 75-100 | 100 | 50-75 |

Pcs I (prototype Pcs I)Isolated from *Paracoccus carotinifaciens* 3K**(5mC)GNNNNN^NN(5mC)G
G(5mC)NN^NNNNNG(5mC)****E505 50 u.a.
E506 250 u.a.****The enzyme cleaves only C5-methylated DNA
and does not cut unmodified DNA ! [1]****Concentration:** 1 000 units/ml**Assayed on** pMHgaI/DriI is a linearized plasmid pMHgaI, which carries genes of DNA-methyltransferases M1.HgaI (recognition sequence 5'-GCGTC-3')

and M2.HgaI (5'-GACGC-3') and includes a unique PcsI canonical site:

5'-W(5mC)GNNNNNNN(5mC)GW-3'/3'-WG(5mC)NNNNNNNG(5mC)W-5' [1].

One unit defined as the amount of enzyme required to digest a unique site

5'-A(5mC)GNNNNNNN(5mC)GT-3'

in 1 µg of DNA pMHgaI/DriI in 1 hour at 37°C in a total reaction volume of 50 µl.

Reagents Supplied with Enzyme:

10×SE Buffer Pcs I

Reaction conditions: 1 x SE-buffer PcsI

Incubate at 37°C.

Storage buffer: 10 mM Tris-HCl (pH 7.5); 200 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 0.1 mg/ml BSA, 50% glycerol.

Store at -20°C.**Overdigestion assay:** No detectable degradation of 1µg of Lambda DNA was observed after incubation with 1 units of enzyme for 16 hours at 37°C in a total reaction volume of 50 µl.**Notes:** When using a buffer other than the optimal (supplied) SEBuffer, it may be necessary to add more enzyme to achieve complete digestion.1. Chernukhin V.A., Nayakshina T.N., Tarasova M.V., Golikova L.N., Akishev A.G., Dedkov V.S., Mikhnenkova N.A., Degtyarev S.Kh. Bacterial strain *Paracoccus carotinifaciens* 3K- producer of PcsI site specific endonuclease. // Russian Federation patent RU 2377294 C1 (2009).**Heat inactivation:****Yes** (65° C for 20 minutes)

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|-------|---|-------|-------|
| Activity in SE Buffers (% of max) | 50-75 | 25-50 | 0 | 10-25 | 50-75 |



New product



Mammalian Genome Qualified



New package

Pkr I (prototype Pkr I)

Isolated from
Planomicrobium koreense 78k



DNA sequence with at least
three 5mC:

G(5mC)N^G(5mC)
(5mC)G^N(5mC)G

E579

50 u.a.

E580

250 u.a.

The enzyme cleaves C5-methylated DNA and doesn't cut unmodified DNA! [1]

Concentration: 1000 units/ml

Assayed on DNA pFsp4HI3/DriI is a linearized plasmid pFsp4HI3, which carries a gene of DNA-methyltransferase M.Fsp4HI and includes three sites:

5'-G(5mC)NG(5mC)-3'/3'-CGN(5mC)G-5' [2].

One unit is defined as the amount of enzyme required to hydrolyze completely 1 µg of linearized plasmid pFsp4HI3 in 1 hour at 37°C in a total reaction volume of 50 µl

Reagents Supplied with Enzyme: 10×SE Buffer Y

Reaction conditions: 1× SE Buffer Y

Incubate at 37°C.

Storage buffer: 20 mM Tris-HCl (pH 7.4); 200 mM NaCl; 0.1 mM EDTA;

7 mM 2-mercaptoethanol; 50% glycerol. **Store at -20°C.**

Overdigestion assay: No detectable degradation of 1 µg of Lambda DNA was observed after incubation with 2 units of enzyme for 16 hours at 37°C in a total reaction volume of 50 µl.

1. V.A. Chernukhin, T.N. Nayakshina, D.A. Gonchar, Ju.E. Tomilova, M.V. Tarasova, V.S. Dedkov, N.A. Mikhnenkova, S.Kh. Degtyarev A new site-specific methyl-directed DNA endonuclease PkrI recognizes and cuts methylated DNA sequence 5'-GCN^GC-3'/3'-CG^NCG-5' carrying at least three 5-methylcytosines. // Bulletin of biotechnology and physico-chemical biology named by Yu.A. Ovchinnikov (Moscow), V.7, No.3, p.35-42 (2011). (In Russian).

2. Chmuzh E.V., Kashirina Yu.G., Tomilova Yu.E., Chernukhin V.A., Okhapkina S.S., Gonchar D. A., Dedkov V.S., Abdurashitov M. A., Degtyarev S. Kh. Gene cloning, comparative analysis of the protein structures from Fsp4HI restriction-modification system and biochemical characterization of the recombinant DNA methyltransferase M.Fsp4HI. // Molecular Biology, V.41, No 1, p. 43-50 (2007)

Heat inactivation:

Yes (65° C for 20 minutes)

| SE Buffers | B | G | O | W | Y |
|-----------------------------------|-------|--------|-------|-------|-----|
| Activity in SE Buffers (% of max) | 50-75 | 75-100 | 10-25 | 25-50 | 100 |



DNA – Methyltransferases

M3.BstF5I

Isolated from *E.coli* strain, that carries the cloned M3.BstF5I gene from *Bacillus stearothermophilus* F5

CH₃
I
5'-GGATG-3'
3'-CCTAC-5'

M007

1000 u

Description: M3.BstF5I Methylase modifies the adenin residue (mA) in the recognition sequence 5'-GGATG-3'

Concentration: 10000 units/ml

Reagents Supplied: 10 × SEBuffer K, SAM

Incubate at 60°C.

Storage Conditions: 10 mM Tris-HCl (pH 7.5); 50 mM NaCl; 50% glycerol. **Store at -20°C.**

Unit Definition: One unit of the enzyme is the amount required to protect 1 µg of Lambda DNA in 1 hour at 60°C in a total reaction volume of 20 µl against cleavage by BstF5 I restriction endonuclease.

I. Golikova L.N., Netesova N.A., Gutorov V.V., Belavin P.A., Abdurashitov M.A., Gonchar D.A., Degtyarev S.Kh. Multiplicity of site-specific DNA-methyltransferases of the BstF5I restriction modification system from *Bacillus stearothermophilus* F5. // *Molecular Biology (Moscow)*, V.34, No.3, p.443-447 (2000) (In Russian).

M.Fsp4HI

Isolated from *E.coli* strain, that carries the cloned M.Fsp4HI gene from *Flavobacterium* species 4H

CH₃
I
5'-GCNGC-3'
3'-CGNCG-5'

M001

100 u

Description: M.Fsp4HI Methylase modifies the internal cytosine residue (C5) in the recognition sequence 5'-GCNGC-3'.

Concentration: 500-1000 units/ml

Reagents Supplied: 10 × SEBuffer B, SAM

Incubate at 30°C.

Storage Conditions: 10 mM Tris-HCl (pH 7.6); 100 mM NaCl; 0.1 mM EDTA;

10 mM 2-mercaptoethanol; 50% glycerol. **Store at -20°C.**

Unit Definition: One unit of the enzyme is the amount required to protect 1 µg of Lambda DNA in 1 hour at 30°C in a total reaction volume of 20 µl against cleavage by Fsp4HI restriction endonuclease.

I. Chmuzh E.V., Kashirina Yu.G., Tomilova Yu.E., Chernukhin V.A., Okhapkina S.S., Gonchar D. A., Dedkov V.S., Abdurashitov M. A., Degtyarev S. Kh. Gene cloning, comparative analysis of the protein structures from Fsp4HI restriction-modification system and biochemical characterization of the recombinant DNA methyltransferase M.Fsp4HI. // *Molecular Biology*, V.41, No 1, p.43-50 (2007) (In Russian).

M.HspAI

Isolated from *E.coli* strain, that carries the cloned M.HspAI gene from *Haemophilus* species AI

CH₃
I
5'-GCGC-3'
3'-CGCG-5'

M003

100 u

Description: M.HspAI Methylase modifies the internal cytosine residue (C5) in the recognition sequence 5'-GCGC-3'.

Concentration: 1000-3000 units/ml

Reagents Supplied: 10 × SEBuffer B, SAM

Incubate at 37°C.

Storage Conditions: 10 mM Tris-HCl (pH 7.6); 100 mM NaCl; 0.1 mM EDTA;

10 mM 2-mercaptoethanol; 50% glycerol. **Store at -20°C.**

Unit Definition: One unit of the enzyme is the amount required to protect 1 µg of Lambda DNA in 1 hour at 37°C in a total reaction volume of 20 µl against cleavage by HspAI restriction endonuclease.

I. Chernukhin V.A., Najakshina T.N., Abdurashitov M.A., Tomilova J.E., Mezentzeva N.V., Dedkov V.S., Mikhnenkova N.A., Gonchar D.A., Degtyarev S. Kh A novel restriction endonuclease GlAI recognizes methylated sequence 5'-G(5mC)⁵GC-3' // *Biotechnologia* V.4. p.31-35 (2006) (In Russian). Online version in English: http://science.sibenzyme.com/article8_article_11_1.phtml

R = A or G

W = A or T

S = G or C

K = G or T

M = A or C

Y = T or C

D = A or G or T

H = A or C or T

B = C or G or T

V = A or C or G

N = A or C or G or T



New product



Mammalian Genome Qualified



New package

Polymerases

DNA Polymerase I, Large (Klenow) Fragment

E325 200 u
E326 1000 u

Isolated from *E. coli* strain that carries the cloned DNA Polymerase I, Large (Klenow) Fragment gene

Description: DNA Polymerase I, Large (Klenow) Fragment is a proteolytic product of *E. coli* DNA Polymerase I which retains polymerization and 3'→5' exonuclease activity, but has lost 5'→3' exonuclease activity.

Applications:

- dideoxy sequencing;
- polishing ends;
- second strand cDNA synthesis

Concentration: 5 000 units/ml

Reagents Supplied with Enzyme: 10 × Klenow Buffer

Storage Conditions: 10 mM KH₂PO₄ (pH 7.5); 50 mM KCl; 0.5 mM EDTA; 7 mM 2-mercaptoethanol; 50% glycerol.

Store at -20°C.

Heat Inactivation: 75°C for 20 minutes.

Unit Definition: One unit is the amount of enzyme required to incorporate 10 nmol of dNTPs into an acid-insoluble material in 30 minutes at 37°C.

Unit Assay Conditions: 1 × Klenow Buffer, 33 μM dNTP including [³H]-dTTP and 70 μg/ml denatured calf thymus DNA.

Quality Control: Purified free of contaminating endonucleases and exonucleases.

M-MuLV Reverse Transcriptase RNase H -

E317 5000 u
E318 25000 u

Isolated from *E. coli* strain that carries the recombinant plasmid

Description: Moloney Murine Leukemia Virus (M-MuLV) Reverse Transcriptase is an RNA directed DNA polymerase. This enzyme can synthesize a complementary DNA strand initiating from a primer using either RNA (cDNA synthesis) or single-stranded DNA as a template. The absence of RNase H activity enhances the synthesis of long cDNAs and therefore the enzyme is recommended for preparing long cDNAs.

Application: first strand cDNA synthesis.

Concentration: 50 000–200 000 units/ml

Reagents Supplied with Enzyme:

10 × M-MuLV Reverse Transcriptase Buffer

Dilution Buffer M-MuLV Reverse Transcriptase

Storage Conditions: 10 mM KH₂PO₄ (pH 7.5); 0,1 mM EDTA; 200 mM NaCl; 7 mM 2-mercaptoethanol; 50% glycerol. **Store at -20°C.**

Unit Definition: One unit is the amount of the enzyme required to incorporate 1 nmol of dTTP into an acid-insoluble material in 10 minutes at 37°C using poly(rA)•oligo(dT).

Unit Assay Conditions: 50 mM Tris-HCl (pH 8.3), 6 mM MgCl₂, 10 mM DTT, 0.5 mM [³H]-dTTP, 0.4 mM poly(rA)•oligo(dT) 12-18.

Quality Control: Purified free of detectable levels of RNase, endonuclease and exonuclease activities.

Note: High enzyme concentration may lead to RT-PCR inhibition. In this case the enzyme preparation should be diluted in 5, 10 or 20 times with M-MuLV Reverse Transcriptase dilution Buffer

T4 DNA Polymerase

E339 200 u.
E340 1000 u

Isolated from *E. coli* strain that carries the recombinant plasmid

Description: T4 DNA Polymerase catalyzes the synthesis of DNA in the 5'→3' direction and requires the presence of template and primer. This enzyme has a 3'→5' exonuclease activity.

Applications:

- polishing ends;
- probe labeling using replacement synthesis.

Concentration: 2 000 – 5 000 units/ml

Reagents Supplied with Enzyme: 10 × T4 DNA Polymerase Buffer

Storage Conditions: 20 mM Tris-HCl (pH 7.5); 50 mM KCl;

10 mM 2-mercaptoethanol; 50% glycerol. **Store at -20°C.**

Heat Inactivation: 75°C for 10 minutes.

Unit Definition: One unit is the amount of enzyme required to incorporate 10 nmol of dNTP into acid-precipitable material in 30 minutes at 37°C.

Unit Assay Conditions: 1 × T4 DNA Polymerase Reaction Buffer, 33 μM dNTPs including [³H]-dTTP and 70 μg/ml denatured calf thymus DNA.

Quality Control: Purified free of contaminating endonucleases and exonucleases.

SP-Taq DNA Polymerase

Isolated from *E.coli* strain that carries the cloned Taq DNA polymerase gene from *Thermus aquaticus*



E333 200 u
E334 1000 u

Description: SP-Taq DNA Polymerase is a fraction of Taq DNA Polymerase which was specially treated and additionally purified. It doesn't contain PCR detected DNA contaminations. Enzyme is suitable for different manipulations in the field of PCR-diagnostics.

Applications:

- nick translation;
- primer extension reaction.

Concentration: 5 000 units/ml

Reagents Supplied with Enzyme: 10 × Taq-DNA-polymerase buffer.

Storage conditions: 20 mM Tris-HCl (pH 7.6); 100 mM KCl; 0.1 mM EDTA; 1 mM DTT; 0.5 % Triton X-100; 50% glycerol.

Store at -20°C.

Unit Definition: One unit is the amount of enzyme required to incorporate 10 nmol of dNTP into acid-insoluble material in 30 minutes at 72°C.

Unit Assay Conditions: 60 mM Tris-HCl (pH 8.5 at 25°C); 25 mM KCl; 1.5 mM MgCl₂, 10 mM 2-mercaptoethanol, 0.1% Triton X-100, 200 μM dATP, dCTP, dGTP, 50 μM H-TTP, 12.5 μg activated Calf Thymus DNA in a total reaction volume of 50 μl.

Quality Control: Purified free of contaminating endonucleases and exonucleases.

Taq DNA Polymerase with Standard Taq Buffer

Isolated from *E.coli* strain that carries the cloned Taq DNA polymerase gene from *Thermus aquaticus*

E331 200 u
E332 1000 u

Description: Taq DNA polymerase is a thermostable DNA polymerase that possesses a 5'→3' polymerase activity and a double strand specific 5'→3' exonuclease activity.

Applications:

- nick translation;
- primer extension reaction.

Concentration: 5 000 units/ml

Reagents Supplied with Enzyme: 10 × Taq-DNA-polymerase buffer.

Storage conditions: 20 mM Tris-HCl (pH 7.6); 100 mM KCl; 0.1 mM EDTA; 1 mM DTT; 0.5 % Triton X-100; 50% glycerol.

Store at -20°C.

Unit Definition: One unit is the amount of enzyme required to incorporate 10 nmol of dNTP into acid-insoluble material in 30 minutes at 72°C.

Unit Assay Conditions: 1 × Taq-DNA-polymerase buffer, 200 μM dNTPs including [³H]-dTTP and 250 μg/ml activated calf thymus DNA.

Quality Control: Purified free of contaminating endonucleases and exonucleases.

Taq DNA Polymerase with Standard Taq Buffer (Mg²⁺ free)

Isolated from *E.coli* strain that carries the cloned Taq DNA polymerase gene from *Thermus aquaticus*

E341 200 u
E342 1000 u

Description: Taq DNA polymerase is a thermostable DNA polymerase that possesses a 5'→3' polymerase activity and a double strand specific 5'→3' exonuclease activity.

Applications:

- nick translation;
- primer extension reaction.

Concentration: 5 000 units/ml

Reagents Supplied with Enzyme:

10 × Taq-DNA-polymerase buffer (Mg²⁺ free).

Storage conditions: 20 mM Tris-HCl (pH 7.6); 100 mM KCl; 0.1 mM EDTA; 1 mM DTT; 0.5 % Triton X-100; 50% glycerol.

Store at -20°C.

Unit Definition: One unit is the amount of enzyme required to incorporate 10 nmol of dNTP into acid-insoluble material in 30 minutes at 72°C.

Unit Assay Conditions: 1 × Taq-DNA-polymerase buffer, 200 μM dNTPs including [³H]-dTTP and 250 μg/ml activated calf thymus DNA.

Quality Control: Purified free of contaminating endonucleases and exonucleases.

Taq DNA Polymerase with AS Buffer

Isolated from *E.coli* strain that carries the cloned Taq DNA polymerase gene from *Thermus aquaticus*

E337 200 u
E338 1000 u

Description: Taq DNA polymerase is a thermostable DNA polymerase that possesses a 5'→3' polymerase activity and a double strand specific 5'→3' exonuclease activity.

Applications:

- nick translation;
- primer extension reaction.

Concentration: 5 000 units/ml

Reagents Supplied with Enzyme:

10 × AS (Ammonium Sulfate) buffer.

Storage conditions: 20 mM Tris-HCl (pH 7.6); 100 mM KCl; 0.1 mM EDTA; 1 mM DTT; 0.5 % Triton X-100; 50% glycerol.

Store at -20°C.

Unit Definition: One unit is the amount of enzyme required to incorporate 10 nmol of dNTP into acid-insoluble material in 30 minutes at 72°C.

Unit Assay Conditions: 1 × AS (Ammonium Sulfate) buffer, 200 μM dNTPs including [³H]-dTTP and 250 μg/ml activated calf thymus DNA.

Quality Control: Purified free of contaminating endonucleases and exonucleases.

TaqSE DNA Polymerase

E313 200 u
E314 1000 u

Isolated from *E.coli* strain that carries the recombinant plasmids.

Description: TaqSE DNA polymerase is a complex mix of thermostable DNA polymerase that possesses a 5'→3' polymerase activity, 3'→5' exonuclease (proofreading) activity and a double strand specific 5'→3' exonuclease activity.

It may increase yield of reaction product compare to Taq DNA polymerase.

Application: long high fidelity primer extension reaction.

Concentration: 5 000 units/ml

Reagents Supplied with Enzyme: 10 × Taq-DNA-polymerase buffer.

Storage conditions: 20 mM Tris-HCl (pH 7.6); 100 mM KCl;

0.1 mM EDTA; 1 mM DTT; 0.5 % Triton X-100; 50% glycerol.

Store at -20°C.

Unit Definition: One unit is the amount of enzyme required to incorporate 10 nmol of dNTP into acid-insoluble material in 30 minutes at 72°C.

Unit Assay Conditions: 1 × Taq-DNA-polymerase buffer, 200 μM dNTPs including [³H]-dTTP and 250 μg/ml activated calf thymus DNA.

Quality Control: Purified free of contaminating endonucleases and exonucleases.

Hot Start Taq-DNA Polymerase

E351 200 u
E352 1000 u

Description: Hot Start Taq DNA Polymerase is complex mixture of a thermostable 94 kD Taq DNA Polymerase purified from *E.coli* recombinant strain expressing *Thermus aquaticus* polymerase gene and specific monoclonal antibodies from mouse.

Hot Start Taq DNA Polymerase is inactive under conditions of amplification reaction preparation. It can eliminate amplification artefacts such as primer-dimer formation and mispriming during preamplification stage and thus may provide improved specificity when compared to standard DNA polymerases.

An advantage of Hot Start Taq DNA Polymerase is the absence of additional heating step for polymerase activation. Heat activation of enzyme occurs during the first denaturation step. An inactive complex of Hot Start Taq DNA Polymerase dissociates automatically over +70°C, allowing activation of DNA polymerase.

Applications: -Highly specific PCR;

-Multiplex PCR (highly recommended);

-High sensitivity applications.

Concentration: 5 000 units/ml

Reagents Supplied with Enzyme:

10 × Hot Start Taq-DNA-polymerase buffer.

Supplementary material is 50 mM MgCl₂.

Storage and dilution buffer: 20 mM Tris-HCl, pH 8.0, 100 mM KCl, 0,1 mM EDTA, 1 mM DTT, 50 % glycerol, 0.5 % Nonidet P-40, 0.5 % Tween-20. **Store at -20°C.**

Unit Definition: One unit of activity is the amount of enzyme required to incorporate 10 nmoles of dNTP into acid-insoluble material in 30 minutes at 74 °C.

Notes: The recommended amount of enzyme is 1 u per 50μl of a total reaction volume.

Pfu DNA Polymerase

E353 200 u
E354 1000 u

Isolated from *E.coli* strain that carries the cloned DNA Polymerase gene from *Pyrococcus furiosus*.

Description: Pfu DNA Polymerase catalyzes the DNA-dependent polymerization of nucleotides into duplex DNA in the 5'→3' direction in the presence of magnesium ions. The enzyme also exhibits 3'→5' exonuclease (proofreading) activity, that enables the polymerase to correct nucleotide incorporation errors. Products of reaction have blunt ends.

Application: Pfu DNA Polymerase useful for high fidelity synthesis and polishing ends.

Concentration: 5 000 units/ml

Reagents Supplied with Enzyme:

10 × Pfu DNA polymerase buffer, 50 × BSA (5 mg/ml BSA).

Storage conditions: 10 mM K₂HPO₄ (pH 7.4); 0,1 mM DTT, 0.1 mM EDTA; 0.5 % Tween 20; 50% glycerol.

Store at -20°C.

To obtain 100% activity, BSA should be added to the 1 × reaction mix to a final concentration of 100 μg/ml.

Unit Definition: One unit of enzyme catalyzes the incorporation of 10 nanomoles of deoxyribonucleotides into a polynucleotide fraction (adsorbed on DE- 81) in 30 minutes at 72°C.

Unit Assay Conditions: 1 × Pfu DNA polymerase buffer, 0.1 mg/ml BSA, 200μg/ml activated calf thymus DNA, 0.2 mM of each dNTP, 0.4MBq/ ml [³ H]- dTTP in 50 μl reaction mix.

Quality Control: Purified free of contaminating endonucleases and exonucleases.

Notes: 1. Do **not** use dU-contained templates.

2. This enzyme is **not** recommended for a number of experiments dealing with very low-annealing temperature amplification approaches (e.g. RAPDs, Random Amplified Polymorphic DNAs).

3. Prepare PCR mix at 0°C.

Set in amplificator preheated to 95°C.

4. The recommended amount of enzyme is 2.5 u per 50μl of a total reaction volume.

T7 RNA Polymerase

E355 5000 u
E356 25000 u

Isolated from *E.coli* strain that carries the cloned T7 gene I

Description: T7 RNA Polymerase catalyzes the synthesis of RNA in the 5'→3' direction in the presence of a DNA template containing a T7 phage promoter.

Application: -Radiolabeled RNA probe preparation
-RNA generation for *in vitro* translation
-RNA generation for studies of RNA structure, processing and catalysis

Concentration: 100 000 units/ml

Reagents Supplied with Enzyme: 10 × T7 RNA polymerase buffer

Storage conditions: 50 mM Tris-HCl (pH 7.5); 100 mM NaCl; 20 mM 2-mercaptoethanol; 1 mM EDTA; 50% Glycerol; 0.1% Triton X-100.

Store at -20°C.

Unit Definition: One unit is the amount of enzyme required to incorporate 1 nmol NTP into an acid-insoluble material in 60 minutes at 37°C.

Quality Control: Purified free of contaminating endonucleases and exonucleases.

DNA Modifying Enzymes

Alkaline Phosphatase, Calf Intestinal

Calf intestinal mucosa

E327 100 u
E328 500 u

Description: Alkaline phosphatase catalyzes the removal of 5' phosphate groups from DNA, RNA and ribo- and deoxyribonucleoside triphosphates.

Applications:

- removing 5' and 3' phosphoryl groups from nucleic acids;
- preparing templates for 5' end labeling;
- preventing fragments from self ligating;
- dephosphorylation of proteins

Concentration: 10 000 units/ml.

Reagents Supplied with Enzyme: 10 × SEBuffer O

Storage conditions: 10 mM Tris-HCl (pH 8.2); 50 mM KCl; 1 mM MgCl₂, 0,1 mM ZnCl₂, 50% glycerol.

Store at -20°C.

Unit Definition: One unit is the amount of enzyme that hydrolyzes 1 μmol of p-nitrophenylphosphate to p-nitrophenol in a total reaction volume of 1 ml in 1 min at 37°C.

Unit Assay Conditions: 1 M diethanolamine-HCl (pH 9.8), 0.5 mM MgCl₂ and 10 mM p-nitrophenylphosphate. These conditions are only used for quantitating enzyme activity.

Quality Control: Purified free of contaminating exonuclease, endonuclease, ribonuclease activities.

Thermolabile Alkaline PhosphataseIsolated from *E.coli* strain that carries the cloned *Alkaline Phosphatase* gene from *Alteromonas undina* P2E365 200 u
E366 1000 u

Description: Thermolabile Alkaline phosphatase catalyzes the removal of 5' phosphate groups from DNA, RNA and ribo- and deoxyribonucleoside triphosphates.

Concentration: 5 000 units/ml.

Reagents Supplied with Enzyme: 10 × SEBuffer W

Heat inactivation: Yes (65° C for 20 minutes)

Storage conditions: 20 mM Tris-HCl (pH 7.6); 0,1 mM ZnCl₂, 50% glycerol.

Store at -20°C.

Unit Definition: One unit is the amount of enzyme that will dephosphorylate 1 μg of pUC19 DNA (linearized with Hind III) in 30 minutes at 25°C.

Dephosphorylation is defined as >95% inhibition of recirculation in a self-ligation reaction that is measured by transformation into *E.coli*.

Quality Control: Purified free of contaminating exonuclease, endonuclease, ribonuclease activities.

Endonuclease IIsolated from *Proteus vulgaris* 84E323 4000 u
E324 20000 u

Description: Endonuclease I hydrolyzes double- and single-stranded nucleic acids to oligonucleotides of 3-5 nucleotide in length with 5'-terminal phosphates.

Application: DNA and RNA degradation.

Concentration: 10 000 – 50 000 units/ml.

Reagents Supplied with Enzyme: 10 × Endonuclease I Buffer.

Storage Conditions: 10 mM Tris-HCl (pH 7.4); 250 mM NaCl; 0,2 mM EDTA; 7 mM 2-mercaptoethanol; 100 μg/ml BSA; 50% glycerol.

Store at -20°C.

Unit Definition: One unit is the amount of enzyme that hydrolyze 1 μg of Lambda DNA in 30 minutes at 37°C.

Unit Assay Conditions: 20 mM Glycine-NaOH (pH 9.5), 100 mM NaCl, 25 mM MgCl₂; 1 mM 2-mercaptoethanol.

Quality Control: Purified free of contaminating phosphatase activity.

Exonuclease III (*E. coli*)

Isolated from a recombinant source

E345 5000 u
E346 25000 u

Description: Exonuclease III (*E. coli*) catalyzes the stepwise removal of mononucleotides from 3'-hydroxyl termini of double stranded DNA. Exonuclease III activity depends partially on helical structure and displays sequence dependence (C>A=T>G). Temperature, salt concentration and the ratio of enzyme to DNA greatly affect enzyme activity, requiring reaction conditions to be tailored to specific applications.

Applications:

- unidirectional nested deletions;
- site-directed mutagenesis;
- preparation of strand-specific probes;
- preparation of single-stranded substrates for dideoxy sequencing.

Concentration: 40 000 – 100 000 units/ml.

Reagents Supplied with Enzyme: 10 × Exonuclease III Buffer .

Storage Conditions: 10 mM Tris-HCl (pH 7.5); 50 mM KCl; 0.5 mM EDTA; 1 mM 2-mercaptoethanol; 50% glycerol. **Store at -20°C.**

Unit Definition: One unit is the amount of enzyme required to produce 1 nmol of acid soluble total nucleotide in 30 minutes at 37°C.

Unit Assay Conditions: 50 mM Tris-HCl (pH 7.5 at 25°C); 1 mM MgCl₂; 1 mM DTT,

0.15 mM sonicated pancreatic DNA.

RNA Ligase T4

Isolated from *E.coli* strain that carries the cloned *RNA Ligase* gene from bacteriophage T4

Description: RNA Ligase catalyzes ligation of a 5' phosphoryl-terminated nucleic acid donor to a 3' hydroxyl-terminated nucleic acid acceptor through the formation of a 3'→5' phosphodiester bond, with hydrolysis of ATP to AMP and PPi. Substrates include single-stranded RNA and DNA as well as dinucleoside pyrophosphates.

Applications:

- labeling of 3'-termini of RNA with 5'-[³²P] pCp;
- inter- and intra-molecular joining of RNA and DNA molecules.

Concentration: 10 000 units/ml.

Reagents Supplied with Enzyme: 10 × T4 RNA Ligase Buffer.

E349 **1000 u**

E350 **5000 u**

Storage conditions: 10 mM Tris-HCl (pH 7.4); 50 mM KCl; 0.1 mM EDTA; 1 mM DTT; 50% glycerol.

Store at -20°C.

Unit Definition: One unit is defined as the amount of enzyme required to convert of 1 pmol of [³H]ATP in AMP-Ligase complex in 15 minutes at 25°C.

Unit Assay Conditions: 50 mM Tris-HCl (pH 7.8), 10 mM MgCl₂, 10 mM DTT, 0.2 μM [³H]ATP.

Quality Control: Free of contaminating single-stranded DNA exonuclease, endonuclease, ribonuclease and phosphatase activities.

T4 Polynucleotide Kinase

Isolated from *E.coli* strain that carries the cloned *Polynucleotide Kinase* gene from bacteriophage T4

Description: T4 polynucleotide Kinase catalyzes the transfer and exchange of P_i from the γ position of ATP to the 5' hydroxyl terminus of polynucleotides (double- and single-stranded DNA and RNA) and nucleoside 3'-monophosphates. The enzyme also catalyzes the removal of 3'- phosphoryl groups from 3'- phosphoryl polynucleotides, deoxynucleoside 3'- monophosphates and deoxynucleoside 3'-diphosphates.

Applications:

- end-labeling DNA or RNA for probes and DNA sequencing;
- addition of 5'- phosphates to oligonucleotides to allow subsequent ligation;
- removal of 3'- phosphoryl groups.

Concentration: 10 000 units/ml.

Reagents Supplied with Enzyme:

10 × T4 Polynucleotide Kinase Buffer.

E311 **500 u**

E312 **2500 u**

Storage conditions: 10 mM Tris-HCl (pH 7.5); 50 mM NaCl; 0.1 mM EDTA; 1 mM DTT; 50% glycerol.

Store at -20°C.

Unit Definition One unit is the amount of enzyme catalyzing the incorporation of 1 nmol of acid-insoluble [³²P] in 30 minutes at 37°C.

Unit Assay Conditions: 1 × T4 Polynucleotide Kinase Buffer, 66 μM [γ- ³²P] ATP (5 × 10⁶ cpm/μmol) and 0.26 mM 5'- hydroxyl-terminated salmon sperm DNA.

Quality Control: Free of exonuclease, phosphatase, endonuclease and RNase activities.

T4 DNA Ligase

Isolated from *E.coli* strain that carries the cloned *DNA Ligase* gene from bacteriophage T4

Description: T4 DNA Ligase catalyzes the formation of a phosphodiester bond between juxtaposed 5' phosphate and 3' hydroxyl termini in duplex DNA or RNA. This enzyme will join blunt-end and cohesive end termini as well as repair single stranded nicks in duplex DNA, RNA, or DNA/RNA hybrids.

Concentration: 50 000 - 200 000 units/ml

Applications:

- cloning of restriction fragments;
- joining linkers and adapters to blunt-ended DNA;

Reagents Supplied with Enzyme: 10 × T4 DNA Ligase Buffer.

E319 **10000 u**

E320 **50000 u**

Storage conditions: 10 mM Tris-HCl (pH 7.5); 50 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 50% glycerol. **Store at -20°C.**

Unit Definition: One unit is defined as the amount of enzyme required to give 50% ligation of Hind III fragments of Lambda DNA (5' DNA termini concentration of 0.12 μM [300 μg/ml]) in 20 μl of 1 × T4 DNA Ligase Reaction Buffer in 30 minutes **at 16°C.**

Quality Control: Free of contaminating exonuclease and endonuclease.

Notes on Use: ATP is an essential cofactor for the reaction.

T4 DNA Ligase Concentrated

Isolated from *E.coli* strain that carries the cloned *DNA Ligase* gene from bacteriophage T4

Description: T4 DNA Ligase catalyzes the formation of a phosphodiester bond between juxtaposed 5' phosphate and 3' hydroxyl termini in duplex DNA or RNA. This enzyme will join blunt-end and cohesive end termini as well as repair single stranded nicks in duplex DNA, RNA, or DNA/RNA hybrids.

Concentration: 2 000 000 units/ml.

Applications:

- cloning of restriction fragments;
- joining linkers and adapters to blunt-ended DNA.

Reagents Supplied with Enzyme: 10 × T4 DNA Ligase Buffer.

E329 **50000 u**

E330 **250000 u**

Storage conditions: 10 mM Tris-HCl (pH 7.5); 50 mM NaCl; 0.1 mM EDTA; 7 mM 2-mercaptoethanol; 50% glycerol. **Store at -20°C.**

Unit Definition: One unit is defined as the amount of enzyme required to give 50% ligation of Hind III fragments of Lambda DNA (5' DNA termini concentration of 0.12 μM [300 μg/ml]) in 20 μl of 1 × T4 DNA Ligase Reaction Buffer in 30 minutes **at 16°C.**

Quality Control: Free of contaminating exonuclease and endonuclease.

Notes on Use: ATP is an essential cofactor for the reaction.

Tli-Inorganic Pyrophosphatase

Isolated from *E.coli* strain that carries the cloned Inorganic pyrophosphatase gene from *Thermococcus litoralis*.

E315 100 u
E316 500 u

Description: Tli-Inorganic Pyrophosphatase catalyzes the hydrolysis of inorganic pyrophosphate to form orthophosphate:

$P_2O_7^{4-} + H_2O \rightarrow 2HPO_4^{2-}$. Enzyme is extremely thermostable.

Concentration: 1 000 units/ml.

Reagents Supplied with Enzyme:

10 × Inorganic Pyrophosphatase Buffer.

Storage Conditions: 20 mM Tris-HCl (pH 7.6); 100 mM KCl; 0.1 mM EDTA; 0.2% Tween-20; 50% glycerol. **Store at -20°C.**

Unit Definition: One unit is the amount of enzyme that will generate 40 nmoles of phosphate per minute from pyrophosphate under standard reaction conditions (a 10 minute reaction at 75°C in 50 mM Tris-HCl (pH 8.5), 1 mM MgCl₂, 0.32 mM PPi, reaction volume of 0.5 ml).

Quality Control: Free of exonuclease and endonuclease activities.

Uracil-DNA Glycosylase (UDG)

Isolated from a recombinant source

E335 1000 u
E336 5000 u

Description: *E.coli* uracil-DNA glycosylase (UDG) catalyzes the release of free uracil from uracil-containing DNA. UDG efficiently hydrolyzes uracil from single-stranded or double-stranded DNA, but not from oligomers (6 or fewer bases).

Applications:

- site-directed mutagenesis;

- as a probe for protein-DNA interaction studies.

Treatment of 1 µg of uracil-containing DNA with 1 unit of UDG for 10 minutes at 37°C renders the DNA incapable of being copied by DNA polymerase. The enzyme can be 95% heat killed by incubation at 95°C for 10 minutes.

Concentration: 20 000 – 50 000 units/ml.

Reagents Supplied with Enzyme: 10 × UDG Buffer .

Storage Conditions: 10 mM Tris-HCl (pH 7.6); 50 mM NaCl; 1 mM EDTA; 1 mM DTT; 50% glycerol. **Store at -20°C.**

Unit Definition: One unit is defined as the amount of enzyme that catalyzes the release of 60 pmol of uracil per minute from double-stranded, uracil-containing DNA. Activity is measured by release of [³H]-uracil in a 50 µl reaction containing 0.2 µg DNA (10⁴ – 10⁵ cpm/µg) in 30 minutes at 37°C.

Notes on Use: UDG is active over a broad pH range with an optimum at pH 8.0, does not require divalent cation, and is inhibited by high ionic strength (>200 mM).

DNA Ladders

1 Kb DNA Ladder

M11 50 mkg
M12 250 mkg

Description: The 1 Kb DNA Ladder has a number of proprietary plasmids which are digested to completion with appropriate restriction enzymes to yield 13 bands suitable for use as molecular weight standards for agarose gel electrophoresis.

| Base Pairs | DNA Mass | Base Pairs | DNA Mass | Base Pairs | DNA Mass |
|------------|----------|-------------|----------|------------|----------|
| 10 000 | 60 ng | 3000 | 200 ng | 750 | 60 ng |
| 8000 | 60 ng | 2500 | 70 ng | 500 | 30 ng |
| 6000 | 60 ng | 2000 | 60 ng | 250 | 20 ng |
| 5000 | 60 ng | 1500 | 50 ng | | |
| 4000 | 60 ng | 1000 | 210 ng | | |

Concentration: 200 mkg/ml

Storage buffer: 10 mM Tris-HCl (pH 8.0), 1.0 mM EDTA, 50 mM NaCl. **Store at -20°C.**

Usage recommendation: The 1 Kb DNA Ladder was not designed for precise quantification of DNA mass, but can be used for approximating the mass of DNA in comparably intense samples of similar size.

We recommend loading 1 mkg of the DNA Ladder per line.

50 Kb DNA Ladder

M29 50 mkg
M30 250 mkg

Description: The 50 Kb DNA Ladder has a number of proprietary plasmids and phages DNA which are digested to completion with appropriate restriction enzymes to yield 17 bands suitable for use as molecular weight standards for agarose pulsed-field gel electrophoresis.

| Base Pairs | DNA Mass | Base Pairs | DNA Mass | Base Pairs | DNA Mass |
|------------|----------|-------------|----------|-------------|----------|
| 48502 | 70 | 6000 | 45 | 1500 | 35 |
| 39936 | 70 | 5000 | 45 | 1000 | 150 |
| 24730 | 70 | 4000 | 45 | 750 | 45 |
| 15206 | 70 | 3000 | 140 | 500 | 20 |
| 10000 | 45 | 2500 | 50 | 250 | 10 |
| 8000 | 45 | 2000 | 45 | | |

Concentration: 200 mkg/ml

Storage buffer: 10 mM Tris-HCl (pH 8.0), 1.0 mM EDTA, 50 mM NaCl. **Store at -20°C.**

Usage recommendation: The 50 Kb DNA Ladder was not designed for precise quantification of DNA mass, but can be used for approximating the mass of DNA in comparably intense samples of similar size.

We recommend loading 0,5-1 mkg of the DNA Ladder per line.

100 bp DNA Ladder

M15 50 mkg
M16 250 mkg

Description: The 100 bp DNA Ladder has a number of proprietary plasmids which are digested to completion with appropriate restriction enzymes to yield 10 bands suitable for use as molecular weight standards for agarose gel electrophoresis.

| Base Pairs | DNA Mass | Base Pairs | DNA Mass | Base Pairs | DNA Mass |
|------------|----------|------------|----------|------------|----------|
| 1 000 | 190 ng | 600 | 110 ng | 200 | 40 ng |
| 900 | 170 ng | 500 | 130 ng | 100 | 40 ng |
| 800 | 150 ng | 400 | 80 ng | | |
| 700 | 90 ng | 300 | 40 ng | | |

Concentration: 200 mkg/ml

Storage buffer: 10 mM Tris-HCl (pH 8.0), 1.0 mM EDTA, 50 mM NaCl. **Store at -20°C.**

Usage recommendation: The 100 bp DNA Ladder was not designed for precise quantification of DNA mass, but can be used for approximating the mass of DNA in comparably intense samples of similar size.

We recommend loading 1 mkg of the DNA Ladder per line.

100 bp + 50 bp DNA Ladder

M33 50 mkg
M34 250 mkg

Description: The 100 bp + 50bp DNA Ladder has a number of proprietary plasmids which are digested to completion with appropriate restriction enzymes to yield 11 bands suitable for use as molecular weight standards for agarose gel electrophoresis.

| Base Pairs | DNA Mass | Base Pairs | DNA Mass | Base Pairs | DNA Mass |
|------------|----------|------------|----------|------------|----------|
| 1 000 | 160 ng | 600 | 100 ng | 200 | 30 ng |
| 900 | 140 ng | 500 | 170 ng | 100 | 30 ng |
| 800 | 120 ng | 400 | 70 ng | 50 | 20 ng |
| 700 | 110 ng | 300 | 50 ng | | |

Concentration: 200 mkg/ml

Storage buffer: 10 mM Tris-HCl (pH 8.0), 1.0 mM EDTA, 50 mM NaCl. **Store at -20°C.**

Usage recommendation: The 100 bp + 50 bp DNA Ladder was not designed for precise quantification of DNA mass, but can be used for approximating the mass of DNA in comparably intense samples of similar size.

We recommend loading 1 mkg of the DNA Ladder per line.

100 bp + 1.5 Kb DNA Ladder

M23 50 mkg
M24 250 mkg

Description: The 100 bp + 1.5 Kb DNA Ladder has a number of proprietary plasmids which are digested to completion with appropriate restriction enzymes to yield 11 bands suitable for use as molecular weight standards for agarose gel electrophoresis.

| Base Pairs | DNA Mass | Base Pairs | DNA Mass | Base Pairs | DNA Mass |
|-------------|----------|------------|----------|------------|----------|
| 1500 | 150 ng | 700 | 70 ng | 300 | 30 ng |
| 1000 | 200 ng | 600 | 80 ng | 200 | 20 ng |
| 900 | 120 ng | 500 | 150 ng | 100 | 20 ng |
| 800 | 110 ng | 400 | 50 ng | | |

Concentration: 200 mkg/ml

Storage buffer: 10 mM Tris-HCl (pH 8.0), 1.0 mM EDTA, 50 mM NaCl. **Store at -20°C.**

Usage recommendation: The 100 bp + 1.5 Kb DNA Ladder was not designed for precise quantification of DNA mass, but can be used for approximating the mass of DNA in comparably intense samples of similar size.

We recommend loading 1 mkg of the DNA Ladder per line.

100 bp + 1.5 Kb + 3 Kb DNA Ladder

M27 50 mkg
M28 250 mkg

Description: The 100 bp + 1.5 Kb + 3 Kb DNA Ladder has a number of proprietary plasmids which are digested to completion with appropriate restriction enzymes to yield 12 bands suitable for use as molecular weight standards for agarose gel electrophoresis.

| Base Pairs | DNA Mass | Base Pairs | DNA Mass | Base Pairs | DNA Mass |
|-------------|----------|------------|----------|------------|----------|
| 3000 | 70 ng | 800 | 120 ng | 400 | 60 ng |
| 1500 | 60 ng | 700 | 80 ng | 300 | 40 ng |
| 1000 | 140 ng | 600 | 90 ng | 200 | 30 ng |
| 900 | 130 ng | 500 | 120 ng | 100 | 30 ng |

Concentration: 200 mkg/ml

Storage buffer: 10 mM Tris-HCl (pH 8.0), 1.0 mM EDTA, 50 mM NaCl. **Store at -20°C.**

Usage recommendation: The 100 bp + 1.5 Kb + 3 Kb DNA Ladder was not designed for precise quantification of DNA mass, but can be used for approximating the mass of DNA in comparably intense samples of similar size.

We recommend loading 1 mkg of the DNA Ladder per line.

100 bp + 2 Kb + 3 Kb DNA Ladder

M25 50 mkg
M26 250 mkg

Description: The 100 bp + 2 Kb + 3 Kb DNA Ladder has a number of proprietary plasmids which are digested to completion with appropriate restriction enzymes to yield 12 bands suitable for use as molecular weight standards for agarose gel electrophoresis.

| Base Pairs | DNA Mass | Base Pairs | DNA Mass | Base Pairs | DNA Mass |
|-------------|----------|------------|----------|------------|----------|
| 3000 | 120 ng | 800 | 110 ng | 400 | 50 ng |
| 2000 | 120 ng | 700 | 70 ng | 300 | 30 ng |
| 1000 | 160 ng | 600 | 80 ng | 200 | 20 ng |
| 900 | 120 ng | 500 | 100 ng | 100 | 20 ng |

Concentration: 200 mkg/ml

Storage buffer: 10 mM Tris-HCl (pH 8.0), 1.0 mM EDTA, 50 mM NaCl. **Store at -20°C.**

Usage recommendation: The 100 bp + 2 Kb + 3 Kb DNA Ladder was not designed for precise quantification of DNA mass, but can be used for approximating the mass of DNA in comparably intense samples of similar size.

We recommend loading 1 mkg of the DNA Ladder per line.

Lambda DNA - Hind III Digest

M01 100 mkg
M02 500 mkg

Description: The Hind III digest of lambda DNA yields 8 fragments suitable for use as molecular weight standards for agarose gel electrophoresis.

| Length of the DNA fragments (Base Pairs) | | | |
|--|------|------|-----|
| 23130 | 6557 | 2322 | 564 |
| 9416 | 4361 | 2027 | 125 |

Concentration: 500 mkg/ml

Storage buffer (TE): 10 mM Tris-HCl (pH 8.0), 1.0 mM EDTA. **Store at -20°C.**

Usage recommendation: Dilute in TE buffer (10 – 20 times) before loading.

We recommend loading 1 mkg of the DNA Ladder per line.

Lambda DNA - BssT1 I (Sty I) Digest

M05 100 mkg
M06 500 mkg

Description: The BssT1 I (Sty I) digest of lambda DNA yields 11 fragments suitable for use as molecular weight standards for agarose gel electrophoresis.

| Length of the DNA fragments (Base Pairs) | | | |
|--|------|------|-----|
| 19329 | 4254 | 1882 | 421 |
| 7743 | 3472 | 1489 | 74 |
| 6223 | 2690 | 925 | |

Concentration: 500 mkg/ml

Storage buffer (TE): 10 mM Tris-HCl (pH 8.0), 1.0 mM EDTA. **Store at -20°C.**

Usage recommendation: Dilute in TE buffer (10 – 20 times) before loading.

We recommend loading 1 mkg of the DNA Ladder per line.

Lambda DNA - Bme18 I (Ava II) Digest

M03 100 mkg
M04 500 mkg

Description: The Bme18 I (Ava II) digest of lambda DNA yields 36 fragments suitable for use as molecular weight standards for agarose gel electrophoresis.

| Length of the DNA fragments (Base Pairs) | | | | | |
|--|------|-----|-----|-----|----|
| 8126 | 2134 | 985 | 511 | 272 | 67 |
| 6555 | 2005 | 974 | 433 | 242 | 45 |
| 6442 | 1951 | 894 | 398 | 215 | 42 |
| 3676 | 1612 | 597 | 345 | 151 | 32 |
| 2605 | 1420 | 590 | 310 | 88 | 28 |
| 2555 | 1284 | 513 | 308 | 73 | 23 |

Concentration: 500 mkg/ml

Storage buffer (TE): 10 mM Tris-HCl (pH 8.0), 1.0 mM EDTA. **Store at -20°C.**

Usage recommendation: Dilute in TE buffer (10 – 20 times) before loading.

We recommend loading 1 mkg of the DNA Ladder per line.

Lambda DNA - Bgl I Digest

M17 100 mkg
M18 500 mkg

Description: The Bgl I digest of lambda DNA yields 30 fragments suitable for use as molecular weight standards for agarose gel electrophoresis.

| Length of the DNA fragments (Base Pairs) | | | | | |
|--|------|------|-----|-----|-----|
| 16179 | 1650 | 1138 | 562 | 366 | 126 |
| 9649 | 1446 | 790 | 499 | 267 | 115 |
| 3009 | 1441 | 773 | 489 | 210 | 91 |
| 2481 | 1249 | 669 | 447 | 186 | 51 |
| 2256 | 1203 | 621 | 404 | 126 | 9 |

Concentration: 500 mkg/ml

Storage buffer (TE): 10 mM Tris-HCl (pH 8.0), 1.0 mM EDTA. **Store at -20°C.**

Usage recommendation: Dilute in TE buffer (10 – 20 times) before loading.

We recommend loading 1 mkg of the DNA Ladder per line.

Lambda DNA – BstE II Digest

M09 100 mkg
M10 500 mkg

Description: The BstEII digest of lambda DNA yields 14 fragments suitable for use as molecular weight standards for agarose gel electrophoresis.

| Length of the DNA fragments (Base Pairs) | | | | | |
|--|------|------|------|-----|--|
| 8454 | 5686 | 3675 | 1371 | 224 | |
| 7242 | 4822 | 2323 | 1264 | 117 | |
| 6369 | 4324 | 1929 | 702 | | |

Concentration: 500 mkg/ml

Storage buffer (TE): 10 mM Tris-HCl (pH 8.0), 1.0 mM EDTA. **Store at -20°C.**

Usage recommendation: Dilute in TE buffer (10 – 20 times) before loading.

We recommend loading 1 mkg of the DNA Ladder per line.

pBR322 DNA - BsuR I (HaeIII) Digest

M21 50 mkg
M22 250 mkg

Description: The BsuR I digest of pBR322 DNA yields 22 fragments suitable for use as molecular weight standards for agarose gel electrophoresis.

| Length of the DNA fragments (Base Pairs) | | | | | |
|--|-----|-----|-----|----|----|
| 587 | 434 | 192 | 104 | 57 | 11 |
| 540 | 267 | 184 | 89 | 51 | 8 |
| 502 | 234 | 124 | 80 | 21 | |
| 458 | 213 | 123 | 64 | 18 | |

Concentration: 200 mkg/ml

Storage buffer: 10 mM Tris-HCl (pH 8.0), 1.0 mM EDTA, 50 mM NaCl. **Store at -20°C.**

Usage recommendation: We recommend loading 1 mkg of the DNA Ladder per line.

pBR322 DNA - Alu I Digest

M19 50 mkg
M20 250 mkg

Description: The Alu I digest of pBR322 DNA yields 9 fragments suitable for use as molecular weight standards for agarose gel electrophoresis.

| Length of the DNA fragments (Base Pairs) | | |
|--|-----|-----|
| 908 | 521 | 257 |
| 659 | 403 | 226 |
| 656 | 281 | 100 |

Concentration: 200 mkg/ml

Storage buffer: 10 mM Tris-HCl (pH 8.0), 1.0 mM EDTA, 50 mM NaCl. **Store at -20°C.**

Usage recommendation: We recommend loading 1 mkg of the DNA Ladder per line.

pUC19 DNA – Kzo9 I (Sau3AI) Digest

M13 50 mkg
M14 250 mkg

Description: The Kzo9 I digest of pUC19 DNA yields 15 fragments suitable for use as molecular weight standards for agarose gel electrophoresis.

| Length of the DNA fragments (Base Pairs) | | | |
|--|-----|----|----|
| 955 | 141 | 46 | 12 |
| 585 | 105 | 36 | 11 |
| 341 | 78 | 18 | 8 |
| 258 | 75 | 17 | |

Concentration: 200 mkg/ml

Storage buffer: 10 mM Tris-HCl (pH 8.0), 1.0 mM EDTA, 50 mM NaCl. **Store at -20°C.**

Usage recommendation: We recommend loading 1 mkg of the DNA Ladder per line.

pUC19 DNA - Msp I Digest

M07 50 mkg
M08 250 mkg

Description: The Msp I digest of pUC19 DNA yields 12 bands suitable for use as molecular weight standards for agarose gel electrophoresis.

| Length of the DNA fragments (Base Pairs) | | | |
|--|-----|-----|----|
| 501 | 331 | 147 | 67 |
| 489 | 242 | 111 | 34 |
| 404 | 190 | 110 | 26 |

Concentration: 200 mkg/ml

Storage buffer: 10 mM Tris-HCl (pH 8.0), 1.0 mM EDTA, 50 mM NaCl. **Store at -20°C.**

Usage recommendation: We recommend loading 1 mkg of the DNA Ladder per line.

Plasmid and Phage DNAs

Lambda DNA (dam⁻, dcm⁻)

Isolated from bacteriophage lambda (c1857ind 1 Sam 7) obtained from heat inducible lysogen *E.coli* strain (dam⁻, dcm⁻)

D10 500 mkg

Description: Duplex DNA is 48502 base pairs in length. The molecular weight is 31.5x10⁶ daltons.

Concentration: 500 mkg/ml

Storage Buffer: 10 mM Tris-HCl (pH8.0) and 1 mM EDTA. **Store at -20°C.**

Lambda DNA

Isolated from bacteriophage lambda (c1857ind 1 Sam 7) obtained from heat inducible lysogen *E.coli* strain (dam⁺, dcm⁺).

D11 500 mkg

Description: Duplex DNA is 48502 base pairs in length. The molecular weight is 31.5x10⁶ daltons.

Concentration: 500 mkg/ml

Storage Buffer: 10 mM Tris-HCl (pH 8.0) and 1 mM EDTA. **Store at -20°C.**

T7 Phage DNA

Isolated from T7 phage obtained from infected *E.coli* strain

D02 500 mkg

Description: Duplex DNA is 39936 base pairs in length. The molecular weight is 26x10⁶ daltons.

Concentration: 500 mkg/ml.

Storage Buffer: 10 mM Tris-HCl (pH 8.0) and 1 mM EDTA. **Store at -20°C.**

pBR322 DNA

Isolated from *E.coli* XL1-Blue (dam⁺, dcm⁺).

D03 50 mkg

D04 250 mkg

Description: pBR322 is commonly used plasmid cloning vector in *E.coli*. The molecule is a double-stranded circle, 4361 base pairs in length. pBR322 contains the genes for resistance to ampicillin and tetracycline, and may be amplified with chloramphenicol. The molecular weight of pBR322 is 2.83x10⁶ daltons.

Concentration: 200 mkg/ml

Storage Buffer: 10 mM Tris-HCl (pH8.0) and 1 mM EDTA. **Store at -20°C.**

pHspAI2/GsaI DNA

pHspAI2 is isolated from *E.coli* (dam⁺,dcm⁺) by a standard plasmid purification procedure.

D09 10 mkg

Description: DNA pHspAI2/GsaI is a linearized plasmid pHspAI2, which includes a gene of DNA-methyltransferase M.HspAI (recognition sequence 5'-GCGC-3') and contains a unique GsaI canonical site 5'-G(5mC)G(5mC)-3'/3'-(5mC)G(5mC)G-5'.

Concentration: 200 mkg/ml

Storage Buffer: 10 mM Tris-HCl (pH8.0) and 1 mM EDTA. **Store at -20°C.**

pUC19 DNA

isolated from *E.coli* XL1-Blue (dam⁺, dcm⁺)

D05 50 mkg

D06 250 mkg

Description: pUC19 is commonly used plasmid cloning vector in *E.coli*. The molecule is a small double-stranded circle, 2686 base pairs in length, and has a high copy number. pUC19 carries a 54 base-pair multiple cloning site polylinker that contains unique sites for 13 different hexanucleotide-specific restriction endonucleases. The molecular weight of pUC19 is 1.75x10⁶ daltons

Concentration: 200 mkg/ml

Storage Buffer: 10 mM Tris-HCl (pH8.0) and 1 mM EDTA. **Store at -20°C.**

Human Genomic DNA

Genomic DNA HeLa

D07 10 mkg

Isolated from HeLa cervix adenocarcinoma cells.

Description: Human female HeLa genomic DNA.

Cells were grown to confluency in DMEM – 90%, plus fetal bovine serum – 10%.

Genomic DNA was isolated by a standard protocol [1].

Concentration: 200 mkg /ml.

Quality Control: Purified free of contaminating proteins and RNA

Storage Buffer: 10 mM Tris-HCl (pH 7.8), 1 mM EDTA. **Store at -20°C.**

Note: Avoid multiple freeze/thaw cycles.

1. Sambrook J. and Russell D. (2001) *Molecular Cloning: A Laboratory Manual*, (3rd ed.), pp. 6.4-6.12. Cold Spring Harbor: Cold Spring Harbor Laboratory Press.

Genomic DNA Jurkat

D08 10 mkg

Isolated from Jurkat acute T-cell leukemia cells.

Description: Human male Jurkat genomic DNA.

Cells were grown to confluency in RPMI-1640 – 90%, plus fetal bovine serum – 10%.

Genomic DNA was isolated by a standard protocol [1].

Concentration: 200 mkg /ml.

Quality Control: Purified free of contaminating proteins and RNA

Storage Buffer: 10 mM Tris-HCl (pH 7.8), 1 mM EDTA. **Store at -20°C.**

Note: Avoid multiple freeze/thaw cycles.

1. Sambrook J. and Russell D. (2001) *Molecular Cloning: A Laboratory Manual*, (3rd ed.), pp. 6.4-6.12. Cold Spring Harbor: Cold Spring Harbor Laboratory Press.

Genomic DNA L-68

D014 10 mkg

Isolated from L-68 fibroblast cells from lung.

Description: Human embryo L-68 genomic DNA.

Cells were grown to confluency in DMEM – 90%, plus fetal bovine serum – 10%.

Genomic DNA was isolated by a standard protocol [1].

Concentration: 200 mkg /ml.

Quality Control: Purified free of contaminating proteins and RNA

Storage Buffer: 10 mM Tris-HCl (pH 7.8), 1 mM EDTA. **Store at -20°C.**

Note: Avoid multiple freeze/thaw cycles.

1. Sambrook J. and Russell D. (2001) *Molecular Cloning: A Laboratory Manual*, (3rd ed.), pp. 6.4-6.12. Cold Spring Harbor: Cold Spring Harbor Laboratory Press.

Genomic DNA Raji

D13 10 mkg

Isolated from Raji Burkitt's lymphoma cells.

Description: Human male Raji genomic DNA.

Cells were grown to confluency in RPMI-1640 – 90%, plus fetal bovine serum – 10%.

Genomic DNA was isolated by a standard protocol [1].

Concentration: 200 mkg /ml.

Quality Control: Purified free of contaminating proteins and RNA

Storage Buffer: 10 mM Tris-HCl (pH 7.8), 1 mM EDTA. **Store at -20°C.**

Note: Avoid multiple freeze/thaw cycles.

1. Sambrook J. and Russell D. (2001) *Molecular Cloning: A Laboratory Manual*, (3rd ed.), pp. 6.4-6.12. Cold Spring Harbor: Cold Spring Harbor Laboratory Press.

Genomic DNA U-937

D12 10 mkg

Isolated from U-937 histiocytic lymphoma cells.

Description: Human male U-937 genomic DNA.

Cells were grown to confluency in RPMI-1640 – 90%, plus fetal bovine serum – 10%.

Genomic DNA was isolated by a standard protocol [1].

Concentration: 200 mkg /ml.

Quality Control: Purified free of contaminating proteins and RNA


Storage Buffer: 10 mM Tris-HCl (pH 7.8), 1 mM EDTA. **Store at -20°C.**

Note: Avoid multiple freeze/thaw cycles.






1. Sambrook J. and Russell D. (2001) *Molecular Cloning: A Laboratory Manual*, (3rd ed.), pp. 6.4-6.12. Cold Spring Harbor: Cold Spring Harbor Laboratory Press.





dNTPs(enzymatic)

| | | | |
|--|---------------------------------|--------------|-----------------------|
| dATP | 40 mM water solution | N011 N012 | 5 µmoles 25 µmoles |
| <p>Description: 2'-deoxyadenosine 5'-triphosphate sodium salt Concentration: 40 mM Method of production: enzymatic synthesis. dATP are Qualified free of nucleases and phosphatases. 15 Kb and longer DNA fragments may be amplified in PCR with this dNTP(enzymatic). Store at -20°C.</p> | | | |
| dCTP | 40 mM water solution | N013 N014 | 5 µmoles 25 µmoles |
| <p>Description: 2'-deoxycytidine 5'-triphosphate sodium salt Concentration: 40 mM Method of production: enzymatic synthesis. dCTP are Qualified free of nucleases and phosphatases. 15 Kb and longer DNA fragments may be amplified in PCR with this dNTP(enzymatic). Store at -20°C.</p> | | | |
| dGTP | 40 mM water solution | N015 N016 | 5 µmoles 25 µmoles |
| <p>Description: 2'-deoxyguanosine 5'-triphosphate sodium salt Concentration: 40 mM Method of production: enzymatic synthesis. dGTP are Qualified free of nucleases and phosphatases. 15 Kb and longer DNA fragments may be amplified in PCR with this dNTP(enzymatic). Store at -20°C.</p> | | | |
| dTTP | 40 mM water solution | N017 N018 | 5 µmoles 25 µmoles |
| <p>Description: 2'-deoxythymidine 5'-triphosphate sodium salt Concentration: 40 mM Method of production: enzymatic synthesis. dTTP are Qualified free of nucleases and phosphatases. 15 Kb and longer DNA fragments may be amplified in PCR with this dNTP(enzymatic). Store at -20°C.</p> | | | |
| dUTP | 100 mM water solution | N031 N032 | 5 µmoles 25 µmoles |
| <p>Description: 2'-deoxyuridine 5'-triphosphate sodium salt Concentration: 100 mM Method of production: enzymatic synthesis. dUTP are Qualified free of nucleases and phosphatases. 15 Kb and longer DNA fragments may be amplified in PCR with this dNTP(enzymatic). Store at -20°C.</p> | | | |
| dNTP Mix | 0,5 µmol of each Water solution | N024 | 1 ml |
| <p>Description: Deoxynucleotide solution mix is an equimolar solution of ultrapure dATP, dCTP, dGTP and dTTP, each at a final concentration of 0,5 mM. Quantity: 1.0 ml contains 0,5 µmol of each dATP, dCTP, dGTP and dTTP. Application: For use in molecular biology and biotechnology, including PCR, long-PCR, and for labeling and sequencing DNA. PCR-diagnostics for medicine. Quality: Molecular Biology Grade. Functional test: The product has been tested for receiving products with length more than 15,000 base pair nucleotides with T7 phage DNA as a template. Store at -20°C.</p> | | | |

| | | | | |
|---|---|---|-------------|-------------------------------------|
| dNTP Mix | | 2,5 μmol of each 40 mM water solution | N026 | 1 ml |
| <p>Description: Deoxynucleotide solution mix is an equimolar solution of ultrapure dATP, dCTP, dGTP and dTTP, each at a final concentration of 2,5 mM.</p> <p>Quantity: 1.0 ml contains 2,5 μmol of each dATP, dCTP, dGTP and dTTP.</p> <p>Application: For use in molecular biology and biotechnology, including PCR, long-PCR, and for labeling and sequencing DNA. PCR-diagnostics for medicine.</p> <p>Quality: Molecular Biology Grade.</p> <p>Functional test: The product has been tested for receiving products with length more than 15,000 base pair nucleotides with T7 phage DNA as a template.</p> <p>Store at -20°C.</p> | | | | |
| dNTP Mix | | 4 μmol of each 40 mM water solution | N027 | 1 ml |
| <p>Description: Deoxynucleotide solution mix is an equimolar solution of ultrapure dATP, dCTP, dGTP and dTTP, each at a final concentration of 4 mM.</p> <p>Quantity: 1.0 ml contains 4 μmol of each dATP, dCTP, dGTP and dTTP.</p> <p>Application: For use in molecular biology and biotechnology, including PCR, long-PCR, and for labeling and sequencing DNA. PCR-diagnostics for medicine.</p> <p>Quality: Molecular Biology Grade.</p> <p>Functional test: The product has been tested for receiving products with length more than 15,000 base pair nucleotides with T7 phage DNA as a template.</p> <p>Store at -20°C.</p> | | | | |
| dNTP Mix | | 10 μmol of each 40 mM water solution | N025 | 1 ml |
| <p>Description: Deoxynucleotide solution mix is an equimolar solution of ultrapure dATP, dCTP, dGTP and dTTP, each at a final concentration of 10 mM.</p> <p>Quantity: 1.0 ml contains 10 μmol of each dATP, dCTP, dGTP and dTTP.</p> <p>Application: For use in molecular biology and biotechnology, including PCR, long-PCR, and for labeling and sequencing DNA. PCR-diagnostics for medicine.</p> <p>Quality: Molecular Biology Grade.</p> <p>Functional test: The product has been tested for receiving products with length more than 15,000 base pair nucleotides with T7 phage DNA as a template.</p> <p>Store at -20°C.</p> | | | | |
| dNTP solution set (enzymatic) |  | 40 mM of each Water solution. | N030 | 4 x 25 μmoles |
| <p>Description: Deoxynucleotide solution set: Four separate solutions of ultrapure nucleotides (dATP, dCTP, dGTP and dTTP). Each nucleotide is supplied at a concentration of 40 mM aqueous solution .</p> <p>Application: For use in molecular biology and biotechnology, including PCR, long-PCR, and for labeling and sequencing DNA. PCR-diagnostics for medicine.</p> <p>Quality: Molecular Biology Grade.</p> <p>Functional test: The product has been tested for receiving products with length more than 15,000 base pairs with T7 phage DNA as a template.</p> <p>Store at -20°C.</p> | | | | |

dNTPs(chemical)

| | | | | |
|---|---|--|--------------|-----------------------|
| dATP |  | 100 mM water solution. | N001 N002 | 5 µmoles 25 µmoles |
| <p>Description: 2'-deoxyadenosine 5'-triphosphate lithium salt Concentration: 100 mM Method of production: chemical synthesis. Purity of dATP is 98-99% according to HPLC analysis (Column ProntoSIL-120-5-C18AQ). The product is free of DNases, RNases and phosphatases. Coefficient of extinction: 15200 M⁻¹ X cm⁻¹ at λ_{max}=259 nm Functional test: The product has been tested in PCR for 0.1-15 kb fragments in length. Store at -20°C.</p> | | | | |
| dCTP |  | 100 mM water solution. | N003 N004 | 5 µmoles 25 µmoles |
| <p>Description: 2'-deoxyadenosine 5'-triphosphate lithium salt Concentration: 100 mM Method of production: chemical synthesis. Purity of dATP is 98-99% according to HPLC analysis (Column ProntoSIL-120-5-C18AQ). The product is free of DNases, RNases and phosphatases. Coefficient of extinction: 9300 M⁻¹ X cm⁻¹ at λ_{max}=271 nm Functional test: The product has been tested in PCR for 0.1-15 kb fragments in length. Store at -20°C.</p> | | | | |
| dGTP |  | 100 mM water solution. | N005 N006 | 5 µmoles 25 µmoles |
| <p>Description: 2'-deoxyadenosine 5'-triphosphate lithium salt Concentration: 100 mM Method of production: chemical synthesis. Purity of dATP is 98-99% according to HPLC analysis (Column ProntoSIL-120-5-C18AQ). The product is free of DNases, RNases and phosphatases. Coefficient of extinction: 13700 M⁻¹ X cm⁻¹ at λ_{max}=253 nm Functional test: The product has been tested in PCR for 0.1-15 kb fragments in length. Store at -20°C.</p> | | | | |
| dTTP |  | 100 mM water solution. | N007 N008 | 5 µmoles 25 µmoles |
| <p>Description: 2'-deoxyadenosine 5'-triphosphate lithium salt Concentration: 100 mM Method of production: chemical synthesis. Purity of dATP is 98-99% according to HPLC analysis (Column ProntoSIL-120-5-C18AQ). The product is free of DNases, RNases and phosphatases. Coefficient of extinction: 9600 M⁻¹ X cm⁻¹ at λ_{max}=267 nm Functional test: The product has been tested in PCR for 0.1-15 kb fragments in length. Store at -20°C.</p> | | | | |
| dNTP Mix |  | 0,5 µmol of each 100 mM water solution. | N020 | 1 ml |
| <p>Description: Deoxynucleotide solution mix is an equimolar solution of ultrapure dATP, dCTP, dGTP, dTTP, each at a final concentration of 0,5 mM. Concentration: 100 mM Method of production: chemical synthesis. Quantity: 1.0 ml contains 0,5 µmol of each dATP, dCTP, dGTP and dTTP. Application: For use in molecular biology and biotechnology, including PCR, long-PCR, and for labeling and sequencing DNA. PCR-diagnostics for medicine. Quality: Molecular Biology Grade. Functional test: The product has been tested for receiving products with length more than 15,000 base pair nucleotides with T7 phage DNA as a template. Store at -20°C.</p> | | | | |

| | | | | |
|--|---|--|-------------|----------------------|
| dNTP Mix |  | 2,5 µmol of each 100 mM water solution. | N021 | 1 ml |
| <p>Description: Deoxynucleotide solution mix is an equimolar solution of ultrapure dATP, dCTP, dGTP and dTTP, each at a final concentration of 2,5 mM.</p> <p>Concentration: 100 mM</p> <p>Method of production: chemical synthesis.</p> <p>Quantity: 1.0 ml contains 2,5 µmol of each dATP, dCTP, dGTP and dTTP.</p> <p>Application: For use in molecular biology and biotechnology, including PCR, long-PCR, and for labeling and sequencing DNA. PCR-diagnostics for medicine.</p> <p>Quality: Molecular Biology Grade.</p> <p>Functional test: The product has been tested for receiving products with length more than 15,000 base pair nucleotides with T7 phage DNA as a template.</p> <p>Store at -20°C.</p> | | | | |
| dNTP Mix |  | 4 µmol of each 100 mM water solution. | N023 | 1 ml |
| <p>Description: Deoxynucleotide solution mix is an equimolar solution of ultrapure dATP, dCTP, dGTP and dTTP, each at a final concentration of 4 mM.</p> <p>Concentration: 100 mM</p> <p>Method of production: chemical synthesis.</p> <p>Quantity: 1.0 ml contains 4 µmol of each dATP, dCTP, dGTP and dTTP.</p> <p>Application: For use in molecular biology and biotechnology, including PCR, long-PCR, and for labeling and sequencing DNA. PCR-diagnostics for medicine.</p> <p>Quality: Molecular Biology Grade.</p> <p>Functional test: The product has been tested for receiving products with length more than 15,000 base pair nucleotides with T7 phage DNA as a template.</p> <p>Store at -20°C.</p> | | | | |
| dNTP Mix |  | 10 µmol of each 100 mM water solution. | N022 | 1 ml |
| <p>Description: Deoxynucleotide solution mix is an equimolar solution of ultrapure dATP, dCTP, dGTP and dTTP, each at a final concentration of 10 mM.</p> <p>Concentration: 100 mM</p> <p>Method of production: chemical synthesis.</p> <p>Quantity: 1.0 ml contains 10 µmol of each dATP, dCTP, dGTP and dTTP.</p> <p>Application: For use in molecular biology and biotechnology, including PCR, long-PCR, and for labeling and sequencing DNA. PCR-diagnostics for medicine.</p> <p>Quality: Molecular Biology Grade.</p> <p>Functional test: The product has been tested for receiving products with length more than 15,000 base pair nucleotides with T7 phage DNA as a template.</p> <p>Store at -20°C.</p> | | | | |
| dNTP solution set (enzymatic) |  | 100 mM of each water solution. | N028 | 4 x 25 µmoles |
| <p>Description: Deoxynucleotide solution set: Four separate solutions of ultrapure nucleotides (dATP, dCTP, dGTP and dTTP). Each nucleotide is supplied at a concentration of 100 mM aqueous solution .</p> <p>Concentration: 100 mM</p> <p>Method of production: chemical synthesis.</p> <p>Application: For use in molecular biology and biotechnology, including PCR, long-PCR, and for labeling and sequencing DNA. PCR-diagnostics for medicine.</p> <p>Quality: Molecular Biology Grade.</p> <p>Functional test: The product has been tested for receiving products with length more than 15,000 base pairs with T7 phage DNA as a template.</p> <p>Store at -20°C.</p> | | | | |

Substrates

S-Adenosyl – L - methionine (SAM)

S005

100 mkl

Source: fermentative synthesis with subsequent purification

Description: S-Adenosyl-L-methionine (SAM) is a substrate for methyltransferases and a cofactor for some restriction endonucleases. SAM provided in 5 mM H₂SO₄ and 10% EtOH.

Formula: C₁₅H₂₃N₆O₅S

Molecular weight: 399,5

Concentration: 32 mM

Use and Storage: Store at -20°C.

The SAM solution should be thawed on ice and added into the reaction mixture immediately before the incubation. For use in a methylation reaction the SAM should be diluted to a final concentration of 80 µM.

Quality Control: the SAM is tested in a reaction of methylation of phage λ DNA with DNA-methyltransferase M.HspAI (SE #003).

Kits

DNA Quick Ligation Kit



K006

20
reactions

Description: The Quick Ligation Kit includes:

1. Quick T4 DNA Ligase (recombinant, 2000 u/μl) - 20μl
2. 2 x Quick Ligation Buffer - 200μl (132 mM Tris-HCl (pH 7.6 at 25°C); 20 mM MgCl₂, 2 mM DTT; 2 mM ATP; 15% PEG 6000).

Mix thoroughly before use.

Ligation reaction protocol (for 20 μl)

1. Mix:
 - H₂O - calculated quantity
 - 2 x Quick Ligation Buffer - 10μl
 - Add DNA
 - Add 1 μl Quick T4 DNA Ligase
2. Incubate at room temperature (+25°C) for 5 min.
3. Cool the mixture on ice, store at -20°C.

Notes on Use: Before the first use the 2 x SE-DNA Quick Ligation Buffer should be divided into small aliquots and store at -20°C.

Avoid defrosting this Buffer more than 2-3 times.

The Buffer can be stored at +4°C during 7 days.

The efficiency of ligation starts to decrease after 2 hours of incubation and is reduced by up to 75% if the reaction is allowed to go overnight at 25°C.

Before using the products of a Quick Ligation reaction for electrotransformation, it is necessary to reduce the PEG concentration.

Store at -20°C.

GLAD-PCR assay Kit



K009

200
reactions

GLAD-PCR assay Kit was developed for 200 GLAD-PCR reactions.

The following reagents are supplied with this product:

1. dH₂O - 3 mL,
2. 10x SE TMN Buffer - 150 μL,
3. DMSO - 30 μL,
4. BSA, 10 mg/mL - 80 μL,
5. GluI (10 u.a./μL) - 10 μL,
6. Universal adapter, double-stranded (10 μM) - 100 μL,
7. ATP, 10 mM - 100 μL,
8. DNA ligase T4 concentrated (2000 u.a./μL) - 10 μL,
9. 10X SE GLAD Buffer - 400 μL,
10. MgCl₂, 50 mM - 110 μL,
11. dNTP Mix, 10 mM each - 80 μL,
12. SP Taq DNA Polymerase, 5 u.a./μL - 15 μL,
13. Control DNA Raji, 9 ng/μL - 70 μL,
14. Control DNA HeLa, 9 ng/μL - 70 μL,
15. Control DNA Lambda, 9 ng/μL - 70 μL,
16. Control URB1 mix (oligonucleotide primers + TaqMan probe, 10 μM each) - 25 μL,
17. Control CEBPD mix (oligonucleotide primers + TaqMan probe, 10 μM each) - 25 μL

Additional services are available with GLAD-PCR assay Kit order:

1. Primers and TaqMan probes design (free for up to 3 RCGY sites);
2. Primes and TaqMan probes synthesis - price upon request;
3. GLAD PCR assay with synthesized oligonucleotides on selected DNA from malignant cell lines (see SE list of products) - price upon request.

Reaction requires:

1. DNA sample,
2. Genome primer and TaqMan probe designed for DNA region of interest,
3. Hybrid primer (includes constant part complimentary to universal adapter and tetranucleotide part which is complimentary to DNA at the point of GluI hydrolysis). Hybrid primers may be ordered by request (see Cat.# K010/1, K010/3, K010/5 and K010/32).

All reagents for the control experiments are included.

The time of assay including DNA hydrolysis, ligation of the adaptor and PCR is less than 4 hours.

Store at -20°C.



GLAD-PCR assay Kit was developed for 200 GLAD-PCR reactions.

The following reagents are supplied with this product:

1. dH₂O - 3 mL,
2. 10x SE TMN Buffer - 150 µL,
3. DMSO - 30 µL,
4. BSA, 10 mg/mL - 80 µL,
5. GluI (10 u.a./µL) - 10 µL,
6. Universal adapter, double-stranded (10 µM) - 100 µL,
7. ATP, 10 mM - 100 µL,
8. DNA ligase T4 concentrated (2000 u.a./µL) - 10 µL,
9. 10X SE GLAD Buffer - 400 µL,
10. MgCl₂, 50 mM - 110 µL,
11. dNTP Mix, 10 mM each - 80 µL,
12. SP Taq DNA Polymerase, 5 u.a./µL - 15 µL,
13. Control DNA Raji, 9 ng/µL - 70 µL,
14. Control DNA HeLa, 9 ng/µL - 70 µL,
15. Control DNA Lambda, 9 ng/µL - 70 µL,
16. Control URB1 mix (oligonucleotide primers + TaqMan probe, 10 µM each) - 25 µL,
17. Control CEBPD mix (oligonucleotide primers + TaqMan probe, 10 µM each) - 25 µL
18. Hybrid primers¹ (10 µM) - 170 µL.

¹ - You can choose 1 (Cat.# K010/1), 3 (Cat.# K010/3), 5 (Cat.# K010/5) of hybrid primers, or order full range of 32 hybrid primers (Cat.# K010/1).

The list of 32 hybrid primers (HP1-HP32) you can find in the Instruction manual

Additional services are available with GLAD-PCR assay Kit order:

1. Primers and TaqMan probes design (free for up to 3 RCGY sites);
2. Primes and TaqMan probes synthesis - price upon request;
3. GLAD PCR assay with synthesized oligonucleotides on selected DNA from malignant cell lines (see SE list of products) - price upon request.

Reaction requires:

1. DNA sample,
2. Genome primer and TaqMan probe designed for DNA region of interest,
3. Hybrid primer (includes constant part complimentary to universal adapter and tetranucleotide part which is complimentary to DNA at the point of GluI hydrolysis). Hybrid primers may be ordered by request. The Kits K010/1, K010/3, K010/5 and K010/32 includes 1, 3, 5 of chosen primers, or full range of 32 hybrid primers, respectively. The concentration of each hybrid primer is 10 µM, quantity - 170 µL. All reagents for the control experiments are included. The time of assay including DNA hydrolysis, ligation of the adaptor and PCR is less than 4 hours.

Store at -20°C.

Buffers composition

| Cat # | Buffer | Contents ×1 |
|-------|---|---|
| B001 | <u>B:</u> | 10 mM Tris-HCl (pH 7.6 at 25° C); 10 mM MgCl ₂ ; 1 mM DTT. |
| B002 | <u>G:</u> | 10 mM Tris-HCl (pH 7.6 at 25° C); 10 mM MgCl ₂ ; 50 mM NaCl; 1 mM DTT. |
| B003 | <u>O:</u> | 50 mM Tris-HCl (pH 7.6 at 25° C); 10 mM MgCl ₂ ; 100 mM NaCl; 1 mM DTT. |
| B004 | <u>W:</u> | 10 mM Tris-HCl (pH 8.5 at 25° C); 10 mM MgCl ₂ ; 100 mM NaCl; 1 mM DTT. |
| B005 | <u>Y:</u> | 33 mM Tris- acetate (pH 7.9 at 25° C); 10 mM magnesium acetate; 66 mM potassium acetate; 1 mM DTT. |
| B006 | <u>2W:</u> | 20 mM Tris-HCl (pH 8.5 at 25° C); 10 mM MgCl ₂ ; 200 mM NaCl; 1 mM DTT. |
| B007 | <u>K:</u> | 10 mM Tris-HCl (pH 7.6 at 25° C); 10 mM MgCl ₂ ; 100 mM KCl; 1 mM DTT. |
| B008 | <u>2K:</u> | 10 mM Tris-HCl (pH 7.6 at 25° C); 10 mM MgCl ₂ ; 200 mM KCl; 1 mM DTT. |
| B009 | <u>FaeI:</u> | 33 mM Tris- acetate (pH 8.3 at 25° C); 10 mM magnesium acetate; 66 mM potassium acetate; 1 mM DTT. |
| B010 | <u>AbsI:</u> | 10 mM Tris-HCl (pH 9.0 at 25° C); 10 mM MgCl ₂ ; 50 mM KCl; 1 mM DTT. |
| B020 | <u>AoxI</u> | 10 mM Tris-HCl (pH 7.5 at 25°C); 200 mM KCl; 0,1 mM EDTA, 7 mM 2-mercaptoethanol; 200 µg/ml BSA, 50% glycerol. |
| B011 | <u>EcoRI:</u> | 100 mM Tris-HCl (pH 7.6 at 25° C); 10 mM MgCl ₂ ; 50 mM NaCl; 1 mM DTT. |
| B012 | <u>BisI</u> | 10 mM Tris-HCl (pH 9.0 at 25° C); 10 mM MgCl ₂ ; 150 mM KCl; 1 mM DTT. |
| B016 | <u>Mall</u> | 20 mM Tris-HCl (pH 9.0 at 25° C); 10 mM MgCl ₂ ; 150 mM NaCl; 1 mM DTT. |
| B017 | <u>N-Bst9I:</u> | 10 mM Tris-HCl (pH 8.5 at 25° C); 10 mM MgCl ₂ ; 150 mM KCl; 1 mM DTT. |
| B019 | <u>PcsI</u> | 10 mM Tris-HCl (pH 8.3 at 25° C); 20 mM NaCl; 3 mM MgCl ₂ ; 1 mM DTT. |
| B018 | <u>RigI:</u> | 10 mM Tris-HCl (pH 8.5 at 25° C); 5 mM MgCl ₂ ; 1 mM DTT |
| B021 | <u>ROSE</u> | |
| B022 | <u>GLAD:</u> | 50 mM Tris-SO ₄ , pH 9.0 at 25°C; 30 mM KCl; 10 mM [NH ₄] ₂ SO ₄ |
| B301 | <u>T4-Polynucleotide Kinase:</u> | 50 mM Tris-HCl (pH 7.6 at 25°C); 10 mM MgCl ₂ ; 5 mM DTT. |
| B302 | <u>T4-DNA Ligase:</u> | 50 mM Tris-HCl (pH 7.8 at 25°C); 10 mM MgCl ₂ ; 10 mM DTT; 1 mM ATP. Storage conditions: by small portions avoiding repeated defrosting to prevent ATP decomposition. |
| B303 | <u>T4-RNA Ligase:</u> | 50 mM Tris-HCl (pH 7.8 at 25°C); 10 mM MgCl ₂ ; 10 mM DTT; 1 mM ATP. |
| B304 | <u>DNA polymerase I E.coli (Klenow Fragment):</u> | 50 mM Tris-HCl (pH 7.6 at 25°C); 10 mM MgCl ₂ ; 5 mM DTT. |
| B309 | <u>Hot Start Taq DNA-polymerase:</u> | 67 mM Tris-HCl (pH 8.8 at 25 C), 16.6 mM (NH ₄) ₂ SO ₄ , 0.01 % Tween-20. Supplementary material is 50 mM MgCl ₂ . |
| B310 | <u>Pfu DNA- polymerase:</u> | 20 mM Tris-HCl (pH 8.8 at 25° C), 10 mM KCl, 10 mM (NH ₄) ₂ SO ₄ , 2 mM MgSO ₄ , 0.1% Triton X-100. |
| B321 | <u>AS (Ammonium Sulfate):</u> | 67 mM Tris-HCl (pH 8.8 at 25 C), 16.6 mM (NH ₄) ₂ SO ₄ , 0.01 % Tween-20. Supplementary material is 50 mM MgCl ₂ . |
| B305 | <u>Taq-DNA- polymerase , TaqSE-DNA- polymerase:</u> | 60 mM Tris-HCl (pH 8.5 at 25°C); 1.5 mM MgCl ₂ ; 25 mM KCl; 10 mM 2- mercaptoethanol; 0.1% Triton X-100. |
| B306 | <u>Taq-DNA- polymerase (Mg²⁺ – free):</u> | 60 mM Tris-HCl (pH 8.5 at 25°C); 25 mM KCl; 10 mM 2- mercaptoethanol; 0.1% Triton X-100. |
| B311 | <u>T4- DNA- polymerase:</u> | 67 mM Tris-HCl (pH 8.8 at 25°C); 6.7 mM MgCl ₂ ; 16.7 mM (NH ₄) ₂ SO ₄ ; 1 mM DTT. |
| B319 | <u>T7- RNA- polymerase:</u> | 50 mM Tris-HCl (pH 7.5 at 25°C); 6 mM MgCl ₂ ; 10 mM DTT; 2 mM spermidine. |
| B312 | <u>M-MuLV reverse transcriptase</u> | 50 mM Tris-HCl (pH 8.3 at 25°C); 75 mM KCl; 3 mM MgCl ₂ , 10 mM DTT. |
| B313 | <u>Inorganic pyrophosphatase</u> | 50 mM Tris-HCl (pH 8.5 at 25° C); 1 mM MgCl ₂ |
| B314 | <u>BAL-31 nuclease:</u> | 20 mM Tris-HCl (pH 8.0 at 25°C); 600 mM NaCl; 12 mM CaCl ₂ ; 12 mM MgCl ₂ ; 1 mM EDTA. |
| B315 | <u>Mung Bean nuclease</u> | 50 mM sodium acetate; 30 mM NaCl; 1 mM ZnSO ₄ ; (pH 5.0 at 25°C). |
| B316 | <u>Exonuclease III:</u> | 50 mM Tris-HCl (pH 7.6 at 25°C); 1 mM MgCl ₂ . |
| B317 | <u>Endonuclease I:</u> | 20 mM Glycin -NaOH (pH 9.5 at 25° C); 25 mM MgCl ₂ ; 100mM NaCl; 1 mM 2- mercaptoethanol. |
| B318 | <u>Uracil-DNA-glycosylase:</u> | 20 mM Tris-HCl (pH 8.0 at 25° C); 1 mM EDTA; 1 mM DTT. |
| B003 | <u>Alkaline phosphatase:</u> | 50 mM Tris-HCl (pH 7.6 at 25° C); 10 mM MgCl ₂ ; 100 mM NaCl; 1 mM DTT. |

| | | |
|-------------|--|--|
| B100 | <u>A (Storage and dilution buffer) :</u> | 10 mM Tris-HCl (pH 7.6 at 25° C); 50 mM KCl; 0,1 mM EDTA; 200 µg/ml BSA; 1 mM DTT; 50% glycerol. |
| B102 | <u>M-MuLV Reverse Transcriptase dilution Buffer</u> | 10 mM KH ₂ PO ₄ (pH 7.5); 0,1 mM EDTA; 200 mM NaCl; 7 mM 2-mercaptoethanol; 50% glycerol |
| B307 | <u>MgCl₂, 50 mM water solution</u> | 50 mM MgCl ₂ ; 500 µl |
| B101 | <u>BSA (for Restrictases)</u> | 10 mg/ml BSA. ; 500 µl |

SE Buffer Activity Chart for Restriction Enzymes

| Enzyme | Recognition sequence | SE Buffer | BSA | Activity (% from maximum) | | | | | Optimum t, °C | Inactivation, 20min |
|-----------|--------------------------------------|-----------|-----|---------------------------|--------|--------|--------|--------|---------------|---------------------|
| | | | | B | G | O | W | Y | | |
| Aat II | GACGT [^] C | Y | - | 10-25 | 25-50 | 10-25 | 25-50 | 100 | 37 | 65°C |
| Abs I | CC [^] TCGAGG | * | - | 75-100 | 10-25 | 0 | 50-75 | 0-10 | 37 | 65°C |
| Acc16 I | TGC [^] GCA | W | - | 50-75 | 75-100 | 25-50 | 100 | 75-100 | 37 | 65°C |
| Acc36 I | ACCTGC(4/8) | Y | - | 25-50 | 25-50 | 50-75 | 50-75 | 100 | 37 | 65°C |
| Acc65 I | G [^] GTACC | W | - | 10-25 | 25-50 | 75-100 | 100 | 10-25 | 37 | 65°C |
| AccB1 I | G [^] GYRCC | K | + | 50-75 | 10-25 | 10-25 | 75-100 | 50-75 | 37 | 65°C |
| AccB7 I | CCANNN [^] NTGG | G | - | 10-25 | 100 | 25-50 | 50-75 | 50-75 | 37 | 65°C |
| AccBS I | GAG [^] CGG | Y | - | 75-100 | 75-100 | 25-50 | 25-50 | 100 | 37 | 65°C |
| Acl I | AA [^] CGTT | Y | + | 0-10 | 0-10 | 0-10 | 0-10 | 100 | 37 | 65°C |
| AclW I | GGATC(4/5) | Y | + | 75-100 | 50-75 | 0-10 | 0-10 | 100 | 37 | 65°C |
| Aco I | Y [^] GGCCR | G | - | 50-75 | 100 | 50-75 | 25-50 | 75-100 | 37 | 65°C |
| Acs I | R [^] AATTY | W | + | 25-50 | 50-75 | 50-75 | 100 | 10-25 | 50 | 80°C |
| Acu I | CTGAAG(16/14) | Y+SAM | + | 25-50 | 50-75 | 50-75 | 75-100 | 100 | 37 | 65°C |
| Afe I | AGC [^] GCT | Y | - | 10-25 | 25-50 | 75-100 | 75-100 | 100 | 37 | 65°C |
| Ags I | TTS [^] AA | Y | + | 75-100 | 50-75 | 10-25 | 10-25 | 100 | 37 | 65°C |
| Ahl I | A [^] CTAGT | B | + | 100 | 75-100 | 25-50 | 25-50 | 75-100 | 37 | No |
| Ajn I | [^] CCWGG | Y | - | 25-50 | 10-25 | 10-25 | 25-50 | 100 | 55 | 65°C |
| Alu I | AG [^] CT | Y | - | 75-100 | 75-100 | 10-25 | 50-75 | 100 | 37 | 65°C |
| AluB I | AG [^] CT | B | + | 100 | 75-100 | 10-25 | 10-25 | 75-100 | 37 | 65°C |
| Ama87 I | C [^] YCGRG | W | + | 10-25 | 50-75 | 75-100 | 100 | 0-10 | 37 | 65°C |
| Aox I | [^] PuG(5mC)Py | * | - | 75-100 | 25-50 | 10-25 | 25-50 | 75-100 | 60 | No |
| Apa I | GGGCC [^] C | Y | + | 50-75 | 25-50 | 0-10 | 0-10 | 100 | 37 | 65°C |
| Ars I | (8/13)GAC(N) ₆ TTYG(11/6) | Y | + | 0 | 0 | 0 | 0 | 100 | 30 | 65°C |
| AsiG I | A [^] CCGGT | O | - | 10-25 | 25-50 | 100 | 75-100 | 10-25 | 37 | 65°C |
| AspA2 I | C [^] CTAGG | W | + | 10-25 | 50-75 | 75-100 | 100 | 75-100 | 37 | 80°C |
| AspLE I | GCG [^] C | O | - | 0-25 | 75-100 | 100 | 50-75 | 25-50 | 37 | 80°C |
| AspS9 I | G [^] GNCC | W | - | 50-75 | 50-75 | 75-100 | 100 | 50-75 | 37 | 65°C |
| AsuC2 I | CC [^] SGG | Y | - | 75-100 | 50-75 | 10-25 | 25-50 | 100 | 37 | 65°C |
| AsuHP I | GGTGA(8/7) | O | - | 10-25 | 50-75 | 100 | 75-100 | 25-50 | 37 | 65°C |
| AsuNH I | G [^] CTAGC | Y | + | 75-100 | 50-75 | 0-10 | 0-10 | 100 | 37 | 65°C |
| BamH I | G [^] GATCC | G | + | 25-50 | 100 | 75-100 | 75-100 | 25-50 | 37 | 65°C |
| Bar I | (7/12)GAAG(N) ₆ TAC(12/7) | 2K | - | 0 | 0-10 | 25-50 | 50-75 | 10-25 | 37 | 65°C |
| Bbv12 I | GWGCW [^] C | O | - | 0-10 | 10-25 | 100 | 75-100 | 10-25 | 37 | 80°C |
| Bgl I | GCCNNNN [^] NGGC | 2W | - | 50-75 | 50-75 | 0-10 | 75-100 | 25-50 | 37 | 65°C |
| Bgl II | A [^] GATCT | O | - | 0-10 | 10-25 | 100 | 25-50 | 10-25 | 37 | 80°C |
| Bis I | G(5mC) [^] NGC | * | - | 10-25 | 25-50 | 50-75 | 75-100 | 50-75 | 37 | 65°C |
| Bls I | PuPyN [^] PuPy | W | - | 10-25 | 10-25 | 50-75 | 100 | 75-100 | 30 | 65°C |
| Bme18 I | G [^] GWCC | O | - | 10-25 | 25-50 | 100 | 75-100 | 10-25 | 37 | 65°C |
| Bmt I | GCTAG [^] C | W | - | 10-25 | 50-75 | 50-75 | 100 | 75-100 | 37 | 65°C |
| BmuI | ACTGGG(5/4) | Y | - | 75-100 | 75-100 | 25-50 | 50-75 | 100 | 37 | 65°C |
| Bpm I | CTGGAG(16/14) | W | + | 25-50 | 50-75 | 75-100 | 100 | 50-75 | 37 | 65°C |
| Bpu10 I | CC [^] TNAGC | K | - | 10-25 | 25-50 | 50-75 | 50-75 | 25-50 | 37 | 80°C |
| Bpu14 I | TT [^] CGAA | G | - | 50-75 | 100 | 25-50 | 25-50 | 75-100 | 37 | 65°C |
| Bsa29 I | AT [^] CGAT | G | + | 25-50 | 100 | 50-75 | 50-75 | 75-100 | 37 | 65°C |
| Bsc4 I | CCNNNN [^] NNGG | W | + | 75-100 | 75-100 | 50-75 | 100 | 25-50 | 55 | 80°C |
| BseI | ACTGG(1/-1) | Y | - | 75-100 | 75-100 | 25-50 | 10-25 | 100 | 65 | 80°C |
| Bse118 I | R [^] CCGGY | O | - | 0-10 | 50-75 | 100 | 75-100 | 25-50 | 65 | 80°C |
| Bse21 I | CC [^] TNAGG | Y | - | 50-75 | 50-75 | 10-25 | 25-50 | 100 | 37 | 80°C |
| Bse3D I | GCAATG(2/0) | G | - | 10-25 | 100 | 25-50 | 50-75 | 75-100 | 60 | 80°C |
| Bse8 I | GATNN [^] NNATC | G | - | 25-50 | 100 | 75-100 | 75-100 | 50-75 | 60 | 80°C |
| BseP I | G [^] CGCGC | G | - | 50-75 | 100 | 75-100 | 50-75 | 50-75 | 50 | 65°C |
| BseX3 I | C [^] GGCCG | O | - | 10-25 | 25-50 | 100 | 50-75 | 10-25 | 50 | 80°C |
| BsIF I | GGGAC(10/14) | Y | + | 25-50 | 25-50 | 10-25 | 25-50 | 100 | 37 | 80°C |
| Bso31 I | GGTCTC(1/5) | O | + | 25-50 | 75-100 | 100 | 75-100 | 25-50 | 55 | 80°C |
| Bsp13 I | T [^] CCGGA | 2K | - | 25-50 | 50-75 | 75-100 | 50-75 | 0-10 | 50 | 65°C |
| Bsp1720 I | GC [^] TNAGC | G | - | 50-75 | 100 | 50-75 | 50-75 | 75-100 | 37 | 80°C |
| Bsp19 I | C [^] CATGG | 2W | + | 0-10 | 10-25 | 50-75 | 75-100 | 10-25 | 37 | 65°C |

* - supplied with its own unique reaction buffer that is different from the five standard SE-Buffers.

SE Buffer Activity Chart for Restriction Enzymes

| Enzyme | Recognition sequence | SE Buffer | BSA | Activity (% from maximum) | | | | | Optimum t,°C | Inactivation, 20min |
|----------|-------------------------------------|-----------|-----|---------------------------|--------|--------|--------|--------|--------------|---------------------|
| | | | | B | G | O | W | Y | | |
| BspAC I | CCGC(-3/-1) | O | + | 10-25 | 25-50 | 100 | 75-100 | 10-25 | 37 | 65°C |
| BspFN I | CG^CG | Y | - | 50-75 | 75-100 | 75-100 | 50-75 | 100 | 37 | 65°C |
| BssEC I | C^CNNGG | Y | - | 50-75 | 50-75 | 50-75 | 75-100 | 100 | 60 | 80°C |
| BssNA I | GTA^TAC | W | + | 50-75 | 50-75 | 75-100 | 100 | 75-100 | 37 | No |
| BssT1 I | C^CWGG | 2K | - | 10-25 | 25-50 | 25-50 | 75-100 | 10-25 | 60 | 80°C |
| Bst2B I | CTCGTG(-5/-1) | Y | + | 75-100 | 25-50 | 10-25 | 25-50 | 100 | 60 | 80°C |
| Bst2U I | CC^WGG | G | + | 75-100 | 100 | 50-75 | 50-75 | 10-25 | 60 | 80°C |
| Bst4C I | ACN^GT | Y | - | 75-100 | 75-100 | 10-25 | 25-50 | 100 | 65 | 80°C |
| Bst6 I | CTCTTC(1/4) | Y | + | 75-100 | 75-100 | 50-75 | 75-100 | 100 | 65 | 80°C |
| BstAC I | GR^CGYC | W | - | 75-100 | 75-100 | 50-75 | 100 | 75-100 | 37 | 80°C |
| BstAF I | C^TTAAG | W | + | 10-25 | 25-50 | 75-100 | 100 | 25-50 | 55 | 80°C |
| BstAP I | GCANNN^NTGC | W | - | 25-50 | 25-50 | 75-100 | 100 | 25-50 | 60 | 80°C |
| BstAU I | T^GTACA | W | - | 10-25 | 50-75 | 25-50 | 100 | 25-50 | 37 | 80°C |
| BstBA I | YAC^GTR | W | + | 25-50 | 25-50 | 75-100 | 100 | 25-50 | 65 | 80°C |
| BstC8 I | GCN^NGC | Y | - | 10-25 | 25-50 | 50-75 | 75-100 | 100 | 55 | 80°C |
| BstDE I | C^TNAG | G | - | 75-100 | 100 | 25-50 | 50-75 | 10-25 | 60 | 80°C |
| BstDS I | C^CRYGG | Y | - | 0-10 | 75-100 | 50-75 | 25-50 | 100 | 65 | 80°C |
| BstEN I | CCTNN^NNNAGG | Y | - | 50-75 | 50-75 | 25-50 | 25-50 | 100 | 65 | 80°C |
| BstF5 I | GGATG(2/0) | Y | - | 75-100 | 50-75 | 25-50 | 50-75 | 100 | 65 | 80°C |
| BstFN I | CG^CG | Y | - | 75-100 | 50-75 | 25-50 | 25-50 | 100 | 60 | 80°C |
| BstH2 I | RGCGC^Y | Y | + | 50-75 | 50-75 | 0-10 | 10-25 | 100 | 65 | 80°C |
| BstHH I | GCG^C | Y | + | 75-100 | 50-75 | 25-50 | 50-75 | 100 | 50 | No |
| BstKI I | GAT^C | W | - | 25-50 | 50-75 | 75-100 | 100 | 50-75 | 37 | 65°C |
| BstMA I | GTCTC(1/5) | W | + | 25-50 | 50-75 | 50-75 | 100 | 75-100 | 55 | 65°C |
| BstMB I | ^GATC | O | - | 10-25 | 25-50 | 100 | 75-100 | 10-25 | 65 | 80°C |
| BstMC I | CGRY^CG | B | + | 100 | 75-100 | 10-25 | 10-25 | 50-75 | 50 | 80°C |
| BstMW I | GCNNNN^NNGC | Y | - | 10-25 | 25-50 | 25-50 | 50-75 | 100 | 55 | 80°C |
| BstNS I | RCATG^Y | B | + | 100 | 50-75 | 10-25 | 10-25 | 75-100 | 37 | 65°C |
| BstPA I | GACNN^NNGTC | Y | - | 50-75 | 25-50 | 50-75 | 50-75 | 100 | 65 | No |
| BstSC I | ^CCNGG | Y | - | 50-75 | 50-75 | 50-75 | 50-75 | 100 | 55 | 80°C |
| BstSF I | C^TRYAG | O | + | 75-100 | 25-50 | 100 | 50-75 | 50-75 | 60 | No |
| BstSL I | GKGCM^C | G | + | 50-75 | 100 | 50-75 | 75-100 | 75-100 | 55 | 65°C |
| BstSN I | TAC^GTA | B | - | 100 | 50-75 | 0-10 | 10-25 | 50-75 | 37 | 80°C |
| BstV1 I | GCAGC(8/12) | G | - | 75-100 | 100 | 75-100 | 75-100 | 75-100 | 55 | 80°C |
| BstV2 I | GAAGAC(2/6) | Y | + | 75-100 | 75-100 | 25-50 | 25-50 | 100 | 55 | 65°C |
| BstX I | CCANNNNN^NTGG | O | - | 10-25 | 10-25 | 100 | 75-100 | 25-50 | 37 | 65°C |
| BstX2 I | R^GATCY | G | - | 75-100 | 100 | 0-10 | 10-25 | 25-50 | 60 | 80°C |
| Bsu I | GTATCC(6/5) | Y | - | 75-100 | 50-75 | 10-25 | 25-50 | 100 | 37 | 65°C |
| BsuR I | GG^CC | G | - | 75-100 | 100 | 25-50 | 50-75 | 50-75 | 37 | 80°C |
| Btr I | CACGTC(-3/-3) | O | + | 75-100 | 75-100 | 100 | 75-100 | 75-100 | 60 | 80°C |
| Cci I | T^CATGA | W | + | 0-10 | 10-25 | 25-50 | 100 | 75-100 | 55 | 80°C |
| CciN I | GC^GGCCGC | Y | - | 25-50 | 50-75 | 75-100 | 75-100 | 100 | 37 | 65°C |
| Dra I | TTT^AAA | G | + | 75-100 | 100 | 25-50 | 75-100 | 75-100 | 37 | 65°C |
| Dra III | CACNNN^GTG | 2K | + | 25-50 | 50-75 | 75-100 | 75-100 | 50-75 | 37 | 65°C |
| Dri I | GACNNN^NNGTC | Y | - | 75-100 | 75-100 | 10-25 | 10-25 | 100 | 37 | 65°C |
| DseD I | GACNNNN^NNGTC | Y | + | 75-100 | 75-100 | 25-50 | 50-75 | 100 | 37 | 80°C |
| EcoICR I | GAG^CTC | G | + | 75-100 | 100 | 0-10 | 0-10 | 75-100 | 37 | 65°C |
| EcoR I | G^AATTC | * | + | 50-75 | 75-100 | 75-100 | 100 | 50-75 | 37 | 65°C |
| EcoR V | GAT^ATC | W | + | 0-10 | 25-50 | 50-75 | 100 | 25-50 | 37 | 80°C |
| Ege I | GGC^GCC | B | + | 100 | 75-100 | 10-25 | 50-75 | 75-100 | 37 | 65°C |
| Erh I | C^CWGG | 2W | + | 10-25 | 25-50 | 25-50 | 75-100 | 10-25 | 37 | 65°C |
| Fae I | CATG^ | * | + | 25-50 | 50-75 | 10-25 | 10-25 | 75-100 | 37 | 65°C |
| Fai I | YA^TR | B | - | 100 | 50-75 | 10-25 | 25-50 | 25-50 | 50 | 80°C |
| Fal I | (8/13)AAG(N) ₅ CTT(13/8) | W+SAM | - | 0-10 | 25-50 | 75-100 | 100 | 50-75 | 37 | 65°C |
| Fat I | ^CATG | G | - | 10-25 | 100 | 25-50 | 10-25 | 50-75 | 55 | 65°C |
| Fau I | CCC GC(4/6) | B | - | 100 | 25-50 | 0-10 | 0-10 | 50-75 | 55 | 65°C |
| FauND I | CA^TATG | Y | + | 50-75 | 75-100 | 10-25 | 50-75 | 100 | 37 | 65°C |
| Fbl I | GT^MKAC | Y | - | 50-75 | 75-100 | 0-10 | 50-75 | 100 | 55 | 80°C |
| Fok I | GGATG(9/13) | Y | - | 50-75 | 50-75 | 25-50 | 25-50 | 100 | 37 | 65°C |

* - supplied with its own unique reaction buffer that is different from the five standard SE-Buffers.

SE Buffer Activity Chart for Restriction Enzymes

| Enzyme | Recognition sequence | SE Buffer | BSA | Activity (% from maximum) | | | | | Optimum t, °C | Inactivation, 20min |
|-----------|--------------------------------------|-----------|-----|---------------------------|--------|--------|--------|--------|---------------|---------------------|
| | | | | B | G | O | W | Y | | |
| FriO I | GRGCY^C | Y | + | 75-100 | 75-100 | 10-25 | 0-10 | 100 | 37 | 65°C |
| Fsp4H I | GC^NGC | Y | - | 50-75 | 75-100 | 10-25 | 25-50 | 100 | 37 | 65°C |
| Fbl I | GT^MKAC | Y | - | 50-75 | 75-100 | 0-10 | 50-75 | 100 | 55 | 80°C |
| Fok I | GGATG(9/13) | Y | - | 50-75 | 50-75 | 25-50 | 25-50 | 100 | 37 | 65°C |
| FriO I | GRGCY^C | Y | + | 75-100 | 75-100 | 10-25 | 0-10 | 100 | 37 | 65°C |
| Fsp4H I | GC^NGC | Y | - | 50-75 | 75-100 | 10-25 | 25-50 | 100 | 37 | 65°C |
| Gla I | Pu(5mC)^GPy | Y | - | 75-100 | 75-100 | 25-50 | 25-50 | 100 | 30 | 65°C |
| Glu I | G(5mC)^NG(5mC) | Y | - | 75-100 | 75-100 | 25-50 | 50-75 | 100 | 37 | 80°C |
| Gsa I | CCCAG^C | W | + | 10-25 | 25-50 | 75-100 | 100 | 75-100 | 70 | No |
| Hae III | GG^CC | G | - | 75-100 | 100 | 25-50 | 50-75 | 50-75 | 37 | 80°C |
| Hga I | GACGC(5/10) | B | - | 100 | 75-100 | 10-25 | 25-50 | 50-75 | 37 | 65°C |
| Hind II | GTY^RAC | G | + | 75-100 | 100 | 25-50 | 25-50 | 75-100 | 37 | 65°C |
| Hind III | A^AGCTT | W | + | 10-25 | 25-50 | 0-10 | 100 | 0-10 | 37 | 80°C |
| Hinf I | G^ANTC | O | - | 25-50 | 75-100 | 100 | 75-100 | 75-100 | 37 | 80°C |
| Hpa I | GTT^AAC | Y | - | 0-10 | 50-75 | 10-25 | 25-50 | 100 | 37 | 65°C |
| Hpa II | C^CGG | B | - | 100 | 50-75 | 10-25 | 25-50 | 50-75 | 37 | 80°C |
| HpySE526I | A^CGT | Y | - | 75-100 | 75-100 | 10-25 | 25-50 | 100 | 37 | 65°C |
| HspA I | G^CGC | Y | - | 50-75 | 50-75 | 25-50 | 25-50 | 100 | 37 | 80°C |
| Kpn I | GGTAC^C | B | + | 100 | 25-50 | 25-50 | 25-50 | 75-100 | 37 | 80°C |
| Kro I | G^C(5mC)GGC | G | - | 50-75 | 100 | 25-50 | 50-75 | 75-100 | 37 | 65°C |
| Ksp22 I | T^GATCA | 2K | + | 50-75 | 100 | 50-75 | 50-75 | 25-50 | 37 | 65°C |
| Kzo9 I | ^GATC | G | - | 50-75 | 100 | 50-75 | 50-75 | 50-75 | 37 | 65°C |
| Lmn I | GCTCCN^ | B | - | 100 | 75-100 | 50-75 | 50-75 | 75-100 | 37 | 65°C |
| Mab I | A^CCWGGT | W | + | 25-50 | 50-75 | 75-100 | 100 | 50-75 | 37 | 65°C |
| Mal I | G(mA)^TC | * | - | 25-50 | 25-50 | 50-75 | 75-100 | 50-75 | 37 | 65°C |
| Mbo II | GAAGA(8/7) | Y | - | 75-100 | 75-100 | 25-50 | 50-75 | 100 | 37 | 65°C |
| Mfe I | C^AATTG | B | + | 100 | 75-100 | 10-25 | 25-50 | 75-100 | 37 | No |
| Mhl I | GDGCH^C | W | - | 10-25 | 25-50 | 75-100 | 100 | 10-25 | 37 | 80°C |
| Mlu I | A^CGCGT | O | - | 0-10 | 10-25 | 100 | 25-50 | 10-25 | 37 | 65°C |
| Mly113 I | GG^CGCC | B | - | 100 | 25-50 | 10-25 | 10-25 | 50-75 | 37 | 65°C |
| Mnl I | CCTC(7/6) | G | + | 75-100 | 100 | 25-50 | 25-50 | 75-100 | 37 | 65°C |
| Mox20 I | TGG^CCA | O | - | 10-25 | 25-50 | 100 | 75-100 | 25-50 | 37 | No |
| MroN I | G^CCGGC | B | - | 100 | 50-75 | 10-25 | 0-10 | 10-25 | 37 | 80°C |
| MroX I | GAANN^NNTTC | W | - | 50-75 | 50-75 | 50-75 | 100 | 25-50 | 37 | 65°C |
| Msp I | C^CGG | B | - | 100 | 75-100 | 50-75 | 75-100 | 75-100 | 37 | 65°C |
| MspA1 I | CMG^CKG | Y | + | 10-25 | 75-100 | 10-25 | 25-50 | 100 | 37 | 65°C |
| MspR9 I | CC^NGG | O | - | 50-75 | 50-75 | 100 | 50-75 | 75-100 | 37 | 80°C |
| Mte I | G(5mC)G(5mC)^NG(5mC)G(5mC) | W | - | 25-50 | 75-100 | 75-100 | 100 | 50-75 | 55 | No |
| Nru I | TCG^CGA | W | - | 0-10 | 10-25 | 75-100 | 100 | 10-25 | 37 | 80°C |
| PalA I | GG^CGCGCC | Y | - | 25-50 | 10-25 | 0 | 0 | 100 | 37 | 65°C |
| Pce I | AGG^CCT | Y | - | 75-100 | 75-100 | 50-75 | 25-50 | 100 | 50 | 80°C |
| Pci I | A^CATGT | O | - | 50-75 | 75-100 | 100 | 75-100 | 50-75 | 37 | 65°C |
| PciS I | GCTCTC(1/4) | B | - | 100 | 50-75 | 0-10 | 0-10 | 75-100 | 37 | 65°C |
| Pcs I | (5mC)GNNNNN^NN(5mC)G | * | - | 50-75 | 25-50 | 0 | 10-25 | 50-75 | 37 | 65°C |
| Pct I | GAATGC(1/-1) | O | - | 25-50 | 50-75 | 100 | 75-100 | 10-25 | 37 | 65°C |
| Pkr I | G(5mC)N^G(5mC) | Y | - | 50-75 | 75-100 | 10-25 | 25-50 | 100 | 37 | 65°C |
| Ple19 I | CGAT^CG | Y | - | 75-100 | 75-100 | 25-50 | 25-50 | 100 | 37 | 65°C |
| Pps I | GAGTC(4/5) | Y | + | 50-75 | 10-25 | 0-10 | 25-50 | 100 | 37 | 65°C |
| Psi I | TTA^TAA | B | - | 100 | 25-50 | 10-25 | 25-50 | 75-100 | 37 | 65°C |
| Psp124B I | GAGCT^C | G | - | 75-100 | 100 | 10-25 | 0-10 | 75-100 | 37 | 80°C |
| Psp6I | ^CCWGG | B | - | 100 | 50-75 | 10-25 | 25-50 | 75-100 | 55 | 80°C |
| PspC I | CAC^GTG | B | + | 100 | 50-75 | 0 | 0 | 50-75 | 37 | 65°C |
| PspE I | G^GTNACC | B | - | 100 | 50-75 | 25-50 | 50-75 | 50-75 | 37 | 65°C |
| PspL I | C^GTACG | Y | + | 75-100 | 75-100 | 25-50 | 10-25 | 100 | 37 | 65°C |
| PspN4 I | GGN^NCC | Y | - | 10-25 | 10-25 | 10-25 | 25-50 | 100 | 37 | 65°C |
| PspOM I | G^GGCCC | Y | - | 75-100 | 10-25 | 0-10 | 0-10 | 100 | 37 | 65°C |
| PspPP I | RG^GWCCY | Y | + | 50-75 | 25-50 | 10-25 | 10-25 | 100 | 37 | 65°C |
| PspX I | VC^TCGAGB | Y | + | 50-75 | 50-75 | 25-50 | 75-100 | 100 | 37 | 80°C |
| Psr I | (7/12)GAAC(N) ₆ TAC(12/7) | Y | + | 10-25 | 10-25 | 0 | 0-10 | 100 | 30 | 65°C |
| Pst I | CTGCA^G | O | + | 10-25 | 25-50 | 100 | 25-50 | 25-50 | 37 | 80°C |

* - supplied with its own unique reaction buffer that is different from the five standard SE-Buffers.

SE Buffer Activity Chart for Restriction Enzymes

| Enzyme | Recognition sequence | SE Buffer | BSA | Activity (% from maximum) | | | | | Optimum t,°C | Inactivation, 20min |
|----------|----------------------|-----------|-----|---------------------------|--------|--------|--------|--------|--------------|---------------------|
| | | | | B | G | O | W | Y | | |
| PstN I | CAGNNN^CTG | Y | - | 75-100 | 50-75 | 10-25 | 25-50 | 100 | 37 | 65°C |
| Pvu II | CAG^CTG | G | + | 25-50 | 100 | 25-50 | 25-50 | 25-50 | 37 | 80°C |
| Rga I | GCGAT^CGC | Y | - | 75-100 | 50-75 | 10-25 | 25-50 | 100 | 55 | 80°C |
| Rig I | GGCCGG^CC | * | + | 75-100 | 50-75 | 0-10 | 10-25 | 50-75 | 37 | 65°C |
| Rsa I | GT^AC | B | - | 100 | 50-75 | 0-10 | 50-75 | 75-100 | 37 | 80°C |
| RsaN I | G^TAC | B | - | 100 | 75-100 | 50-75 | 50-75 | 75-100 | 37 | 80°C |
| Rsr2 I | CG^GWCCG | Y | + | 50-75 | 75-100 | 0-10 | 10-25 | 100 | 37 | 65°C |
| Sal I | G^TCGAC | O | - | 0-10 | 10-25 | 100 | 25-50 | 0-10 | 37 | 65°C |
| Sbf I | CCTGCA^GG | Y | - | 75-100 | 50-75 | 0-10 | 0-10 | 100 | 37 | 80°C |
| Set I | ASST^ | Y | - | 25-50 | 25-50 | 75-100 | 75-100 | 100 | 55 | 80°C |
| SfaN I | GCATC(5/9) | O | - | 10-25 | 25-50 | 100 | 75-100 | 0-10 | 37 | 80°C |
| Sfi I | GGCCNNNN^NGGCC | G | + | 75-100 | 100 | 25-50 | 25-50 | 25-50 | 50 | 65°C |
| Sfr274 I | C^TCGAG | B | - | 100 | 75-100 | 50-75 | 50-75 | 75-100 | 50 | 65°C |
| Sfr303 I | CCG^GG | B | - | 100 | 50-75 | 10-25 | 10-25 | 75-100 | 37 | 65°C |
| Sma I | CCC^GGG | Y | - | 0-10 | 0-10 | 0-10 | 0-10 | 100 | 25 | 65°C |
| Smi I | ATTT^AAAT | O | + | 25-50 | 25-50 | 100 | 75-100 | 25-50 | 37 | 65°C |
| SmiM I | CAYNN^NNRTG | W | - | 10-25 | 10-25 | 75-100 | 100 | 10-25 | 37 | 65°C |
| Sph I | GCATG^C | G | + | 25-50 | 100 | 75-100 | 75-100 | 50-75 | 37 | 65°C |
| Sse9 I | ^AATT | B | + | 100 | 75-100 | 50-75 | 50-75 | 75-100 | 55 | 65°C |
| Ssp I | AAT^ATT | K | + | 75-100 | 50-75 | 25-50 | 50-75 | 75-100 | 37 | 65°C |
| SspM I | C^TAG | Y | - | 50-75 | 25-50 | 10-25 | 50-75 | 100 | 55 | No |
| Taq I | T^CGA | Y | + | 50-75 | 75-100 | 75-100 | 50-75 | 100 | 65 | 80°C |
| Tru9 I | T^TAA | W | - | 75-100 | 25-50 | 25-50 | 100 | 50-75 | 65 | 80°C |
| TseF I | ^GTSAC | B | - | 100 | 50-75 | 0-10 | 25-50 | 50-75 | 65 | No |
| Tth111 I | GACN^NNGTC | Y | - | 75-100 | 50-75 | 10-25 | 10-25 | 100 | 65 | 80°C |
| Vne I | G^TGCAC | O | - | 10-25 | 25-50 | 100 | 25-50 | 25-50 | 37 | 65°C |
| Vsp I | AT^TAAT | W | - | 0-10 | 10-25 | 50-75 | 100 | 25-50 | 37 | 65°C |
| Xba I | T^CTAGA | O | + | 75-100 | 75-100 | 100 | 50-75 | 75-100 | 37 | 65°C |
| Xma I | C^CCGGG | Y | - | 75-100 | 50-75 | 0 | 0-10 | 100 | 37 | 65°C |
| Zra I | GAC^GTC | B | - | 100 | 50-75 | 25-50 | 25-50 | 75-100 | 37 | 80°C |
| Zrm I | AGT^ACT | Y | + | 50-75 | 25-50 | 0-10 | 0-10 | 100 | 37 | 65°C |
| Zsp2 I | ATGCA^T | B | + | 100 | 50-75 | 25-50 | 25-50 | 25-50 | 37 | 65°C |
| N-Bst9 I | GAGTC(4/-) | * | - | 10-25 | 75-100 | 100 | 100 | 50-75 | 55 | 80°C |

* - supplied with its own unique reaction buffer that is different from the five standard SE-Buffers.

Isoschizomers

| Enzyme | SE Enzyme | Enzyme | SE Enzyme | Enzyme | SE Enzyme | Enzyme | SE Enzyme | Enzyme | SE Enzyme |
|---------|---------------------|----------|-----------------------|----------|---------------------|----------|---------------------|----------|-----------------------|
| AanI | PsiI | AsuHPI | AsuHPI | BseBI | AjnI [^] | BspTNI | Bso31I | BstX2I | BstX2I |
| AasI | DseDI | AsuNHI | AsuNHI | BseBI | Bst2UI | BspXI | Bsa29I | BstYI | BstX2I |
| AatI | PceI | | BmtI [^] | | Psp6I [^] | BsrI | Bse1I | BstZI | BseX3I |
| AatII | AatII | AvaI | Ama87I | BseCI | Bsa29I | BsrBI | AccBSI | BstZ17I | BssNAI |
| AbsI | AbsI | AvaII | Bme18I | BseDI | BssECI | BsrDI | Bse3DI | Bsul | Bsul |
| AccI | FbII | AvaIII | Zsp2I [^] | Bse3DI | Bse3DI | BsrFI | Bse118I | Bsu15I | Bsa29I |
| AccII | BstFNI | AviII | Acc16I | BseGI | BstF5I | BsrGI | BstAUI | Bsu36I | Bse21I |
| AccIII | Bsp13I | AvrII | AspA2I | BseGI | FokI [^] | BsrSI | Bse1I | BsuRI | BsuRI |
| Acc16I | Acc16I | AxyI | Bse21I | BseJI | Bse8I | BssAI | Bse118I | | HaeIII |
| Acc36I | Acc36I | Ball | Mox20I | BseLI | Bsc4I | BssECI | BssECI | BsuTUI | Bsa29I |
| Acc65I | Acc65I | BamHI | BamHI | BseMI | Bse3DI | BssHI | Sfr274I | BtgI | BstDSI |
| | KpnI [^] | BanI | AccB1I | BseNI | Bse1I | BssHII | BsePI | BtrI | BtrI |
| AccB1I | AccB1I | BanII | FriOI | BsePI | BsePI | BssKI | MspR9I [^] | BveI | Acc36I |
| AccB7I | AccB7I | BanIII | Bsa29I | BseSI | BstSLI | | BstSCI | Cac8I | BstC8I |
| AccBSI | AccBSI | BarI | BarI | BseXI | BstV1I | BssNAI | BssNAI | CaiI | PstNI |
| Acil | BspACI | Bbel | Egel [^] | BseX3I | BseX3I | BssSI | Bst2BI | CauII | AsuC2I |
| AcII | AcII | | Mly113I ^{^^} | BseYI | GsaI | BssT1I | BssT1I | CciI | CciI |
| AcIWI | AcIWI | BbrPI | PspCI | Bsh1236I | BstFNI | | ErhI | CciNI | CciNI |
| AcoI | AcoI | BbsI | BstV2I | Bsh1285I | BstMCI | Bst6I | Bst6I | CellI | Bsp1720I |
| AcsI | AcsI | Bbul | SphI | BshFI | BsuRI | Bst98I | BstAFI | CfoI | AspLEI |
| AcuI | AcuI | BbvI | BstV1I | | HaeIII | Bst1107I | BssNAI | | BstHHI |
| AcvI | PspCI | BbvII | BstV2I | BshNI | AccB1I | BstACI | BstACI | | HspAI [^] |
| Acyl | BstACI | Bbv12I | Bbv12I | BshTI | AsiGI | BstAFI | BstAFI | CfrI | AcoI |
| Adel | DraIII | BclI | Ksp22I | Bsil | Bst2BI | BstAPI | BstAPI | Cfr9I | XmaI |
| AfaI | RsaI | BciVI | Bsul | BsiEI | BstMCI | BstAUI | BstAUI | | SmaI [^] |
| AfeI | AfeI | BcnI | AsuC2I | BsiHKA1 | Bbv12I | BstBI | Bpu14I | Cfr10I | Bse118I |
| AfIII | BstAFI | Bcul | AhlI | BsiHKCI | Ama87I | Bst2BI | Bst2BI | Cfr13I | AspS9I |
| AgeI | AsiGI | Bfil | Bmul | BsiSI | HpaI | BstBAI | BstBAI | Cfr42I | Sfr303I |
| AgsI | AgsI | BfimI | BstSFI | | MspI | Bst4CI | Bst4CI | Clal | Bsa29I |
| AhaIII | DraI | BfriI | BstAFI | BsiWI | PspLI | BstC8I | BstC8I | CpoI | Rsr2I |
| AhdI | Dril | BfribI | Zsp2I [^] | BsiYI | Bsc4I | BstDEI | BstDEI | CspI | Rsr2I |
| AhII | AhII | Bful | Bsul | Bsil | Bsc4I | BstDSI | BstDSI | Csp6I | Rsal [^] |
| AjnI | AjnI | BfuAI | Acc36I | BslFI | BslFI | BstEII | PspEI | Csp45I | Bpu14I |
| | Bst2UI [^] | BfuCI | BstMBI | BsmFI | BslFI [^] | BstENI | BstENI | CspAI | AsiGI |
| | Psp6I | | BstKTI [^] | BsmI | PctI | BstF5I | BstF5I | CviAII | FaeI [^] |
| AluI | AluI | | Kzo9I | BsmAI | BstMAI | | FokI [^] | | FatI [^] |
| | AluBI | BglI | BglI | Bso31I | Bso31I | BstFNI | BstFNI | DdeI | BstDEI |
| AluBI | AluBI | BglII | BglII | BsoBI | Ama87I | BstH2I | BstH2I | DpnI | MaiI |
| AlwI | AcIWI | BimI | AcIWI | BsoMAI | BstMAI | BstHHI | BstHHI | DpnII | BstMBI |
| Alw21I | Bbv12I | Bisl | Bisl | Bsp13I | Bsp13I | | AspLEI | | Kzo9I |
| Alw26I | BstMAI | BlnI | AspA2I | Bsp19I | Bsp19I | | HspAI [^] | | BstKTI [^] |
| Alw44I | Vnel | Blpl | Bsp1720I | Bsp68I | NruI | BstKTI | BstKTI | DraI | DraI |
| AlwNI | PstNI | BlsI | BlsI | Bsp106I | Bsa29I | | BstMBI [^] | DraIII | DraIII |
| Ama87I | Ama87I | Bme18I | Bme18I | Bsp119I | Bpu14I | | Kzo9I [^] | DrdI | DseDI |
| Aor51HI | AfeI | Bme1390I | MspR9I | Bsp120I | PspOMI | BstMBI | BstKTI [^] | Dril | DseDI |
| AoxI | AoxI | | BstSCI [^] | | ApaI [^] | | BstMBI | DsaI | BstDSI |
| Apal | Apal | BmgBI | BtrI | Bsp143I | BstMBI | | Kzo9I | DseDI | DseDI |
| | PspOMI [^] | BmtI | BmtI | | Kzo9I | BstMCI | BstMCI | EaeI | AcoI |
| ApaBI | BstAPI [^] | | AsuNHI [^] | | BstKTI [^] | BstMWI | BstMWI | Eagl | BseX3I |
| ApalI | Vnel | BmyI | MhII | Bsp143II | BstH2I | BstNI | AjnI [^] | Eam1104I | Bst6I |
| ApoI | AcsI | BoxI | BstPAI | Bsp1286I | MhII | | Bst2UI | Eam1105I | Dril |
| ArsI | ArsI | Bpil | BstV2I | Bsp1407I | BstAUI | | Psp6I [^] | EarI | Bst6I |
| Ascl | PalAI | Bpml | Bpml | Bsp1720I | Bsp1720I | BstNSI | BstNSI | Ecl136II | EcoICRI |
| AseI | VspI | Bpu10I | Bpu10I | BspACI | BspACI | BstOI | AjnI [^] | | Psp124BI [^] |
| AsiGI | AsiGI | Bpu14I | Bpu14I | BspANI | BsuRI | | Bst2UI | EclHKI | Dril |
| AsiSI | AsiSI | Bpu1102I | Bsp1720I | | HaeIII | | Psp6I [^] | EclXI | BseX3I |
| AspI | Tth111I | BpuAI | BstV2I | BspCI | Ple19I | BstPI | PspEI | Eco24I | FriOI |
| Asp700I | MroXI | Bsal | Bso31I | BspDI | Bsa29I | BstPAI | BstPAI | Eco31I | Bso31I |
| Asp718I | Acc65I | Bsa29I | Bsa29I | BspEI | Bsp13I | BstSCI | BstSCI | Eco32I | EcoRV |
| | KpnI [^] | BsaAI | BstBAI | BspFNI | BspFNI | | MspR9I [^] | Eco47I | Bme18I |
| AspA2I | AspA2I | BsaBI | Bse8I | BspHI | Ceil | BstSFI | BstSFI | Eco47III | AfeI |
| AspEI | Dril | BsaHI | BstACI | BspLL | PspN4I | BstSLI | BstSLI | Eco52I | BseX3I |
| AspHI | Bbv12I | BsaJI | BssECI | BspLU11I | PciI | BstSNI | BstSNI | Eco57I | AcuI |
| AspLEI | AspLEI | BsaMI | PctI | BspMI | Acc36I | BstUI | BstFNI | Eco72I | PspCI |
| | BstHHI | Bsc4I | Bsc4I | BspMII | Bsp13I | Bst2UI | AjnI [^] | Eco81I | Bse21I |
| | HspAI [^] | Bse1I | Bse1I | BspMAI | PstI | | Bst2UI | Eco88I | Ama87I |
| AspS9I | AspS9I | Bse8I | Bse8I | BspPI | AcIWI | | Psp6I [^] | Eco91I | PspEI |
| AsuI | AspS9I | Bse21I | Bse21I | BspTI | BstAFI | BstV1I | BstV1I | Eco105I | BstSNI |
| AsuII | Bpu14I | Bse118I | Bse118I | BspT104I | Bpu14I | BstV2I | BstV2I | Eco130I | BstT1I |
| AsuC2I | AsuC2I | BseAI | Bsp13I | BspT107I | AccB1I | BstXI | BstXI | | ErhI |

Isoschizomers

| Enzyme | SE Enzyme | Enzyme | SE Enzyme | Enzyme | SE Enzyme | Enzyme | SE Enzyme | Enzyme | SE Enzyme |
|---------|-----------------------|-----------|-----------------------|----------|---------------------|----------|----------------------|----------|----------------------|
| Eco147I | PceI | HincII | HindII | MvaI | Ajnl [^] | Psp124BI | EcoICRI [^] | SnaI | BssNAI |
| EcoICRI | EcoICRI | HindII | HindII | | Bst2UI | | Psp124BI | SnaBI | BstSNI |
| | Psp124BI [^] | HindIII | HindIII | | Psp6I [^] | PspCI | PspCI | SpaHI | SphI |
| EcoNI | BstENI | Hinfl | Hinfl | Mva1269I | PctI | PspEI | PspEI | SpeI | AhII |
| EcoO65I | PspEI | HpaI | HpaI | MvnI | BstFNI | PspGI | Ajnl | SphI | SphI |
| EcoRI | EcoRI | HpaII | HpaII | MwoI | BstMWI | | Bst2UI [^] | | |
| EcoRII | Ajnl | | MspI | NaeI | MroNI [^] | | Psp6I | SplI | PspLI |
| | Bst2UI [^] | HphI | AsuHPI | NarI | Egel [^] | PspLI | PspLI | Sse9I | Sse9I |
| | Psp6I | HpyCH4III | Bst4CI | | Mly113I | PspN4I | PspN4I | Sse8387I | SbfI |
| EcoRV | EcoRV | HpyF10VI | BstMWI | | | PspOMI | ApaI [^] | SseBI | PceI |
| EcoT14I | BssT1I | HpySE526I | HpySE526I | NciI | AsuC2I | | PspOMI | Ssil | BspACI |
| | ErhI | Hsp92I | BstACI | NcoI | Bsp19I | PspPI | AspS9I | SspI | SspI |
| | Zsp2I | Hsp92II | FatI [^] | NdeI | FauNDI | PspPPI | PspPPI | SspBI | BstAUI |
| EcoT22I | FriOI | HspAI | AspLEI [^] | NdeII | BstMBI | PspXI | PspXI | SspMI | SspMI |
| EcoT38I | EgeI | | BstHHI [^] | | Kzo9I | Psrl | Psrl | SstI | EcoICRI [^] |
| EgeI | Mly113I [^] | | HspAI | | BstKTI [^] | PstI | PstI | | Psp124BI |
| EheI | EgeI | ItaI | Fsp4HI | NgoMIV | MroNI | Psul | BstX2I | StuI | PceI |
| | Mly113I [^] | KasI | Egel [^] | NheI | AsuNHI | PsyI | Tth111I | StyI | BstT1I |
| ErhI | BssT1I | | Mly113I ^{^^} | | BmtI [^] | PvuI | Ple19I | | ErhI |
| | ErhI | KpnI | KpnI | NlaIII | FaeI | PvuII | PvuII | StyD4I | BstSCI |
| EspI | Bsp1720I | | Acc65I [^] | | FatI [^] | RgaI | RgaI | | MspR9I [^] |
| FaeI | FaeI | Kpn2I | Bsp13I | NlaIV | PspN4I | RigI | RigI | SunI | PspLI |
| Fail | Fail | KroI | KroI | NotI | CciNI | RsaI | RsaI | SwaI | SmiI |
| Fall | Fall | KspI | Sfr303I | NruI | NruI | | RsaNI [^] | Taal | Bst4CI |
| FatI | FaeI [^] | Ksp22I | Ksp22I | NsbI | Acc16I | RsaNI | RsaNI | Taq I | TaqI |
| | FatI | Ksp632I | Bst6I | NsiI | Zsp2I | RsrII | Rsr2I | TasI | Sse9I |
| FauI | FauI | KspAI | HpaI | NspI | BstNSI | Rsr2I | Rsr2I | TelI | Tth111I |
| FauNDI | FauNDI | Kzo9I | BstMBI | NspIII | Ama87I | SacI | EcoICRI [^] | ThII | Sfr274I |
| FbaI | Ksp22I | | Kzo9I | NspV | Bpu14I | | Psp124BI | TruI I | Tru9I |
| FblI | FblI | | BstKTI [^] | NspBII | MspAII | SacII | Sfr303I | Tru9I | Tru9I |
| FinI | BsFI | Lmn I | Lmn I | PaeI | SphI | Sall | Sall | TseFI | TseFI |
| FnuDII | BspFNI | LweI | SfaNI | PaeR7I | Sfr274I | SapI | PciSI | Tsp45I | TseFI |
| | BstFNI | MabI | MabI | Pall | BsuRI | SatI | Fsp4HI | Tsp4CI | Bst4CI |
| Fnu4HI | Fsp4HI | Mae I | SspMI | | HaeIII | SauI | Bse21I | Tsp509I | Sse9I |
| FokI | FokI | MaeII | HpySE526I | Paul | BsePI | Sau96I | AspS9I | TspEI | Sse9I |
| | BstF5I [^] | Mall | Mall | PceI | PceI | Sau3AI | BstMBI | Tth111I | Tth111I |
| FriOI | FriOI | MamI | Bse8I | PciI | PciI | | Kzo9I | Van9II | AccB7I |
| FseI | RigI | Mbil | AccBSI | PciSI | PciSI | | BstKTI [^] | VneI | VneI |
| FspI | Acc16I | MboI | BstMBI | PesI | PesI | Sbfl | Sbfl | VpaK11BI | Bme18I |
| Fsp4HI | Fsp4HI | | Kzo9I | PctI | PctI | ScaI | Zrml | VspI | VspI |
| FunI | AfeI | | BstKTI [^] | PdII | MroNI [^] | SchI | PpsI [^] | XagI | BstENI |
| FunII | EcoRI | MboII | MboII | Pdml | MroXI | ScrFI | BstSCI [^] | XapI | AcsI |
| Gla I | Gla I | McrI | BstMCI | PfI23II | PspLI | | MspR9I | XbaI | XbaI |
| GluI | GluI | MfeI | MfeI | PfIBI | AccB7I | SdaI | Sbfl | XceI | BstNSI |
| GsaI | GsaI | MflI | BstX2I | PfIFI | Tth111I | SduI | MhII | XhoI | Sfr274I |
| GsuI | BpmI | MhII | MhII | PfIMI | AccB7I | SecI | BssECI | XhoII | BstX2I |
| HaeII | BstH2I | MlsI | Msp20I | PhoI | BsuRI | SetI | SetI | XmaI | SmaI [^] |
| HaeIII | BsuRI | MluI | MluI | | HaeIII | | | | XmaI |
| | HaeIII | MluNI | Msp20I | PinAI | AsiGI | SexAI | MabI | XmaIII | BseX3I |
| HapII | HpaII | MlyI | PpsI [^] | PkrI | PkrI | SfaNI | SfaNI | XmaCI | SmaI [^] |
| | MspI | Mly113I | Egel [^] | PleI | PpsI | SfcI | BstSFI | | XmaI |
| Hgal | Hgal | Mly113I | Mly113I | Ple19I | Ple19I | Sfel | BstSFI | XmaJI | AspA2I |
| HgiAI | Bbv12I | MnII | MnII | PmaCI | PspCI | Sfil | Sfil | Xmil | FbII |
| HgiCI | AccB1I | Mox20I | Mox20I | PmlI | PspCI | Sfol | Egel | Xmnl | MroXI |
| HgiJII | FriOI | Mph1103I | Zsp2I | PpsI | PpsI | | Mly113I [^] | ZhoI | Bsa29I |
| HhaI | AspLEI | MroI | Bsp13I | PpuMI | PspPPI | Sfr274I | Sfr274I | ZraI | AatII [^] |
| | BstHHI | MroNI | MroNI | PpuXI | PspPPI | Sfr303I | Sfr303I | | ZraI |
| | HspAI [^] | MroXI | MroXI | PshAI | BstPAI | Sful | Bpu14I | Zrml | Zrml |
| HinII | BstACI | MscI | Msp20I | PshBI | VspI | Sgfl | AsiSI | Zsp2I | Zsp2I |
| Hin6I | AspLEI [^] | MseI | Tru9I | Psil | Psil | | RgaI | | |
| | BstHHI [^] | MslI | SmiMI | Psp5II | PspPPI | SgrBI | Sfr303I | | |
| | HspAI | MspI | HpaII | Psp6I | Ajnl | SinI | Bme18I | | |
| HinPII | AspLEI [^] | | MspI | | Bst2UI [^] | SlaI | Sfr274I | | |
| | BstHHI [^] | MspA1I | MspA1I | | Psp6I | SmaI | SmaI | | |
| | HspAI | MspCI | BstAFI | Psp1406I | AclI | | XmaI [^] | | |
| | | MspR9I | BstSCI [^] | PspAI | SmaI [^] | SmiI | SmiI | | |
| | | | MspR9I | | XmaI | SmiMI | SmiMI | | |
| | | MstI | Acc16I | | | Smul | Faul | | |
| | | MteI | MteI | | | | | | |
| | | MunI | MfeI | | | | | | |

Alphabetized List of SE Recognition Sequences

| | | | | | |
|-----------------------------------|-----------|--------------------------------------|----------|--------------------------------------|-----------|
| AA^CGTT | Acl I | C^CGG | Hpa II | GACNN^NNGTC | BstPA I |
| A^AGCTT | Hind III | C^CGG | Msp I | GACNNN^NNGTC | Dri I |
| (8/13)AAGN ₅ CTT(13/8) | Fal I | CC^NGG | MspR9 I | GACNNNN^NNGTC | DseD I |
| AAT^ATT | Ssp I | ^CCNGG | BstSC I | (8/13)GAC(N) ₆ TTYG(11/6) | Ars I |
| ^AATT | Sse9 I | C^CNNGG | BssEC I | (5/4)GACTC | Pps I |
| A^CATGT | Pci I | CCNNNN^NNGG | Bsc4 I | (5/1)GAGAC | BstMA I |
| A^CCGGT | AsiG I | C^CRYGG | BstDS I | (5/1)GAGACC | Bso31 I |
| ACCTGC(4/8) | Acc36 I | CC^SGG | AsuC2 I | GAG^CGG | AccBS I |
| A^CCWGGT | Mab I | C^CTAGG | AspA2 I | GAG^CTC | EcoICR I |
| A^CGCGT | Mlu I | CCTC(7/6) | MnlI | GAGCT^C | Psp124B I |
| A^CGT | HpySE526I | CC^TCGAGG | Abs I | (6/7)GAGG | Mnl I |
| ACN^GT | Bst4C I | CCTGCA^GG | Sbf I | GAGTC(4/5) | Pps I |
| A^CTAGT | Ahl I | CC^TNAGC | Bpu10 I | G^ANTC | Hinf I |
| ACTGG(1/-1) | Bse1 I | CC^TNAGG | Bse21 I | GAT^ATC | EcoR V |
| ACTGGG(5/4) | Bmu I | CCTNN^NNNAGG | BstEN I | G(mA)^TC | Mal I |
| A^GATCT | Bgl II | ^CCWGG | Psp6 I | ^GATC | BstMB I |
| AGC^GCT | Afe I | CC^WGG | Bst2U I | ^GATC | Kzo9 I |
| AG^CT | Alu I | C^CWWGG | BssT1 I | GAT^C | BstKTI |
| AG^CT | AluB I | C^CWWGG | Erh I | (5/4)GATCC | AclW I |
| AGG^CCT | Pce I | CGAT^CG | Ple19 I | (9/5)GATGC | SfaN I |
| AGG^CGG | AccBS I | CG^CG | BspFN I | GATNN^NNATC | Bse8 I |
| AGT^ACT | Zrm I | CG^CG | BstFN I | GCAATG(2/0) | Bse3D I |
| ASST^ | Set I | C^GGCCG | BseX3 I | GCAGC(8/12) | BstV1 I |
| AT^CGAT | Bsa29 I | CG^GWCCG | Rsr2 I | (8/4)GCAGGT | Acc36 I |
| ATGCA^T | Zsp2 I | (5mC)GNNNNN^NN(5mC)G | Pes I | GCANNNN^NTGC | BstAP I |
| AT^TAAT | Vsp I | CGRY^CG | BstMC I | GCATC(5/9) | SfaN I |
| ATTT^AAAT | Smi I | C^GTACG | PspL I | GCATG^C | Sph I |
| C^AATTG | Mfe I | CMG^CKG | MspA I | (-1/1)GCATTC | Pct I |
| C^ACGAG | Bst2B I | C^TAG | SspM I | G^CCGGC | MroN I |
| CAC^GTC | Btr I | (14/16)CTCCAG | Bpm I | GCCNNNN^NGGC | Bgl I |
| CAC^GTG | PspC I | C^TCGAG | Sfr274 I | GCGAT^CGC | AsiS I |
| CACNNN^GTG | Dra III | C^TCGTG | Bst2B I | GCGAT^CGC | Rga I |
| CAG^CTG | Pvu II | CTCTTC(1/4) | Bst6 I | G^CGC | HspA I |
| CAGNNN^CTG | PstN I | CTGAAG(16/14) | Acu I | GCG^C | AspLE I |
| CA^TATG | FauND I | CTGCA^G | Pst I | GCG^C | BstHH I |
| (13/9)CATCC | Fok I | CTGGAG(16/14) | Bpm I | G^CGCGC | BseP I |
| (0/2)CATCC | BstF5 I | C^TNAG | BstDE I | G(5mC)G(5mC)^NG(5mC)G(5mC) | Mte I |
| CATG^ | Fae I | C^TRYAG | BstSF I | (-1/-3)GCGG | BspAC I |
| ^CATG | Fat I | C^TTAAG | BstAF I | GC^GGCCGC | CciN I |
| (0/2)CATTGC | Bse3D I | (14/16)CTTCAG | Acu I | (6/4)GCGGG | Fau I |
| CAYNN^NNRTG | SmiM I | C^YCGRG | Ama87 I | (10/5)GCGTC | Hga I |
| (-1/1)CCAGT | Bse1 I | (7/12)GAAC(N) ₆ TAC(12/7) | Psr I | GC^NGC | Fsp4H I |
| CCANNN^NTGG | AccB7 I | GAAGAC(2/6) | BstV2 I | G^C(5mC)GGC | Kro I |
| CCANNNNN^NTGG | BstX I | (4/1)GAAGAG | Bst6 I | G(5mC)^NGC | Bis I |
| C^CATGG | Bsp19 I | (4/1)GAAGAGC | PciS I | G(5mC)N^G(5mC) | Pkr I |
| CCCAG^C | Gsa I | (7/12)GAAG(N) ₆ TAC(12/7) | Bar I | G(5mC)^NG(5mC) | GluI |
| (4/5)CCCAGT | Bmu I | GAANN^NNTTC | MroX I | GCN^NGC | BstC8 I |
| CCCGC(4/6) | Fau I | GAATGC(1/-1) | Pct I | GNNNNN^NNGC | BstMW I |
| CCC^GGG | Sma I | G^AATTC | EcoR I | G^CTAGC | AsuNH I |
| C^CCGGG | Xma I | GACGC(5/10) | Hga I | GCTAG^C | Bmt I |
| CCGC(-3/-1) | BspAC I | GACGT^C | Aat II | GCTCCN^ | Lmn |
| CCG^CCT | AccBS I | GAC^GTC | Zra I | GCTCTTC(1/4) | PciS I |
| CCGC^GG | Sfr303 I | GAC^GTG | Btr I | (12/8)GCTGC | BstV1 I |
| | | GACN^NNGTC | Tth111 I | | |

R = A or G

K = G or T

D = A or G or T

W = A or T

M = A or C

H = A or C or T

S = G or C

Y = T or C

V = A or C or G

B = C or G or T

N = A or C or G or T

Alphabetized List of SE Recognition Sequences

| | | | | | |
|-----------------------------|-----------|--------------------------------------|----------|------------------------|---------|
| G [^] CTGGG | Gsa I | G [^] GYRCC | AccB1 I | R [^] GATCY | BstX2 I |
| GC [^] TNAGC | Bst1720 I | GKGC [^] M [^] C | BstSL I | RGCGC [^] Y | BstH2 I |
| GC [^] TNAGG | Bpu10 I | GR [^] CGYC | BstAC I | RG [^] GWCCY | PspPP I |
| GDGCH [^] C | Mhl I | GRGCY [^] C | FriO I | TAC [^] GTA | BstSN I |
| GGATC(4/5) | AclW I | GT [^] AC | Rsa I | (7/8)TCACC | AsuHP I |
| G [^] GATCC | BamH I | G [^] TAC | RsaN I | T [^] CATGA | Cci I |
| GGATG(2/0) | BstF5 I | (7/12)GTA(N) ₆ CTTC(12/7) | Bar I | T [^] CCGGA | Bsp13 I |
| GGATG(9/13) | Fok I | (7/12)GTA(N) ₆ GTTC(12/7) | Psr I | T [^] CGA | Taq I |
| GG [^] CC | BsuR I | GTA [^] TAC | BssNA I | TCG [^] CGA | Nru I |
| GG [^] CC | Hae III | GTATCC(6/5) | Bsu I | T [^] CTAGA | Xba I |
| GGCCGG [^] CC | Rig I | (14/10)GTCCC | BsI F I | T [^] GATCA | Ksp22 I |
| GGCCNNNN [^] NGGCC | Sfi I | G [^] TCGAC | Sal I | TGC [^] GCA | Acc16 I |
| GG [^] CGCC | Mly113 I | GTCTC(1/5) | BstMA I | TGG [^] CCA | Mox20 I |
| GGC [^] GCC | Ege I | (6/2)GTCTTC | BstV2 I | T [^] GTACA | BstAU I |
| GG [^] CGCGCC | PalA I | G [^] TGCAC | Vne I | T [^] TAA | Tru9 I |
| GGGAC(10/14) | BsI F I | GT [^] MKAC | Fbl I | TTA [^] TAA | Psi I |
| G [^] GGCCC | PspOM I | [^] GTSAC | TseF I | TT [^] CGAA | Bpu14 I |
| GGGCC [^] C | Apa I | GTT [^] AAC | Hpa I | TTS [^] AA | Ags I |
| G [^] GNCC | AspS9 I | GTY [^] RAC | Hind II | TTT [^] AAA | Dra I |
| GGN [^] NCC | PspN4 I | GWGCW [^] C | Bbv12 I | VC [^] TCGAGB | PspX I |
| G [^] GTACC | Acc65 I | Pu(5mC) [^] GPy | Gla I | YAC [^] GTR | BstBA I |
| GGTAC [^] C | Kpn I | [^] PuG(5mC)Py | Aox I | YA [^] TR | Fai I |
| GGTCTC(1/5) | Bso31 I | PuPyN [^] PuPy | Bls I | Y [^] GGCCR | Aco I |
| GGTGA(8/7) | AsuHP I | R [^] AATTY | Acs I | | |
| G [^] GTNACC | PspE I | RCATG [^] Y | BstNS I | | |
| G [^] GWCC | Bme18 I | R [^] CCGGY | Bse118 I | | |

R = A or G

K = G or T

D = A or G or T

W = A or T

M = A or C

H = A or C or T

S = G or C

Y = T or C

V = A or C or G

B = C or G or T

N = A or C or G or T